



[Manuals.plus](#) /

› [Golden Technologies](#) /

› Golden Technologies GB120 Buzzaround Carry-On Folding Travel Scooter User Manual

Golden Technologies GB120

Golden Technologies GB120 Buzzaround Carry-On Folding Travel Scooter User Manual

Model: GB120

Brand: Golden Technologies

INTRODUCTION

This manual provides essential information for the safe and efficient operation, maintenance, and care of your Golden Technologies GB120 Buzzaround Carry-On Folding Travel Scooter. Please read this manual thoroughly before using your scooter to ensure proper understanding of its features and functions. The GB120 Buzzaround Carry-On is a compact, portable, and stylish mobility scooter designed for ease of transport and convenience, offering a 300 lb weight capacity, an adjustable tiller, and integrated lighting.

SAFETY INFORMATION

Your safety is paramount. Adhere to all warnings and instructions to prevent injury or damage.

- **Maximum User Weight:** Do not exceed the maximum user weight of 300 lbs (136 kg). Exceeding this limit can compromise stability and performance.
- **Ground Clearance:** The scooter has a ground clearance of 4 inches (10 cm). Avoid operating on uneven or rough terrain that may cause the scooter to become unstable or stuck.
- **Battery Safety:** Always follow the battery charging and handling instructions provided. Use only the approved charger and battery pack. Do not expose the battery to extreme temperatures or puncture it. The included 12.5 AH lithium-ion battery is airline-approved.
- **Night Operation:** When operating the scooter at night or in low-light conditions, always use the integrated front and rear lights for enhanced visibility to others and to illuminate your path.
- **Braking:** The scooter is equipped with electro-mechanical brakes. Familiarize yourself with their operation before use.
- **Disassembly/Assembly:** Exercise caution when disassembling or assembling the scooter. Ensure all components are securely locked in place before operation.

PRODUCT FEATURES

The Golden Technologies GB120 Buzzaround Carry-On scooter is designed with several features to enhance your mobility and convenience:

- **Compact Storage:** Ideal for space-saving storage due to its fold-flat design.
- **Folding Design:** Features a folding tiller and a self-locking seat frame for a compact size, folding to 12 inches tall (without the seat).
- **Portability:** Equipped with four carry handles for easy lifting and transport.
- **Weight Capacity:** Supports users up to 300 lbs (136 kg).
- **Adjustable Tiller:** Infinitely adjustable tiller for personalized comfort, includes a USB port and a water bottle holder.
- **Enhanced Visibility:** Full front and rear lights for safety during operation.
- **Convenient Armrests:** Flip-up armrests for easier access and exit.
- **Battery:** Standard 12.5 AH lithium-ion battery pack with easy access, offering an operating range of 14 miles. This battery is airline-approved. An optional long-range 15 AH battery (18-mile range) is available separately.
- **Disassembly:** Easily disassembles into two main pieces for transport.

SETUP

1. **Unpacking:** Carefully remove all components from the packaging. Inspect for any shipping damage.
2. **Unfolding the Scooter:**
 - Place the scooter on a flat, stable surface.
 - Release any locking mechanisms that secure the scooter in its folded position.
 - Gently unfold the main frame until it locks into the operational position.
 - Raise the tiller to its upright position and secure it.
3. **Attaching the Seat:**
 - Align the seat post with the receiver on the scooter's frame.
 - Insert the seat post firmly until it clicks into place.
 - Ensure the seat is stable and does not wobble.
4. **Battery Installation/Charging:**
 - Locate the battery compartment, typically under the floorboard.
 - Insert the 12.5 AH lithium-ion battery pack, ensuring it is properly seated and connected.
 - Before first use, fully charge the battery using the provided off-board charger. Refer to the charger's instructions for specific charging times and indicators.
5. **Tiller Adjustment:** Adjust the tiller to a comfortable height and angle using the adjustment lever. Ensure it is securely locked before operating.
6. **Armrest Adjustment:** The armrests can be flipped up for easier access. Adjust them to your preferred position if they offer height or width adjustments.

OPERATING INSTRUCTIONS

1. **Power On/Off:** Turn the key in the ignition to the "ON" position to power on the scooter. Turn to "OFF" to power off.
2. **Speed Control:** Use the speed adjustment dial or lever on the tiller to set your desired maximum speed. Start at a low speed until you are comfortable with the scooter's handling.

3. **Driving:**

- To move forward, gently push the throttle lever forward.
- To move backward, gently pull the throttle lever backward.
- Release the throttle lever to engage the automatic electro-mechanical brakes and stop the scooter.

4. **Steering:** Use the tiller to steer the front wheels. Turn gradually, especially at higher speeds, to maintain stability.

5. **Lights:** Activate the front and rear lights using the designated switch on the tiller for improved visibility.

6. **Horn:** Use the horn button on the tiller to alert others when necessary.

7. **USB Port:** The tiller includes a USB port for charging small electronic devices.

8. **Freewheel Mode:** If equipped, use the freewheel lever (usually located near the rear wheels) to disengage the motor and allow the scooter to be pushed manually. Always ensure the scooter is in drive mode before operating.

MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your scooter.

- **Cleaning:** Clean the scooter regularly with a damp cloth. Avoid using abrasive cleaners or direct water spray.
- **Battery Care:**
 - Charge the battery after each use, even if only for a short period.
 - If storing the scooter for an extended period, charge the battery fully and then every month to maintain its health.
 - Do not allow the battery to fully discharge.
- **Tire Inspection:** Periodically check tire pressure (if applicable) and for any signs of wear or damage.
- **Fasteners:** Ensure all nuts, bolts, and fasteners are tight. Do not overtighten.
- **Professional Service:** For any complex repairs or issues, contact an authorized Golden Technologies service center.

TROUBLESHOOTING

If you encounter issues with your scooter, refer to the following common troubleshooting steps:

- **Scooter Does Not Power On:**
 - Check if the key is fully inserted and turned to the "ON" position.
 - Ensure the battery is properly installed and charged.
 - Verify that all electrical connections are secure.
- **Scooter Moves Slowly or Not at All:**
 - Check the battery charge level. A low battery will reduce performance.
 - Ensure the speed control dial is not set to the lowest setting.
 - Verify the freewheel lever is in the "drive" position (motor engaged).
- **Unusual Noises:**
 - Stop the scooter immediately.
 - Inspect for any loose parts or obstructions.
 - If the noise persists, contact customer support.
- **Lights Not Working:**

- Check the light switch position.
- Inspect for any visible damage to the light bulbs or wiring.

For issues not covered here, please contact Golden Technologies customer support.

SPECIFICATIONS

Feature	Detail
Model Name	Buzzaround CarryOn
Manufacturer Part Number	GB120
Maximum User Weight	300 lbs (136 kg)
Maximum Speed	3.6 Miles per Hour (5.8 km/h)
Operating Range (Standard Battery)	14 miles (22.5 km)
Battery Type	12.5 AH Lithium-Ion (Airline Approved)
Folded Height (without seat)	12 inches (30.5 cm)
Ground Clearance	4 inches (10 cm)
Turn Radius	48.75 Inches (123.8 cm)
Brake Style	Electromagnetic
Wheel Type	Four-Wheeled
Material Type	Aluminum, Vinyl (Frame: Alloy Steel)
Additional Features	50 Amp Controller, Electronic Speed Control, Freewheel Mode, USB Port, Water Bottle Holder

VISUAL GUIDE

Product Images



Side view of the Golden Tech GB120 Buzzaround Carry-On Folding Travel Scooter in Crisp White, showcasing its compact design and comfortable seat.



Another side view of the scooter, highlighting the front basket and the adjustable tiller.



Front view of the scooter, showing the integrated front lights and the tiller controls.



Rear view of the scooter, displaying the rear lights and the back of the seat.



The scooter in its folded upright position, demonstrating its compact storage capability.



The scooter folded flat, illustrating how it can be stored under a bed or in tight spaces.



View of the scooter disassembled into its main components: seat, battery, basket, and main frame, for easy transport.



Standard Colors:



Velvet Red
(Installed)



Coastal Blue
(Included in Box)

NEW Sporty Optional Colors:



Crisp White



Canary Yellow



Stone Grey



Orchid Pink



Army Green
(Matte)



Desert Tan

Buzzaround GB120



The Buzzaround CarryOn fold-flat scooter is small in stature and big on fun! This super portable, sporty scooter will take you where you want to go in comfort and style! It features a 300 lb. weight capacity, an infinite adjustable tiller with a convenient USB port and drink holder, full front and rear light package, and flip up armrests. Enjoy long range coverage up to 14 miles on the airline friendly 12.5 AH lithium-ion battery, which can be removed and taken as a carry on while traveling by air, sea or train in a stylish carry case. A new rear basket is mounted on the deck to hold an optional spare 12.5 AH battery for even longer range and serves as a kickstand for storing upright. The heaviest piece is 47.5 lbs. (without seat and batteries) and the height when folded is only 12". It disassembles easily into two pieces and is designed to store in the smallest compartments. It comes standard with red color shrouds installed and an extra set in blue. New optional color shroud kits let you express a style all your own! Pack your dreams and go anywhere with the CarryOn Scooter!



A lifestyle image showing the Buzzaround CarryOn scooter in use, along with various color options available for customization.

**NEW & IMPROVED
Sporty Fold-Flat Scooter!
Airline Friendly!
USB Port Standard!**

Standard Features:

- 4-point carry handles for easy transport
- 12" height when folded (without seat)
- 4" ground clearance
- 300 lb. weight capacity
- Flip up armrests
- Infinite adjustable tiller
- USB port & drink holder
- Disassembles easily
- Up to 14 mile range
- Easy access, airline friendly 12.5 AH lithium ion battery pack

Official Product Videos

No official seller videos are available for embedding at this time.

WARRANTY AND SUPPORT

For warranty information, please refer to the documentation included with your purchase or visit the official Golden Technologies website. For technical support, spare parts, or service inquiries, please contact Golden Technologies customer service directly. Keep your purchase receipt and model number handy when contacting support.

Contact Information: Please refer to the contact details provided in your product packaging or on the manufacturer's official website.