

## VERPAL Compact 4 Wheel Mobility Scooter

# VERPAL Compact 4-Wheel Foldable Mobility Scooter User Manual

Model: Compact 4 Wheel Mobility Scooter | Brand: VERPAL

## 1. INTRODUCTION

This manual provides essential information for the safe and efficient operation, assembly, and maintenance of your VERPAL Compact 4-Wheel Foldable Mobility Scooter. Please read it thoroughly before first use and keep it for future reference.

## 2. SAFETY INFORMATION

- The scooter is ISO safety-certified.
- Features anti-rollover rear wheels for enhanced safety.
- Always read this instruction manual carefully before use to ensure proper operation and safety.
- Limit usage to one person at a time.
- Avoid steep slopes and uneven terrain.
- Use rearview mirrors for improved visibility and safety.
- Park the scooter in a stable and upright position when not in use.

## 3. SETUP AND ASSEMBLY

1. **Assemble the Scooter:** Attach the four lightweight pieces together. Ensure the folding tiller and removable seat are securely fastened. Regularly check the tightness of all screws.
2. **Adjust the Scooter:** Set the seat height using the simple lever for a comfortable fit. The padded armrests can be flipped up for easy access.

*Video: Product Installation Instructions. This video demonstrates the assembly process of the mobility scooter, including attaching the main components.*

*Video: ENGWE EASE1 Mobility Scooter Unboxing & Assembly. This video provides a visual guide to unboxing and assembling a similar mobility scooter model.*



*Image: Intuitive Control Panel. This image highlights the user-friendly control panel with clearly labeled buttons for speed, forward/reverse, headlights, and horn.*

## 4. OPERATING INSTRUCTIONS

---

1. **Power On:** Insert the key into the ignition and turn it to the 'ON' position.
2. **Adjust Speed:** Use the variable speed dial on the control panel to set your desired maximum speed.
3. **Direction Control:** Use the throttle input to move forward or backward. The intuitive control panel allows for smooth transitions.
4. **Headlights and Horn:** Turn on the LED headlights for visibility, especially in low-light conditions. Use the horn as needed for safety.



## Headlights

Illuminate Your Path



## Warning Tail Lights

Stay Visible, Stay Safe

*Image: Headlights and Warning Tail Lights. This image shows the scooter's LED headlights for illuminating your path and warning tail lights for visibility.*

*Video: Mobility Scooter for Adults and Seniors. This video demonstrates the scooter in use, showcasing its maneuverability and features.*

## 5. MAINTENANCE

---

- **Cleaning:** Regularly clean the scooter with a damp cloth and mild detergent.
- **Tire Pressure:** Check tire pressure before each use.
- **Battery Care:** Check battery life before each use. Charge the battery after each use to maintain optimal performance and longevity.

## 6. TROUBLESHOOTING

---

If you encounter any issues with your mobility scooter, please refer to the following video for common alarm

issues and potential solutions.

*Video: Common Alarm Issues with Mobility Scooters. This video addresses frequently encountered alarm problems and provides guidance on how to resolve them.*

## 7. SPECIFICATIONS

<b>Product Dimensions</b>	38"D x 25"W x 37"H
<b>Item Weight</b>	92 Pounds
<b>Weight Capacity</b>	Up to 265 pounds
<b>Top Speed</b>	Up to 3.7 mph
<b>Range Per Charge</b>	Up to 12.4 miles
<b>Motor Power</b>	250W Peak Power
<b>Battery</b>	Long-Lasting Extended
<b>Style</b>	Compact; Foldable; Motorized

# Product Dimension



*Image: Product Dimensions. This image provides a visual representation of the scooter's dimensions and weight capacity.*

## 8. WARRANTY AND SUPPORT

VERPAL offers lifetime support for your peace of mind. For any inquiries or assistance, please contact our customer service.