

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

Manuals.plus /

› [TROPOW](#) /

› TROPOW Scientific Lab Incubator 20L, RT-300°C Digital PID Control - Instruction Manual

TROPOW 20L, RT-300°C

TROPOW Scientific Lab Incubator 20L, RT-300°C Digital PID Control - Instruction Manual

Model: 20L, RT-300°C

1. INTRODUCTION

Thank you for choosing the TROPOW Scientific Lab Incubator. This manual provides essential information for the safe and efficient operation, installation, and maintenance of your new incubator. Please read this manual thoroughly before operating the device and retain it for future reference. This incubator is designed for precise temperature control in various laboratory applications, including agricultural, animal husbandry research, and microbiological studies.

2. SAFETY INSTRUCTIONS

To ensure safe operation and prevent damage to the equipment or injury to personnel, observe the following safety precautions:

- Always connect the incubator to a grounded power outlet.
- Do not operate the incubator in environments with flammable or explosive gases.
- Ensure adequate ventilation around the unit to prevent overheating.
- Avoid placing objects on top of the incubator that could obstruct ventilation.
- Do not open the door frequently during operation to maintain stable internal temperature.
- Wear appropriate personal protective equipment (PPE) when handling samples or performing maintenance.
- Disconnect power before cleaning or performing any maintenance.
- In case of malfunction, switch off the unit and consult the troubleshooting section or contact customer support.



Figure 2.1: Front view of the incubator, showing general appearance and location of safety warning labels.

3. PRODUCT OVERVIEW

3.1 Components

The TROPOW Scientific Lab Incubator includes the following main components:

- **Incubator Unit:** The main body of the incubator.
- **Adjustable Shelves:** For organizing samples within the chamber.
- **Power Cord:** For connecting to an electrical outlet.
- **Instruction Manual:** This document.



Figure 3.1: Incubator dimensions. External dimensions are approximately 20.08 x 18.9 x 24.8 inches. Internal chamber dimensions are approximately 11.8 x 11.8 x 9.8 inches.

3.2 Key Features

The incubator is equipped with several features designed for optimal performance and user convenience:

- **High-Quality Construction:** Features a cold-rolled steel plate exterior with electrostatic spraying and a 304 stainless steel inner tank for durability and corrosion resistance.
- **Intelligent PID Control:** Ensures precise and stable temperature regulation, with an overtemperature alarm and temperature deviation correction capabilities.
- **Efficient Hot Air Circulation:** A fan-assisted system provides uniform heating throughout the chamber, reaching temperatures up to 300°C.
- **Observation Window:** A tempered glass window allows for visual monitoring of samples without opening the door.
- **Enhanced Sealing:** Ultra-thickened structure and double silicone seals minimize heat loss and maintain temperature stability.
- **Exhaust Holes:** Located on the top, these allow for controlled thermal discharge.

- **Ergonomic Handle:** A built-in handle facilitates easy opening and closing of the door.



Intelligent PID Control

Ensure temperature stability in the incubator, precise temperature control

Figure 3.2: The intelligent PID control panel for precise temperature management.



Multiple Application Scenarios

A perfect thermostat for the cultivation of bacteria and microorganisms

Figure 3.3: Labeled components illustrating the transparent door, exhaust holes, silicone sealing strips, and built-in handle.

4. SETUP

1. **Unpacking:** Carefully remove the incubator from its packaging. Inspect for any shipping damage.
2. **Placement:** Place the incubator on a stable, level surface away from direct sunlight, heat sources, and strong vibrations. Ensure there is adequate clearance (at least 10 cm) around all sides for proper ventilation.
3. **Power Connection:** Connect the power cord to the incubator's power inlet and then to a grounded electrical outlet with the correct voltage (refer to specifications).
4. **Shelf Installation:** Install the adjustable shelves at the desired height within the chamber.
5. **Initial Cleaning:** Before first use, wipe the interior of the chamber with a mild disinfectant solution and a soft cloth.

5. OPERATING INSTRUCTIONS

1. **Power On:** Flip the main power switch to the "ON" position. The digital display will illuminate.

2. Setting Temperature:

- Press the "SET" button on the PID controller. The set temperature value will flash.
- Use the "UP" and "DOWN" arrow buttons to adjust the desired temperature.
- Press "SET" again to confirm the new temperature. The incubator will begin heating or cooling to reach the set point.

3. **Loading Samples:** Place samples evenly on the shelves, ensuring not to overload the chamber or block the air circulation vents.

4. **Monitoring:** Observe the current temperature on the digital display. The PID controller will maintain the set temperature with minimal fluctuation.

5. **Overtemperature Alarm:** If the temperature exceeds a preset safety limit, an alarm will sound, and the heating element may shut off automatically. Address the cause of the overtemperature condition.

6. **Power Off:** When incubation is complete, remove samples, and then turn off the main power switch.



Thermal Convection System

Provide a constant temperature incubation environment

Figure 5.1: Illustration of the thermal convection system providing uniform heating within the chamber.

6. MAINTENANCE

Regular maintenance ensures the longevity and optimal performance of your incubator.

- **Cleaning Exterior:** Wipe the exterior surfaces with a soft cloth dampened with a mild detergent. Avoid abrasive cleaners.
- **Cleaning Interior:** For the interior, use a laboratory-grade disinfectant. Ensure the chamber is dry before resuming operation. Remove shelves for thorough cleaning if necessary.
- **Shelf Maintenance:** Clean shelves regularly to prevent contamination.
- **Door Seal Inspection:** Periodically check the silicone door seals for cracks or damage. Damaged seals can affect temperature stability.
- **Ventilation Check:** Ensure that the ventilation openings are clear of dust and debris.

7. TROUBLESHOOTING

Refer to the table below for common issues and their potential solutions.

Problem	Possible Cause	Solution
Incubator does not power on.	No power supply; power cord loose; fuse blown.	Check power connection; ensure outlet is functional; inspect fuse.
Temperature not reaching set point.	Door not sealed properly; ambient temperature too low/high; heating element malfunction.	Check door seal; ensure proper room temperature; contact support if heating element is suspected.
Temperature fluctuates significantly.	Frequent door opening; poor door seal; PID parameters require adjustment.	Minimize door openings; inspect door seal; consult advanced settings in the full manual or contact support.
Overtemperature alarm sounds.	Sensor error; control system malfunction; external heat source.	Power off and restart; ensure proper ventilation; contact support.

8. SPECIFICATIONS

Parameter	Value
Model	20L, RT-300°C
Capacity	20 Liters
Temperature Range	Room Temperature to 300°C
Control System	Digital PID Control
Inner Tank Material	304 Stainless Steel
Exterior Material	Cold-rolled Steel Plate with Electrostatic Spraying
Product Dimensions (L x W x H)	20.08 x 18.9 x 24.8 inches
Weight	56.27 Pounds
Manufacturer	TROPOW

9. WARRANTY AND SUPPORT

For warranty information, technical support, or service inquiries, please contact your vendor or the manufacturer directly. Keep your purchase receipt and product serial number handy when contacting support.

Manufacturer: TROPOW

ASIN: B0DF53SFZK