



Manuals.plus /

› ByteTecpeak /

› ByteTecpeak M4 Lightweight Folding Wheelchair User Manual

## ByteTecpeak M4

# ByteTecpeak M4 Lightweight Folding Wheelchair User Manual

Model: M4

[Introduction](#)

[Setup](#)

[Operation](#)

[Maintenance](#)

[Troubleshooting](#)

[Specifications](#)

[Support](#)

## 1. INTRODUCTION

---

Thank you for choosing the ByteTecpeak M4 Lightweight Folding Wheelchair. This manual provides essential information for the safe and effective use, assembly, operation, and maintenance of your new wheelchair. Please read it thoroughly before first use and keep it for future reference.

The M4 wheelchair is designed for adults, offering a lightweight and foldable solution for mobility. Constructed from high-quality aluminum alloy, it combines durability with ease of transport, making it suitable for both home use and travel. It features a maximum capacity of 100kg.



The ByteTecpeak M4 Lightweight Folding Wheelchair in its open configuration.

## 2. SETUP AND UNFOLDING

---

The ByteTecpeak M4 wheelchair is designed for quick and easy setup. Follow these steps to unfold and prepare your wheelchair for use.

1. **Unfolding the Frame:** Gently pull the backrest upwards while pushing down on the seat rails until the frame locks into place. Ensure all locking mechanisms are secure.
2. **Attaching Footrests:** The footrests are adjustable and removable. Align the footrest pins with the corresponding holes on the front frame and push down until they click into place. Adjust the height as needed.
3. **Securing Cushions:** Place the double-layer breathable seat cushion onto the seat. Ensure it is properly aligned and secured, typically with Velcro straps on the underside.
4. **Adjusting Backrest:** Ensure the backrest is fully extended and secured.

Your browser does not support the video tag.

Video demonstrating the unfolding, footrest attachment, and cushion placement for the ByteTecpeak M4 wheelchair.

## Facile da piegare, facile da trasportare

Dimensioni adatte per il bagagliaio  
Sufficientemente leggero da essere portato su scale  
da una persona



Pezzo di bagaglio



Mani in alto

## Meccanismo di piegatura rapida

Si piega in meno di 30 secondi con un  
sollevamento facile, movimento insieme



Quick folding mechanism for easy storage and transport.

## 3. OPERATING INSTRUCTIONS

### 3.1 User Propulsion

The M4 wheelchair features a self-propelled design with smooth-rotating push rims for effortless control. The ergonomic push wheels are optimized for easy self-propulsion across various surfaces. The lightweight frame reduces the energy required for autonomous navigation.

## Flessibile: Guidare da soli o essere spinti

Spingere a mano



Altri spingono



User propelling the wheelchair independently.

### 3.2 Assisted Propulsion

For assisted use, the wheelchair can be easily pushed by a caregiver. The design ensures comfortable handling for both the user and the assistant.

# Sedia Più Larga Per Maggiore Comfort

+ Ideale Per Gli Anziani

+ Design Ergonomico

+ Comodo Per Diversi Tipi Di Corpi



Caregiver assisting with wheelchair propulsion.

## 3.3 Braking System

The M4 wheelchair is equipped with a dual front and rear braking system for maximum safety and control. Engage the brakes when stationary to prevent unintended movement.

## 3.4 Comfort Features

- **Wide Seat:** The extra-wide seat (46cm) provides enhanced comfort and freedom of movement, reducing discomfort during prolonged use.
- **Breathable Cushion:** The double-layer breathable cushion, made from mesh fabric and synthetic leather, ensures comfort in all seasons.
- **Removable Armrests:** Armrests can be easily removed for easier transfers or to fit under tables.
- **Storage Pocket:** A practical storage pocket is located on the backrest for water bottles, magazines, and small items.

# Conservazione pratica sulla schiena

Tasca per bottiglie d'acqua & riviste & piccoli oggetti



Convenient storage pocket on the back of the wheelchair.

## 4. MAINTENANCE

---

Regular maintenance ensures the longevity and safe operation of your ByteTecpeak M4 wheelchair.

- **Cleaning:** Clean the frame with a damp cloth and mild soap. Avoid abrasive cleaners. The cushions can be wiped clean.
- **Wheels:** The solid rubber wheels are maintenance-free and do not require inflation. Regularly check for debris or damage.
- **Brakes:** Periodically check the brake system for proper function. Ensure the brakes engage securely.
- **Fasteners:** Inspect all bolts, screws, and fasteners regularly to ensure they are tight. Tighten any loose components.
- **Storage:** When not in use, store the wheelchair in a dry, clean place, ideally folded to save space.

## 5. TROUBLESHOOTING

---

This section addresses common issues you might encounter with your M4 wheelchair.

Problem	Possible Cause	Solution
Wheelchair does not unfold or fold smoothly.	Frame components are misaligned or obstructed.	Ensure no objects are obstructing the folding mechanism. Apply gentle, even pressure when unfolding/folding. Check for bent parts.
Brakes are not holding securely.	Brake mechanism needs adjustment or is worn.	Inspect the brake levers and calipers. If adjustment is needed or parts are worn, contact customer support.
Footrests are loose or difficult to attach.	Pins or receiving holes are dirty or damaged.	Clean the pins and holes. Ensure correct alignment before pushing to lock. If damaged, contact customer support.

If you encounter an issue not listed here, or if the suggested solutions do not resolve the problem, please contact ByteTecpeak customer support.

## 6. SPECIFICATIONS

Feature	Detail
Brand	ByteTecpeak
Model	M4
Product Dimensions (Open)	93 cm (L) x 63 cm (W) x 101 cm (H)
Product Dimensions (Folded)	Approx. 72 cm (L) x 36 cm (W) x 70 cm (H)
Item Weight	12 kg
Maximum Load Capacity	100 kg
Frame Material	High-quality Aluminum Alloy
Seat Width	46 cm
Rear Wheel Size	20 inches (approx. 51 cm)
Front Wheel Size	Approx. 16 cm
Wheel Type	Solid Rubber (Maintenance-free)
Braking System	Dual Front and Rear Brakes
Country of Origin	China

# Dimensioni pratiche per ogni esigenza



Detailed dimensions of the M4 wheelchair.

## 7. WARRANTY AND SUPPORT

ByteTecpeak is committed to customer satisfaction. Your M4 wheelchair comes with a standard return policy allowing returns and exchanges within 30 days of purchase.

For any questions, concerns, or assistance, please contact ByteTecpeak customer support. We offer 24/7 customer support and multi-language assistance to ensure you receive the help you need.

Please refer to the product packaging for any additional warnings or instructions specific to your purchase.



ByteTecpeak customer support is available to assist you.