

[Manuals.plus](#) /

› [Beurer](#) /

› Beurer IH 17 Nebulizer User Manual

Beurer Beurer IH 17 ph

Beurer IH 17 Nebulizer User Manual

Model: Beurer IH 17 ph

Brand: Beurer

PRODUCT OVERVIEW

The Beurer IH 17 Nebulizer is designed for the treatment of the upper and lower airways, including conditions such as colds, asthma, and various respiratory diseases. It features a high proportion of respirable particles, ensuring effective delivery of medication.

This device is effective in relieving symptoms associated with Asthma, COPD, Bronchitis, and Pneumonia, and is designed for short inhalation times.



Figure 1: The Beurer IH 17 Nebulizer main unit, showing the compact compressor, air tube, and nebulizer cup with mouthpiece attached.

INCLUDED COMPONENTS

The Beurer IH 17 Nebulizer package includes the following essential components:

- Nebulizer compressor unit
- Air tube
- Atomizer/Nebulizer cup
- Adult mask

- Child/Baby mask
- Mouthpiece
- Nosepiece
- Replacement filters (typically 10)



Figure 2: A comprehensive view of the accessories provided with the nebulizer, including adult and baby masks with straps, the atomizer, air hose, and a set of filters.

SETUP INSTRUCTIONS

1. **Unpack Components:** Carefully remove all components from the packaging and inspect them for any damage.
2. **Place the Unit:** Position the nebulizer compressor unit on a stable, flat surface, ensuring adequate ventilation around the device.
3. **Connect Air Tube:** Attach one end of the air tube securely to the air outlet on the compressor unit.
4. **Prepare Atomizer:** Open the atomizer/nebulizer cup and pour the prescribed medication into the medication chamber. Close the atomizer securely.
5. **Connect Atomizer:** Attach the other end of the air tube to the inlet port on the atomizer.
6. **Attach Accessory:** Connect the desired accessory (mouthpiece, adult mask, or child mask) to the top of the atomizer.
7. **Connect Power:** Plug the power cord into a suitable electrical outlet.

OPERATING THE NEBULIZER

1. **Power On:** Press the power button on the nebulizer unit to turn it on. The compressor will begin to generate airflow, and the medication will start to nebulize into a fine mist.
2. **Inhale Medication:**
 - If using the mouthpiece: Place the mouthpiece between your teeth and seal your lips around it. Breathe in slowly and deeply through your mouth until all medication is consumed.
 - If using a mask: Place the mask over your nose and mouth, ensuring a snug fit. Breathe normally through the mask.
3. **Duration:** Continue the inhalation until the medication chamber is empty, or as directed by your healthcare professional. The Beurer IH 17 is designed for short inhalation times.
4. **Power Off:** Once the treatment is complete, press the power button to turn off the nebulizer. Disconnect the power cord from the outlet.



Figure 3: A user demonstrating the use of the nebulizer with a mask, emphasizing the convenience of short inhalation times.



Figure 4: A user holding the nebulizer atomizer, illustrating the product's effectiveness in relieving symptoms of Asthma, COPD, Bronchitis, and Pneumonia.

MAINTENANCE AND CLEANING

Regular cleaning and maintenance are crucial for the optimal performance and longevity of your nebulizer. Always disconnect the device from the power supply before cleaning.

After Each Use:

- **Atomizer and Accessories:** Disassemble the atomizer, mouthpiece, and masks. Wash these components thoroughly with warm water and a mild detergent. Rinse them well under running water.
- **Disinfection:** For disinfection, boil the washed components in water for 5 minutes, or use a cold sterilization solution as recommended by your healthcare provider. Allow all parts to air dry completely on a clean cloth or paper towel.
- **Air Tube:** The air tube does not require cleaning unless visibly soiled. If cleaning is necessary, wipe the exterior with a damp cloth. Do not immerse the air tube in water.

Compressor Unit:

- Wipe the exterior of the compressor unit with a clean, damp cloth. Do not use abrasive cleaners or immerse the unit in water.

Filter Replacement:

- Check the air filter regularly. Replace it when it appears discolored or every 100-120 operating hours, or as indicated in the full product manual.

Storage:

- Store the nebulizer and its accessories in a clean, dry place, away from direct sunlight and extreme temperatures.

TROUBLESHOOTING

If you encounter issues with your Beurer IH 17 Nebulizer, refer to the following common troubleshooting steps:

Problem	Possible Cause	Solution
No mist/low mist production	Atomizer not assembled correctly Air tube kinked or disconnected Medication chamber empty or insufficient medication Air filter clogged	Reassemble atomizer ensuring all parts are secure. Check air tube for kinks or disconnections; reattach securely. Add prescribed medication. Replace air filter.
Unit not turning on	Power cord not plugged in Power outlet not working	Ensure power cord is securely plugged into the unit and outlet. Test the outlet with another device.
Unusual noise from compressor	Unit not on a stable surface Internal issue	Place the unit on a firm, level surface. If noise persists, contact customer support.

If the problem persists after attempting these solutions, please contact Beurer customer support for further assistance.

TECHNICAL SPECIFICATIONS

Specification	Detail
Model Number	Beurer IH 17 ph
Product Dimensions (LxWxH)	160 x 155 x 84 Centimeters
Item Weight	1 kg 300 g
Country of Origin	Germany
Manufacturer	Beurer GmbH Soflinger Str. 218, 89077 Ulm, Germany
Importer (India)	Beurer India Pvt. Ltd., 523-524, JMD Megapolis, Sohna Road, Sector- 48, Gurugram- 122018 India

Specification	Detail
Included Components	Nebulizer compressor unit, air tube, atomizer, adult mask, child mask, mouthpiece, nosepiece, filters.
First Available Date	26 August 2024

WARRANTY AND SUPPORT

For warranty information and customer support, please contact Beurer directly using the details below. It is recommended to retain your purchase receipt for warranty claims.

- **Manufacturer Contact:** Beurer GMBH, Soflinger, Ulm, Germany
- **Email:** info@beurer.com
- **Phone:** +49 (731) 3989-0
- **Importer (India):** Beurer India Pvt. Ltd., 523-524, JMD Megapolis, Sohna Road, Sector- 48, Gurugram-122018 India

Please refer to the official Beurer website or the full product manual for detailed warranty terms and conditions specific to your region.

© 2024 Beurer. All rights reserved.

This manual provides general instructions. Always consult a healthcare professional for medical advice.