

Beurer IH 49

Beurer IH 49 Portable Inhaler User Manual

Model: IH 49

1. INTRODUCTION

The Beurer IH 49 Inhaler is a portable nebulizer designed for the effective treatment of the upper and lower respiratory tracts. It is suitable for both adults and children, providing relief for conditions such as colds, asthma, and other respiratory ailments. Its compact size, quiet operation, and rechargeable battery make it ideal for use at home or while traveling.

This manual provides essential information for the safe and correct use of your Beurer IH 49 Inhaler. Please read it thoroughly before first use and keep it for future reference.

2. IMPORTANT SAFETY INFORMATION

- Always consult your doctor or pharmacist before using the inhaler, especially regarding medication types and dosages.
- Do not use tap water for nebulization. Only use distilled water or saline solution as recommended by your healthcare professional.
- Keep the device out of reach of children when not in use.
- Ensure all components are clean and dry before and after each use to maintain hygiene and prevent contamination.
- Do not immerse the main unit in water or other liquids.
- If you experience any discomfort or adverse reactions during inhalation, stop use immediately and consult a healthcare professional.

3. PACKAGE CONTENTS

Upon unpacking, please verify that all items listed below are present and undamaged:

- Beurer IH 49 Inhaler main unit
- Mouthpiece
- Adult mask
- Child mask
- USB-C charging cable
- Storage bag
- Instruction Manual (this document)



Image 3.1: Beurer IH 49 Inhaler and included accessories.

4. PRODUCT OVERVIEW

Familiarize yourself with the components of your Beurer IH 49 Inhaler:

1. **Main Unit:** Contains the power button, indicator lights, and battery.
2. **Medication Chamber:** Where the liquid medication is placed for nebulization.
3. **Mouthpiece:** Used for direct inhalation.
4. **Masks (Adult/Child):** Provide a comfortable fit for inhalation, especially for children or those who prefer mask use.
5. **USB-C Port:** For charging the device.



Image 4.1: Diagram illustrating the main components of the Beurer IH 49 Inhaler.

5. SETUP AND INITIAL USE

5.1 Charging the Device

Before first use, fully charge the Beurer IH 49 Inhaler. Connect the provided USB-C cable to the device and a suitable USB power adapter (not included). The indicator light will show charging status and turn off when fully charged. A full charge typically provides power for up to 4 treatments.

5.2 Assembling for Use

1. Ensure the medication chamber is clean and dry.
2. Open the medication chamber and carefully pour the prescribed liquid medication into it. Do not overfill.
3. Close the medication chamber securely.
4. Attach either the mouthpiece or the appropriate mask (adult or child) firmly to the top of the medication

chamber.

6. OPERATING THE INHALER

6.1 Starting Inhalation

1. Press the power button to turn on the device. The device will start nebulizing at the strongest setting by default.
2. Place the mouthpiece in your mouth or position the mask comfortably over your nose and mouth, ensuring a good seal.
3. Breathe slowly and deeply through the mouthpiece or mask.
4. The device operates quietly, making it suitable for use without disturbance.

6.2 Adjusting Nebulization Strength

The Beurer IH 49 offers two nebulization strength settings. Press the power button once to switch between the two settings. The strongest setting is activated by default upon power-on.

6.3 Ending Inhalation

When the medication chamber is empty, or your prescribed treatment time is complete, press and hold the power button to turn off the device. Always ensure the device is turned off after use.



Image 6.1: Proper use of the Beurer IH 49 Inhaler with a mouthpiece.



Vrij ademen – overal en altijd

Door de fijne verneveling ontstaan bijzonder kleine deeltjes ($3,89\ \mu\text{m}$) die werkzame stoffen diep in de longen transporteren. Ideaal bij verkoudheid, astma of andere aandoeningen van de luchtwegen.

Image 6.2: The Beurer IH 49 Inhaler is suitable for children, shown here with a child mask.

7. CLEANING AND MAINTENANCE

Regular cleaning is crucial for hygienic and effective operation of your inhaler.

7.1 Daily Cleaning

1. After each use, disassemble the mouthpiece/mask and medication chamber from the main unit.
2. Rinse the mouthpiece/mask and medication chamber thoroughly with warm water.
3. The device features an integrated self-cleaning function. Refer to the detailed instructions in the full manual for activating this feature.
4. Allow all parts to air dry completely before reassembling or storing.
5. Wipe the main unit with a soft, damp cloth. Do not submerge the main unit in water.

7.2 Storage

Store the clean and dry inhaler and its accessories in the provided storage bag in a cool, dry place, away from direct sunlight and extreme temperatures.

8. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Battery is depleted.	Charge the device using the USB-C cable.
No mist or weak mist production.	Medication chamber is empty or clogged.	Refill medication chamber or clean thoroughly. Ensure correct liquid is used.
Device makes unusual noise.	Improper assembly or damaged component.	Check all connections. If problem persists, contact customer support.

If you encounter issues not listed here, or if the problem persists after attempting the suggested solutions, please contact Beurer customer support.

9. TECHNICAL SPECIFICATIONS

Model: IH 49

Dimensions: 11.2 x 6.2 x 4.6 cm

Weight: Approximately 410 g

Battery: 1 Lithium-ion battery (included)

Nebulization Performance: ≥ 0.25 ml/min

Particle Size (MMAD): 3.89 μm (Micrometers)

Operating Principle: Vibrating mesh technology

Country of Origin: China

Manufacturer: Beurer

**German
Engineering**



Klein & licht ■
(ca. 93 g)

Bijzonder stil ■

Oplaadbare accu ■
(met USB-C)

Zelfreinigingsfunctie ■

**Schwingmembran-
technologie (deeltjes-
grootte: 3,89 µm)** ■

**Hoge vernevelings-
prestaties en
korte inhalatietijd** ■

Medisch apparaat ■

Veel toebehoren ■
(met opbergzakje)

Image 9.1: Key specifications and dimensions of the Beurer IH 49 Inhaler.

10. WARRANTY AND CUSTOMER SUPPORT

The Beurer IH 49 Inhaler comes with a **5-year warranty**, ensuring quality and reliability. For warranty claims, technical assistance, or any questions regarding your device, please contact Beurer customer support through their official website or the retailer where you purchased the product.

Please retain your proof of purchase for warranty purposes.

beurer

Adem vrij!
Thuis & onderweg



Image 10.1: Beurer's commitment to quality, featuring a 5-year warranty.