

[Manuals.plus](#) /

› [Furlove](#) /

› Furlove AFO Drop Foot Brace Instruction Manual

Furlove B0DKK5DBCD

Furlove AFO Drop Foot Brace Instruction Manual

Model: B0DKK5DBCD

Brand: Furlove

PRODUCT OVERVIEW

The Furlove AFO Drop Foot Brace is a medical ankle foot orthosis designed to support individuals experiencing foot drop. This brace aids in postural correction and improves stability and gait confidence for users with conditions such as stroke, neuropathy (including diabetic neuropathy), MS, CP, peroneal nerve injuries, or those in post-surgical recovery.

It provides dynamic foot drop lift during the swing phase and strong lateral foot stability, facilitating secure heel-to-toe transitions and reducing tripping risks associated with foot drop and weak ankles.



Image: The Furlove AFO Drop Foot Brace, a black orthosis designed to support the ankle and foot.

SETUP AND FITTING

Proper fitting of the AFO Drop Foot Brace is crucial for optimal support and comfort. Please follow these steps:

1. **Size Selection:** Refer to the Furlove AFO detailed size chart to ensure you select the correct size. Measure your foot length and ankle circumference accurately before ordering.
2. **Custom Trimming:** The thin, flexible footplate can be custom-trimmed to achieve a personalized fit for your foot and shoe. Use standard cast scissors for this adjustment. This allows the brace to be adapted for both men's and women's footwear, for either the left or right foot.
3. **Strap Adjustment:** The brace features three high-quality straps for secure fixation. Adjust these straps to ensure a snug yet comfortable fit around your calf and ankle.
4. **Wearing with Shoes:** If wearing the brace with shoes, the front foot straps are removable. Please remove these straps to ensure a proper fit inside your footwear.

New release

Recommended by orthopedic surgeon

Furlove has been committed to the rehabilitation of the foot, while advancing with the time by time, according to the needs and feedback of the masses to improve step by step, we have been a long-term cooperation with the hospital, also got the orthopedic surgeon highly recommended.

2012



2015



Upgraded



Designed by Authoritative Germany Doctor

Updated to 3 high-quality straps to fix it separately instead of the complicated original 8-character strap. Better fixing effect and simpler wearing steps

Image: Illustration demonstrating the custom design feature, allowing users to trim the footplate for a unique fit.

Removable Foot Straps

Please remove the front foot straps while wear it with shoes



Image: Instructions on removing the front foot straps for use with footwear, showing the brace both with and without the straps.

OPERATING INSTRUCTIONS AND INTENDED USE

The Furlove AFO Drop Foot Brace is designed to assist individuals with various conditions causing foot drop. It is intended for daily use to improve mobility and reduce the risk of falls.

- **Application:** Securely fasten the brace around your ankle and calf using the adjustable straps. Ensure the brace provides a comfortable lift to your foot during walking.
- **Dynamic Support:** The brace provides essential dorsiflexion assistance during the swing phase of gait, preventing the foot from dragging.
- **Lateral Stability:** It offers strong lateral ankle and foot orthosis support, aiding in stable heel-to-toe

transitions and enhancing overall balance.

- **Targeted Conditions:** This brace is suitable for individuals experiencing foot drop due to hemiplegia, peroneal nerve injury, sciatica nerve lesion, talipes valgus, or varus foot.

Your Life Companion Comfortable & Effective



Image: An elderly man demonstrating the use of the AFO brace while walking with support, highlighting its role as a comfortable and effective companion.

Suitable for people



- **Foot drag due to hemiplegia**
- **Drop foot due to peroneal nerve injury or sciatica nerve lesion**
- **Talipes Valgus**
- **Varus Foot**



1



2



3



4

Helped over 100,000+ patients

Image: Visual representation of the types of foot conditions the AFO brace is designed to assist, including foot drag and drop foot.

MAINTENANCE AND CARE

To ensure the longevity and hygiene of your Furlove AFO Drop Foot Brace, please follow these maintenance guidelines:

- **Cleaning Liners:** The brace features removable terry cloth liners. These liners are machine-washable for easy cleaning and to maintain lasting hygiene.
- **Brace Cleaning:** Wipe the polypropylene thermoplastic frame of the brace with a damp cloth and mild

soap. Do not use harsh chemicals or abrasive cleaners.

- **Drying:** Allow all components to air dry completely before reassembling or wearing. Do not expose to direct heat sources.
- **Storage:** Store the brace in a clean, dry place away from direct sunlight when not in use.

TROUBLESHOOTING

While the Furlove AFO Drop Foot Brace is designed for durability and comfort, you may encounter minor issues. Here are some common concerns and suggestions:

- **Discomfort or Rubbing:** Ensure the brace is correctly sized and trimmed. Adjust the straps to be snug but not overly tight. Check for any rough edges on the trimmed footplate and smooth them if necessary. If discomfort persists, consult a medical professional.
- **Brace Slipping:** Re-tighten the straps to ensure a secure fit around your calf and ankle. Ensure the brace is positioned correctly on your leg.
- **Reduced Effectiveness:** Verify that the brace is worn consistently and correctly. If the brace feels less supportive over time, inspect it for any signs of damage or wear.
- **Odor:** Regularly wash the removable terry cloth liners as per the maintenance instructions to prevent odor buildup.

If you experience persistent issues or significant discomfort, discontinue use and consult with a healthcare professional or the manufacturer for further assistance.

PRODUCT SPECIFICATIONS

Feature	Detail
Product Dimensions	2 x 2 x 2 inches
Item Weight	5.61 ounces
Manufacturer	Furlove
ASIN	B0DKK5DBCD
Brand	Furlove
Size	Medium (Left)
Specific Uses	Ankle Sprain, Foot Drop, Fracture, MS, Stroke foot, Hemiplegia Foot Drop, assist Walking
Body Part Use	Ankles, Foot
Age Range	Adult

Feature	Detail
Material	Durable Polypropylene Thermoplastic

WARRANTY AND SUPPORT

For specific warranty information regarding your Furlove AFO Drop Foot Brace, please refer to the product packaging or contact the manufacturer directly. Typically, products purchased through Amazon may be eligible for returns or replacements within 30 days of receipt, subject to Amazon's return policy.

For any product-related inquiries, technical support, or assistance, please contact Furlove customer service. You can often find contact information on the product packaging or by visiting the official Furlove store on Amazon:

[Visit the Furlove Store on Amazon](#)

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is a medical device and should be used as directed by a healthcare professional.