

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Flamaker](#) /
- › [Flamaker Power Lift Recliner Chair User Manual](#)

## Flamaker LR84LP

# Flamaker Power Lift Recliner Chair User Manual

Model: LR84LP

## PRODUCT OVERVIEW

---

The Flamaker Power Lift Recliner Chair is designed to provide ultimate comfort and convenience, especially for the elderly or individuals needing assistance with mobility. It combines power lift and recline functionalities with a soothing massage system, all within a classic and ergonomic design.



*Image: The Flamaker Power Lift Recliner Chair, showcasing its beige corduroy fabric and integrated massage points.*

## **SAFETY INFORMATION**

This chair has a maximum weight capacity of 300 pounds. Please ensure that the user's weight does not exceed this limit to ensure safe and proper use of the chair. Always ensure the power cord is not a tripping hazard and is connected to a suitable power outlet. Keep children and pets away from the moving parts of the chair during operation.

## **ASSEMBLY INSTRUCTIONS**

The Flamaker Power Lift Recliner Chair is designed for easy assembly, requiring only a few steps and no

additional tools. The chair is typically delivered in two separate boxes, which may arrive at different times. Please wait for all packages before beginning assembly.

1. Unpack all components from the boxes and verify that all parts are present according to the included packing list.
2. Carefully align the backrest with the base of the chair and slide it into place until it securely locks.
3. Connect the electrical cables as indicated in the detailed diagram provided in the physical instruction manual. Ensure all connections are firm.
4. Connect the power adapter to the chair and then plug it into a wall outlet.

For detailed visual guidance, refer to the assembly diagrams in your physical manual.

## OPERATING THE CHAIR

---

### Power Lift and Recline Function

The chair features a smooth and quiet electric motor for both lifting and reclining. Use the wired remote control to adjust the chair's position.

- **To Recline:** Press and hold the 'recline' button on the remote. The footrest will extend, and the backrest will smoothly move to your desired reclining angle. Release the button to stop at any position.
- **To Lift/Stand:** Press and hold the 'lift' button on the remote. The chair will gently tilt forward and lift, assisting you to a standing position without strain on your back or knees. Release the button to stop.
- **To Return to Upright Position:** Press the 'upright' button or adjust using the recline/lift buttons until the chair is in the desired sitting position.

# POWER LIFT AND RECLINER



Lift



Recliner



## Multifunction to meet different needs



Image: Visual guide illustrating the various positions of the power lift and recline functions, from standing assistance to full recline for sleeping.



*Image: An elderly woman demonstrating the ease of standing up from the chair using its power lift function.*

## Massage Function

The chair is equipped with a four-point massage system targeting the back, lumbar, thigh, and leg areas. It offers three adjustable modes and three intensity options for a customizable massage experience.

- **Activating Massage:** Use the remote control to select your desired massage mode and intensity.
- **Timer Function:** Set the massage duration for 15, 30, or 60 minutes using the timer function on the remote.
- **Heat Function:** The chair also includes a heat function for added comfort, particularly beneficial for waist pressure relief and blood circulation.



Image: The remote control for the chair, detailing buttons for massage modes, intensity, and timing, alongside a diagram of the massage points.

## FEATURES AND DESIGN

The Flamaker Power Lift Recliner Chair is designed with user comfort and convenience in mind:

- **Material:** Constructed with thick cushions padded with high-density sponges and covered in soft, skin-friendly corduroy fabric.
- **Support:** Features a wide seating space, thick sponge padded back, armrests, and bottom seat. The slightly raised headrest provides neck support.
- **Convenience:** Includes two cup holders and side pockets for storing small items, books, or beverages.
- **Durability:** Built with a sturdy metal frame for long-lasting quality and safety.



*Image: Close-up views highlighting the padded backrest, convenient side pocket, and integrated cup holder, along with a texture sample of the corduroy fabric.*



*Image: A person enjoying the comfort of the recliner chair in a relaxed, reclined position.*

## MAINTENANCE

Proper care will extend the life and appearance of your Flamaker Power Lift Recliner Chair:

- **Cleaning:** For general cleaning, wipe the corduroy fabric with a soft, damp cloth. For spills, blot immediately with a clean, dry cloth. Avoid harsh chemicals or abrasive cleaners.
- **Vacuuming:** Regularly vacuum the chair using a soft brush attachment to remove dust and debris from the fabric.
- **Inspections:** Periodically check all connections and moving parts to ensure they are secure and functioning correctly. Do not attempt to lubricate any internal mechanisms unless specifically instructed by the manufacturer.
- **Power Cord:** Ensure the power cord is not pinched or damaged. If damage occurs, discontinue use and contact customer service.

## TROUBLESHOOTING

If you encounter any issues with your Flamaker Power Lift Recliner Chair, please refer to the following common solutions:

- **Chair Not Responding:**
  - Check if the power cord is securely plugged into both the chair and the wall outlet.
  - Ensure the wall outlet is functional by plugging in another appliance.
  - Verify that all cable connections between the chair components are firm.
- **Massage/Heat Not Working:**
  - Confirm the chair is powered on and the remote control is properly connected.
  - Check the remote for any error indicators or low battery (if applicable).
- **Unusual Noises:**
  - Ensure the chair is on a flat, stable surface.
  - Check for any loose components or obstructions under the chair.

If the problem persists after attempting these steps, please contact customer service for further assistance.

## SPECIFICATIONS

Feature	Detail
Brand	Flamaker
Model Name	Power Lift Chair (LR84LP)
Color	Beige
Product Dimensions (D x W x H)	33.8"D x 38.2"W x 41.7"H
Maximum Weight Recommendation	300 Pounds
Special Features	Power Recliner, Arm Rest, Ergonomic, Power Lift, Cup Holders, Massage, Heat
Number of Cup Holders	2
Fill Material	Foam, Sponge
Item Weight	98 pounds

## PRODUCT DIMENSION

350  
LBS



Image: Detailed product dimensions for the Flamaker Power Lift Recliner Chair.

## WARRANTY AND SUPPORT

Your Flamaker Power Lift Recliner Chair comes with a manufacturer's warranty. Please refer to the warranty card included with your product for specific terms and conditions. Generally, the product is covered for 12 months from the date of purchase against manufacturing defects.

For any questions, concerns, or to report issues, please contact Flamaker customer support. Keep your purchase receipt and model number (LR84LP) ready for faster service.

- **Returns:** This product is eligible for 30-day easy returns.
- **Protection Plans:** Extended protection plans may be available for purchase. Refer to your retailer for details.

