

## ispuoocti Wooden Walking Cane

# ispuoocti Wooden Walking Cane Instruction Manual

## 1. PRODUCT OVERVIEW

---

The ispuoocti Wooden Walking Cane is designed to provide stable support and comfort. Crafted from high-quality ash wood, it features an ergonomic handle and a durable construction suitable for daily use. The cane is enhanced with a non-fading printed aluminum alloy tube for a distinctive appearance.



Figure 1: Overview of the ispuoociti Wooden Walking Cane in use.

## 2. SETUP AND HEIGHT ADJUSTMENT

The ispuoociti Wooden Walking Cane comes at a standard length of 37.4 inches (net length 36.6 inches) and can be customized to your specific height for optimal comfort and support. It weighs approximately 0.99 lbs.

### Determining Correct Cane Height:

1. Remove the non-slip foot pad from the bottom of the cane.
2. Stand upright with your arms hanging naturally at your sides.
3. Position the cane upright on a flat surface next to you.
4. The top of the cane handle should align with the crease of your wrist.

### Adjusting Cane Length:

- Once the correct height is determined, mark the shaft of the cane at the wrist crease level.
- Carefully saw off the excess portion of the cane shaft using a suitable saw. **Important:** When cutting, it is recommended to reserve 1-3 inches of length before the final cut to avoid cutting too short. The

minimum recommended length is 27.56 inches.

- After cutting, firmly re-install the non-slip foot pad onto the bottom of the cane. Tap the cane firmly on the ground until the foot pad is fully seated.



Figure 2: Guide for measuring and adjusting cane height.

### 3. OPERATING INSTRUCTIONS

The ispuoociti Wooden Walking Cane is designed for comfortable and stable support during walking. Its ergonomic handle, shaped like a bird's back, fits the palm's curve to reduce hand fatigue.

#### Proper Usage:

- Hold the cane in the hand opposite to your weaker leg.
- Move the cane forward simultaneously with your weaker leg.
- Ensure the cane's tip is firmly on the ground before applying weight.
- Maintain an upright posture, using the cane for balance and support, not to lean heavily.

#### Safety Considerations:

- The cane is designed to support up to 500 lbs. Do not exceed this weight limit.

- Before each use, check that the handle is securely connected to the shaft. The connection uses 3.15-inch internal and external thread self-tapping screw nuts.
- Ensure the non-slip foot pad is intact and free from wear or damage. Replace worn tips immediately.
- Exercise caution on uneven, wet, or slippery surfaces.



Figure 3: Demonstrating proper walking technique with the cane.

## 4. MAINTENANCE

Proper maintenance ensures the longevity and safety of your ispuoociti Wooden Walking Cane.

- **Cleaning:** Wipe the cane regularly with a soft, damp cloth. Avoid harsh chemicals or abrasive cleaners that could damage the wood or finish.
- **Wood Protection:** The cane is coated with wood wax oil, which protects the ash wood and enhances its natural texture. To maintain its appearance and durability, periodically reapply a suitable wood wax or polish designed for natural wood products.
- **Tip Replacement:** Inspect the rubber foot pad regularly for wear and tear. A worn tip can compromise stability and safety. Replace it promptly with a new, compatible non-slip tip. The product includes two

extra cane tips.

- **Storage:** Store the cane in a dry place away from extreme temperatures and direct sunlight to prevent wood warping or cracking.

## 5. TROUBLESHOOTING

---

If you encounter any issues with your ispuoociti Wooden Walking Cane, refer to the following common solutions:

- **Loose Handle:** If the handle feels loose, gently twist it clockwise to ensure the internal screw nuts are fully tightened. Do not overtighten.
- **Cane Feels Unstable:**
  - Check if the non-slip foot pad is securely attached and not worn out. Replace if necessary.
  - Verify that the cane length is appropriate for your height. An incorrectly sized cane can affect stability.
- **Squeaking Noise:** A squeaking noise might indicate a slightly loose connection. Check the handle and ensure it is tightened.

If problems persist after attempting these solutions, please contact customer support.

## 6. SPECIFICATIONS

---

Detailed specifications for the ispuoociti Wooden Walking Cane:



Figure 4: Component breakdown and key features.

Feature	Detail
Brand	ispuoocti
Model Name	Wooden Walking Cane
Material	Ash Wood
Shaft Material	Ash
Handle Material	Ash Wood
Color	Black brown
Extended Length	37.4 Inches (customizable)
Item Weight	0.99 Pounds
Maximum Height Recommendation	37.4 Inches

Feature	Detail
Maximum Support Capacity	500 lbs
Shaft Diameter	0.83"-0.94"
Included Components	Cane tip (comes with two extra)

## 7. WARRANTY AND SUPPORT

The ispuooci Wooden Walking Cane comes with a **1-year warranty** against non-man-made damage quality problems from the date of receipt.

For warranty claims, technical support, or any inquiries regarding your product, please contact ispuooci customer service through the retailer where the product was purchased. Please have your purchase details and product model information available when contacting support.

© 2025 ispuooci. All rights reserved.