

BEBIRD YD PINK CS

Bebird Earsight Plus R37R Ear Cleaning Camera Wax Removal Tool Kit User Manual

Model: YD PINK CS

1. PRODUCT OVERVIEW

The Bebird Earsight Plus R37R is an advanced ear cleaning camera and wax removal tool kit designed for visual inspection and gentle cleaning of the ear canal. It features a high-definition camera, flexible otoscope, and precision tweezers, providing a comprehensive solution for ear care for both adults and children. This manual provides detailed instructions for the safe and effective use of your device.

2. PRODUCT FEATURES

- **Flexible Otoscope Camera:** Equipped with a bendable otoscope designed to navigate the natural curves of the ear canal. Features a 3P waterproof, dustproof, and anti-fog HD camera for clear internal views.
- **Enhanced Safety and Comfort:** Features a soft silicone coating on tips to prevent scratches and discomfort. The spiral tip design ensures secure attachment during use. Includes a built-in safety scale and protective stopper to prevent over-insertion and protect the eardrum.
- **Precision Visual Ear Tweezers:** Specialized tweezers with an integrated camera allow for gentle and precise removal of ear debris without direct skin contact, suitable for sensitive ear canals.
- **Comprehensive Ear Care Kit:** Comes with 4 types of cleaning tips suitable for both wet and dry earwax. The kit also includes components for ear, nose, and throat checkups.
- **Portable Design:** Supplied with a compact storage case, making it convenient for home use and travel.
- **Durable Construction:** Manufactured with high-quality materials for reliable and long-term performance.

3. PACKAGE CONTENTS

Please verify that all items listed below are included in your package. If any items are missing or damaged, please contact customer support.



Image: The Bebird Earsight Plus R37R kit, showing the main camera unit, ear tweezers, charging case, various ear tips, charging cable, and cleaning brush. Two smartphone screens display the camera's view inside an ear canal.

- Bebird Earsight Plus R37R Main Camera Unit (Flexible Otoscope)
- Visual Ear Tweezers with Camera
- Charging/Storage Case
- Assorted Ear Cleaning Tips (e.g., silicone ear spoons, brush tips, spiral tips, cotton swab tips)
- USB-C Charging Cable
- Cleaning Brush
- User Manual (this document)

4. SAFETY INFORMATION

Read all safety warnings and instructions carefully before using this product. Failure to follow these instructions may result in injury or damage to the device.

- Do not insert the device too deeply into the ear canal. Always observe the safety scale and protective stopper.
- Ensure the ear tips are securely attached before use.

- Clean the ear tips thoroughly before and after each use to prevent infection.
- Do not use the device if you experience pain, discomfort, or bleeding. Consult a medical professional.
- Keep out of reach of children unless under direct adult supervision.
- Avoid extreme temperatures and direct sunlight.
- Do not attempt to disassemble or repair the device yourself.

5. SETUP GUIDE

5.1 Charging the Device

1. Connect the USB-C charging cable to the charging port on the main camera unit or the charging case.
2. Plug the other end of the cable into a standard USB power adapter (not included) or a computer USB port.
3. The indicator light will show charging status (refer to device specific light indicators for full/charging status).
4. Allow the device to fully charge before first use, typically 1-2 hours.

5.2 App Installation and Connection

1. Download the official Bebird app from your smartphone's app store (e.g., Apple App Store or Google Play Store). Search for "Bebird".
2. Turn on the Bebird Earsight Plus R37R main camera unit by pressing and holding the power button.
3. Open the Bebird app on your smartphone.
4. Follow the in-app instructions to connect the device via Wi-Fi. The device will broadcast its own Wi-Fi network (e.g., "Bebird-XXXX"). Connect your phone to this network.
5. Once connected, the live view from the camera will appear on your phone screen.

6. OPERATING INSTRUCTIONS

6.1 Using the Flexible Otoscope Camera

1. Select the appropriate ear tip for your needs and securely attach it to the flexible otoscope camera.
2. Ensure the device is connected to the app and the live view is active.
3. Gently insert the otoscope into the ear canal, observing the live feed on your smartphone.
4. Use the safety scale on the otoscope as a guide to prevent over-insertion.
5. Carefully maneuver the otoscope to inspect the ear canal and identify earwax or debris.
6. Use the attached ear spoon or other cleaning tips to gently scoop or dislodge earwax. Always move slowly and deliberately.
7. Avoid touching the eardrum. If you feel any resistance or pain, stop immediately.
8. Once cleaning is complete, gently remove the otoscope.

6.2 Using the Visual Ear Tweezers

1. Ensure the visual ear tweezers are charged and connected to the app.
2. Gently insert the tweezer tip into the ear canal, guided by the live camera view on your smartphone.
3. Carefully position the tweezers to grasp ear debris.
4. Gently close the tweezers to secure the debris and slowly pull it out.

- 5. Ensure the tweezers do not touch the ear canal skin during removal to prevent discomfort.
- 6. Clean the tweezers thoroughly after each use.

7. MAINTENANCE AND CARE

7.1 Cleaning the Device

- After each use, detach the ear tips from the otoscope and tweezers.
- Wash the ear tips with mild soap and warm water. Rinse thoroughly and allow them to air dry completely before storing.
- Use the provided cleaning brush to remove any wax or debris from the camera lens and tip area.
- Wipe the main body of the device with a soft, damp cloth. Do not immerse the main unit in water.
- Ensure all components are dry before placing them back into the charging/storage case.

7.2 Storage

Store the Bebird Earsight Plus R37R and its accessories in the provided compact storage case when not in use. Store in a cool, dry place away from direct sunlight and extreme temperatures.

8. TROUBLESHOOTING

If you encounter any issues with your Bebird Earsight Plus R37R, please refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Device does not turn on.	Low battery.	Charge the device fully.
Cannot connect to the app.	Wi-Fi not enabled or incorrect network selected.	Ensure device Wi-Fi is on and connect your phone to the correct "Bebird-XXXX" network. Restart both device and app.
Image is blurry or unclear.	Camera lens is dirty or fogged.	Clean the camera lens with the provided brush or a soft, dry cloth. Allow device to acclimate to room temperature if fogging occurs.
Ear tip falls off during use.	Ear tip not securely attached.	Ensure the ear tip is firmly clicked into place before use.

9. SPECIFICATIONS

Feature	Detail
Model Number	YD PINK CS
Camera Resolution	HD (High Definition)
Waterproof Rating	3P (Waterproof, Dustproof, Anti-fog)

Feature	Detail
Connectivity	Wi-Fi (App-controlled)
Power Source	Rechargeable Battery (via USB-C)
Package Dimensions	0.2 x 0.2 x 0.2 inches
Item Weight	4.8 ounces
Manufacturer	Heifeng Zhizao (Shenzhen) Technology Co., Ltd

10. WARRANTY AND SUPPORT

Specific warranty information for the Bebird Earsight Plus R37R is not provided in this manual. For details regarding warranty coverage, technical support, or replacement parts, please refer to the product packaging, the official Bebird website, or contact the retailer from whom you purchased the product. Please retain your proof of purchase for warranty claims.