

## PEPE Mobility P10037

# PEPE Mobility Wheelchair Tray with Cup Holder Instruction Manual

Model: P10037

## 1. PRODUCT OVERVIEW

The PEPE Mobility Wheelchair Tray with Cup Holder is designed to enhance comfort and functionality for wheelchair users. It provides a stable surface for various activities such as eating, reading, or using electronic devices. Its lightweight yet durable construction ensures ease of use and portability.

### Key Features:

- **Functional Design:** Offers a stable and practical surface with an integrated cup holder for beverages.
- **Easy Installation:** Attaches quickly and securely with adjustable Velcro straps.
- **Lightweight and Durable:** Constructed from high-quality ABS plastic, weighing approximately 2.53 lb.
- **Versatile Use:** Suitable for various environments including home, hospitals, and outdoors.

## 2. SAFETY INFORMATION

- Ensure the tray is securely attached to the wheelchair armrests before use.
- Do not overload the tray with excessive weight.
- Keep hot liquids and sharp objects away from the edges to prevent spills or injury.
- Regularly inspect the Velcro straps for wear and tear and replace if necessary.
- Always ensure the user has adequate space and is comfortable with the tray in place.

## 3. PACKAGE CONTENTS

- 1 x PEPE Mobility Wheelchair Tray
- 2 x Adjustable Velcro Straps (pre-attached)

## 4. SETUP AND INSTALLATION

1. Position the wheelchair tray over the armrests of the wheelchair, ensuring the cut-out section aligns with the user's body.
2. Locate the pre-attached Velcro straps on the underside of the tray.

3. Wrap each Velcro strap securely around the corresponding armrest of the wheelchair.
4. Thread the strap through the buckle and pull tightly to ensure a snug fit. Press the Velcro surfaces together to secure.
5. Verify that the tray is stable and does not wobble before use. Adjust the straps as needed for optimal security.



Image: Attaching the Velcro strap to the wheelchair armrest.



Image: Underside of the tray with Velcro straps.

### Installation Video:

Video: Demonstrates the installation process of the PEPE Mobility Wheelchair Tray with Cup Holder.

## 5. OPERATING INSTRUCTIONS

---

### Using the Tray:

- The tray provides a flat, stable surface for meals, reading, writing, or using laptops and tablets.
- Utilize the integrated cup holder to securely place beverages, reducing the risk of spills.
- The raised edges of the tray help prevent items from sliding off.



Image: Demonstrating versatile uses of the wheelchair tray.



Image: User eating comfortably with the tray.





Image: An elderly user enjoying the tray for daily activities.

## 6. MAINTENANCE AND CARE

---

- Clean the tray regularly with a damp cloth and mild soap.
- Avoid using abrasive cleaners or harsh chemicals, as these may damage the surface.
- Wipe dry immediately after cleaning to prevent water spots.
- Store the tray in a clean, dry place when not in use.



Image: Cleaning the tray surface.

## 7. TROUBLESHOOTING

---

- **Tray feels unstable:** Ensure both Velcro straps are pulled tight and securely fastened around the armrests. Re-adjust if necessary.
- **Tray does not fit:** Verify that your wheelchair armrests are compatible with the tray's design and strap length. The tray is designed for most standard wheelchairs.
- **Difficulty attaching straps:** Ensure the straps are threaded correctly through the slots on the tray and around the armrests.

## 8. PRODUCT SPECIFICATIONS

---

Specification	Detail
Model Number	P10037
Material	High-quality ABS plastic
Dimensions (L x W x H)	22.8 x 20.7 x 1.4 inches (58 x 52.5 x 3.5 cm)

Specification	Detail
Weight	2.53 lb (1.15 kg)
Color	Black



Image: Wheelchair tray dimensions (top view).



Image: Wheelchair tray thickness (side view).

9. WARRANTY AND CUSTOMER SUPPORT

For any questions, concerns, or warranty information regarding your PEPE Mobility Wheelchair Tray, please contact PEPE Mobility customer support through the retailer where the product was purchased. Please retain your



proof of purchase for warranty claims.