

[Manuals.plus](#) /

› [Rhevoy](#) /

› Rhevoy PL003 Oversized Power Lift Recliner Chair with Massage and Heat for Elderly

## Rhevoy PL003

# Rhevoy PL003 Oversized Power Lift Recliner Chair with Massage and Heat for Elderly

Model: PL003

## INTRODUCTION

This manual provides comprehensive instructions for the assembly, operation, maintenance, and troubleshooting of your Rhevoy PL003 Oversized Power Lift Recliner Chair. Designed for comfort and ease of use, this chair features power lift assistance, reclining positions, massage, and heat functions, along with convenient USB charging ports.

## IMPORTANT SAFETY INFORMATION

### WARNING:

- Do not exceed the weight capacity of 350 lbs. Exceeding this limit may cause the chair to malfunction or break, leading to injury.
- Keep children and pets away from the chair's moving parts.
- Ensure all electrical connections are secure and properly plugged in before operation.
- Avoid placing the chair near heat sources or in direct sunlight to prevent material damage.
- Do not stand on the footrest or armrests.

## PACKAGE CONTENTS

Your Rhevoy PL003 Power Lift Recliner Chair is shipped in two boxes. Please ensure all components are present before beginning assembly.

- Chair Base (1)
- Chair Backrest (1)
- Armrests (2)
- Power Adapter (1)
- Remote Control for Massage/Heat (1)

- Assembly Hardware (bolts, nuts, etc.)
- Instruction Manual (this document)



*Image: The Rhevoy PL003 Power Lift Recliner Chair in beige, showcasing its oversized design.*

## SETUP AND ASSEMBLY

Assembly is designed to be straightforward. It is recommended to have two people for easier assembly.

1. **Unpack Components:** Carefully remove all parts from the two shipping boxes. Ensure no parts are missing or damaged.
2. **Attach Armrests:** Locate the metal slots on the chair base. Slide each armrest onto the corresponding metal slot until it is securely in place. Ensure the armrest cushion covers are zipped and velcroed properly.
3. **Install Backrest:** Align the backrest with the receiving brackets on the chair base. Slide the backrest down firmly until it locks into place.
4. **Connect Wiring:** Underneath the chair, connect all electrical wires. This typically involves connecting the power adapter to the chair's main power cable, and the massage/heat remote control cable. Ensure all connections are snug.

5. **Attach Rear Rollers (Optional/If Included):** If your model includes rear rollers for mobility, attach them to the designated points on the chair's base using the provided hardware.
6. **Power On:** Plug the power adapter into a wall outlet.

Your browser does not support the video tag.

*Video: Product Installation Guide. This video demonstrates the step-by-step assembly process for the Rhevo Power Lift Recliner Chair, including attaching armrests, backrest, and connecting electrical components.*

## Lift Chairs Provide Any Position for You



*Image: Detailed product dimensions of the Rhevo PL003 Power Lift Recliner Chair, including seat width, depth, and overall height.*

## OPERATING INSTRUCTIONS

### Power Lift and Recline Functions

The chair's power lift and recline functions are controlled by the attached remote. This allows for smooth transitions between sitting, reclining, and standing positions.

- **To Recline:** Press and hold the "recline" button on the remote control. The chair will smoothly move into a reclined position. Release the button when the desired angle is reached (up to 150 degrees).
- **To Sit Up:** Press and hold the "sit up" button on the remote control. The chair will return to an upright

sitting position.

- **To Lift/Assist Standing:** Continue holding the "sit up" button past the upright position. The chair will gently tilt forward and lift, assisting the user to a standing position.

Your browser does not support the video tag.

*Video: Demonstration of the 26" Extra Wide Seat Power Lift Recliner Chair. This video highlights the chair's power lift and reclining capabilities, showing various comfortable positions.*

## Massage and Heat Functions

The chair features 8 built-in massage points and lumbar heating, controlled by a separate remote.

- **Power On/Off:** Press the "POWER" button to activate or deactivate the massage and heat functions.
- **Heat Function:** Press the "HEAT" button to turn on the lumbar heating. Press again to turn off.
- **Massage Modes:** Use the "MODE" button to cycle through 5 different massage modes (e.g., Pulse, Press, Wave, Auto, Normal).
- **Intensity Adjustment:** Adjust massage intensity using the "HIGH" and "LOW" buttons.
- **Targeted Massage:** Select specific massage zones (Back, Lumbar, Thigh, Leg) using the corresponding buttons.
- **Timer:** Set the massage duration to 15, 30, or 60 minutes using the "TIMER" button.



Image: Close-up of the remote control for the massage and heat functions, illustrating the various buttons for mode, intensity, and zone selection.

## Additional Features

- **USB Charging Port:** A built-in USB port is located on the side of the chair for convenient charging of your electronic devices.
- **Side Pockets:** Handy side pockets are available for storing remote controls, books, or other small items.
- **Removable Armrest Cover:** The armrest cover is easily removable for cleaning.

## MAINTENANCE

Proper care will extend the life and appearance of your recliner.

- **Cleaning:** Spot clean the chenille fabric with a mild detergent and a damp cloth. Avoid harsh chemicals or abrasive cleaners.
- **Spills:** Blot spills immediately with a clean, dry cloth.
- **General Care:** Regularly vacuum the chair to remove dust and debris.
- **Electrical Components:** Ensure all electrical cords are free from kinks or damage. Do not attempt to repair electrical components yourself.

## TROUBLESHOOTING

If you encounter any issues with your recliner, please refer to the following common solutions:

Problem	Possible Cause	Solution
Chair does not operate (no lift/recline)	Power cord disconnected; outlet not working; remote control issue.	Check power connections. Test outlet with another device. Ensure remote control batteries (if applicable) are functional.
Massage/Heat not working	Remote control not connected; function not activated.	Ensure massage/heat remote is properly connected. Press the "POWER" and "HEAT" buttons to activate.
Chair makes unusual noises	Loose hardware; obstruction in mechanism.	Inspect for loose bolts and tighten if necessary. Check for any objects obstructing the reclining mechanism.

If the problem persists, please contact customer support.

## SPECIFICATIONS

Feature	Detail
Brand	Rhevoy
Model Name	PL003
Color	Beige

Feature	Detail
Product Dimensions	29"D x 40"W x 44"H
Seat Depth	21 inches
Maximum Weight Recommendation	350 Pounds
Material	Chenille Fabric, Foam, Fiber Batting, Metal Frame
Special Features	Power Lift, Massage, Heat, USB Port, Removable Armrest Cover, Oversized
Reclining Position Count	3 (Upright, Reclined, Lifted)
Product Care Instructions	Spot Clean

## WARRANTY AND SUPPORT

---

Your Rhevoy PL003 Power Lift Recliner Chair comes with a **2-year warranty** covering faulty, damaged, or missing parts.

For any questions regarding your purchase, assembly, or operation, please contact Rhevoy customer service. We are committed to providing support and ensuring your satisfaction.

Please refer to the contact information provided with your purchase documentation or visit the official Rhevoy website for assistance.