



Manuals.plus /

› FACULX /

› FACULX 360° Swivel Car Seat Cushion User Manual

FACULX FACULX

FACULX 360° Swivel Car Seat Cushion User Manual

Model: FACULX

1. INTRODUCTION

The FACULX 360° Swivel Car Seat Cushion is designed to provide effortless rotation, assisting individuals with mobility challenges, such as pregnant women, the elderly, and those with hip or back pain, in entering and exiting vehicles or other seating arrangements. This portable cushion aims to reduce strain and enhance comfort during use.

Key features include:

- 360-degree rotating design for easy movement.
- Non-slip base for stability.
- Comfortable and breathable materials (pleated cloth, soft sponge, durable nylon, memory foam).
- Portable design with a built-in handle.
- Suitable for various seating environments beyond cars.



Image: The FACULX 360° Swivel Car Seat Cushion in gray, shown from a top-down perspective, highlighting its circular design and integrated handle.

2. SETUP AND INSTALLATION

Proper installation ensures optimal performance and safety of your swivel seat cushion. This product is suitable for most car models, provided the seating surface is flat and stable.

Installation Steps:

1. **Prepare the Seat:** Ensure the car seat surface where the cushion will be placed is clean, flat, and stable. Soft or uneven surfaces may hinder the rotation.
2. **Position the Cushion:** Place the FACULX swivel seat cushion firmly on the desired seat.
3. **Secure the Cushion:** If applicable, stuff any excess material or sponge from the cushion into the gap of the car seat to ensure it is packed tightly and remains close to the seat. This prevents tilting or hanging that could affect rotation.
4. **Verify Stability:** Before use, gently test the cushion to ensure it is stable and does not slide or move unexpectedly.

How to install and use 360° swivel seat cushion

1. Stuff the sponge into the gap of the car seat
2. Make sure the seat cushion is packed tightly and can be used normally



Image: A visual guide illustrating the four steps for installing the 360° swivel seat cushion into a car seat, emphasizing proper placement and securing.

3. OPERATING INSTRUCTIONS

The FACULX 360° Swivel Car Seat Cushion is designed for ease of use in various scenarios.

Using the Swivel Function:

- Sit centrally on the cushion.
- To rotate, gently shift your weight and turn your body in the desired direction. The cushion will rotate smoothly with you.
- The 360-degree rotation allows for easy entry and exit from vehicles, chairs, or other seating, reducing the need for twisting your body.



Image: The swivel seat cushion positioned on a car seat, demonstrating its 360-degree rotation capability with directional arrows.

Recommended Use:

- The cushion is recommended for users weighing less than 220 lbs (90 kg) to ensure optimal rotation flexibility and product longevity. Exceeding this weight may increase rotation resistance.
- This versatile cushion can be used in cars, office chairs, home chairs, for travel, or even as a floor mat for yoga or meditation.



Image: A collage showcasing the FACULX swivel seat cushion's adaptability for use in diverse environments, including office chairs, car seats, during travel, and with wheelchairs.

4. MAINTENANCE

Proper care and maintenance will extend the life and performance of your FACULX swivel seat cushion.

Cleaning Instructions:

- The cushion is designed with a double-layer structure, allowing for easy removal of the top layer for cleaning.
- **Hand Wash Only:** The fabric cover should be hand washed with mild soap and cold water. Do not machine wash, bleach, or tumble dry.
- Allow the cover to air dry completely before reassembling or using.



Image: A close-up view demonstrating the cushion's double-layer design, where the top layer can be easily removed for cleaning.

General Maintenance:

- Regularly clean the rotating base to prevent dust or debris accumulation, which can affect smoothness.
- If the rotation smoothness decreases after long-term use, slightly adjust the cushion's position on the seat to restore it.

5. TROUBLESHOOTING

If you encounter issues with your FACULX swivel seat cushion, refer to the following common problems and solutions:

Issue: Cushion does not rotate smoothly or is difficult to turn.

- **Cause 1: Uneven or Soft Seating Surface.**
Solution: Ensure the cushion is placed on a flat and stable surface. Car seats that are too soft or uneven can impede rotation. Adjust the cushion's position or consider adding a firm base underneath if possible.

- **Cause 2: Exceeding Weight Limit.**

Solution: The cushion is designed for users weighing less than 220 lbs (90 kg). Exceeding this weight may cause increased resistance.

- **Cause 3: Debris in Rotating Mechanism.**

Solution: Regularly inspect and clean the rotating base to remove any dust or debris that might be obstructing smooth movement.

- **Cause 4: Cushion Not Tightly Packed.**

Solution: Ensure the cushion is fully settled and any excess material is tucked into the seat gap to prevent tilting, which can affect rotation.

6. SPECIFICATIONS

Detailed specifications for the FACULX 360° Swivel Car Seat Cushion:

- **Brand:** FACULX
- **Model Name:** 360° Rotating Function
- **Color:** Gray
- **Dimensions:** 17.32"L x 11.42"W x 0.8"H (44cm x 29cm x 2cm)
- **Item Weight:** 9.9 ounces (0.28 Kilograms)
- **Fill Material:** Pleated cloth, Soft sponge, Memory Foam
- **Fabric Type:** Combination of high-quality pleated cloth, soft sponge, and durable nylon
- **Special Feature:** 360° Rotating Function, Non-slip bottom
- **Product Care Instructions:** Hand Wash Only
- **Recommended Uses:** Driving, Traveling, Rocking chair, Office chairs, Home chairs, Yoga, Meditation
- **Age Range (Description):** Adult
- **Number of Pieces:** 1

Product Size



Image: A visual representation of the product's dimensions, indicating its length, width, and thickness for clarity.

7. WARRANTY AND SUPPORT

For information regarding warranty coverage, returns, or technical support, please refer to the retailer's policy or contact FACULX customer service directly. Specific warranty details are not provided within this manual.

For further assistance, please visit the official FACULX store or contact their support channels.