

## HOMCOM 713-194V80CF

# HOMCOM Power Lift Recliner Chair User Manual

Model: 713-194V80CF

### INTRODUCTION

Thank you for choosing the HOMCOM Power Lift Recliner Chair. This manual provides essential information for the safe assembly, operation, and maintenance of your new recliner. Please read it thoroughly before use and retain it for future reference.

### SAFETY INFORMATION

Ensure proper assembly and stability before use. Do not exceed the maximum weight capacity of 330 pounds. Follow instructions for safe operation of the power lift mechanism.

### PRODUCT OVERVIEW

The HOMCOM Power Lift Recliner Chair is designed for comfort and ease of use, particularly for the elderly. It features a power lift mechanism, stepless reclining up to 150°, and a remote control for effortless adjustments.

#### Key Features:

- Effortless Lift and Stepless Recline: Smoothly elevates to a near-standing position and reclines up to 150° with remote control.
- Comfort: Features pocket spring support and plush foam padding, upholstered in a soft, breathable, linen-feel fabric.
- Elegant Design: Tufted back and flared arms complement various home decors.
- Convenient Storage: Includes a side pocket for storing remotes, magazines, or books.

#### Product Dimensions:



Image: Detailed diagram illustrating the various dimensions of the recliner chair in both upright and reclined positions.

- Overall Dimensions: 34.6"D x 26.8"W x 40.9"H
- Seat Height: 19.3 Inches
- Seat Depth: 20.1 inches
- Seat Length: 18.9 Inches
- Backrest Width: 19.7 Inches
- Seat Back Interior Height: 21.7 Inches
- Arm Height: 3.9 Inches
- Maximum Weight Recommendation: 330 Pounds

#### Included Components:

- Arm Pad
- Remote Controller
- Side Pocket
- Tufted Back

#### SETUP AND ASSEMBLY

The HOMCOM Power Lift Recliner Chair is designed for quick assembly, requiring no tools and taking

approximately 5 minutes. Follow the steps below for proper setup.

### Assembly Steps:

1. Unpack all components from the packaging.
2. Attach the backrest to the seat base by aligning the metal brackets and pushing down until it clicks into place.
3. Connect the power cables for the lift and recline mechanism.
4. Ensure all connections are secure and the chair is stable before use.



Image: The recliner chair in a living room setting, highlighting its quick and easy assembly process.

### Assembly Video:

Video: Installation guide for a similar power lift recliner chair, demonstrating the assembly process.

Video: Detailed installation instructions for a power lift recliner chair, showing step-by-step assembly.

## OPERATING INSTRUCTIONS

Your HOMCOM Power Lift Recliner Chair is operated using a wired remote control, allowing for easy adjustment of the lift and recline functions.

### Remote Control Functions:



Image: The remote control and its functions for adjusting the chair's position.

- **Lift Function:** Press and hold the 'Up' button to gently raise the chair to a near-standing position, assisting users in getting up easily.
- **Recline Function:** Press and hold the 'Down' button to smoothly recline the chair to any angle up to 150°, extending the footrest for full relaxation.
- Release the button at any desired position to stop the movement.

### Operation Videos:

Video: Demonstration of the HOMCOM 713-194V80CF Lift Recliner Chair in operation, showing its lift and recline capabilities.

Video: Overview of the HOMCOM Power Lift Assist Recliner Chair for Elderly, demonstrating its remote-controlled functions.

Video: Demonstration of the HOMCOM 713-042V80BK Power Lift Chair, showcasing its movement and features.

Video: Short preview of the HOMCOM 713-042V81BG Power Lift Chair, demonstrating its basic lift function.

Video: Demonstration of a power lift recliner chair, showing its recline and lift functions.

### MAINTENANCE



To ensure the longevity and optimal performance of your recliner, follow these simple maintenance guidelines:

- **Cleaning:** Wipe clean with a damp cloth. Avoid harsh chemicals that may damage the fabric.
- **Inspection:** Periodically check all connections and moving parts for any signs of wear or damage.
- **Power Cord:** Ensure the power cord is not pinched or damaged.

### TROUBLESHOOTING

If you encounter any issues with your HOMCOM Power Lift Recliner Chair, please refer to the following common troubleshooting tips:

- **Chair not moving:** Check if the power cord is securely plugged into both the chair and a working electrical outlet. Ensure all internal connections are firm.
- **Remote not responding:** Verify that the remote control cable is properly connected to the chair.
- **Unusual noises:** If you hear any unusual noises, stop using the chair immediately and inspect the mechanism for any obstructions or loose parts.

If the problem persists, please contact customer support for assistance.

### TECHNICAL SPECIFICATIONS

Feature	Specification
Brand	HOMCOM
Model Name	Lift Recliner Chair
Item Model Number	713-194V80CF
Color	Coffee
Material	100% Polyester Fabric, Foam Fill
Product Dimensions	34.6"D x 26.8"W x 40.9"H
Item Weight	84 Pounds
Maximum Weight Recommendation	330 Pounds
Reclining Position Count	Infinite Position (up to 150°)
Special Features	Adjustable Backrest, Arm Rest, Back Support, Reclining, Side Pocket
Product Care Instructions	Wipe Clean
UPC	196393230176

### WARRANTY AND SUPPORT

Your HOMCOM Power Lift Recliner Chair comes with a 30-day easy return policy. For any product inquiries, technical support, or warranty claims, please contact HOMCOM customer support.

- **Returns:** 30-day easy returns.
- **Customer Support:** Refer to your purchase documentation for contact details or visit the HOMCOM store on Amazon.

