

## TRIO PICCOLO

# Mega Motion MM-4001-PE Trio Piccolo Lift/Recline Chair User Manual

Model: TRIO PICCOLO

## 1. INTRODUCTION

---

Thank you for choosing the Mega Motion MM-4001-PE Trio Piccolo Single Motor 3 Position Lift/Recline Chair. This manual provides essential information for the safe and efficient operation, setup, and maintenance of your new chair. Please read this manual thoroughly before using the product and retain it for future reference.

## 2. SAFETY INFORMATION

---

Always observe the following safety precautions to prevent injury or damage to the chair:

- **Weight Capacity:** Do not exceed the maximum weight capacity of 375 lbs (170 kg).
- **Clearance:** Ensure adequate clearance around the chair during operation to avoid obstructions.
- **Children and Pets:** Keep children and pets away from moving parts of the chair.
- **Power Cord:** Do not place the power cord in high-traffic areas where it could be tripped over or damaged.
- **Liquid Spills:** Avoid spilling liquids on the chair, especially near electrical components.
- **Emergency Stop:** In case of emergency, disconnect the power cord from the wall outlet.
- **Professional Repair:** Do not attempt to repair the chair yourself. Contact qualified service personnel for any repairs.

## 3. COMPONENTS OVERVIEW

---

Your Mega Motion Trio Piccolo chair consists of the following main components:

- Chair Base with Motor and Lift Mechanism
- Backrest
- Handwand (Remote Control) with USB Port
- Power Cord
- Emergency Battery Backup (Integrated)



**Figure 1:** The Mega Motion Trio Piccolo Lift/Recline Chair shown in its upright, seated position. This view highlights the overall design and the handwand on the side.

## 4. SETUP AND ASSEMBLY

---

The Mega Motion MM-4001-PE Trio Piccolo chair requires minimal assembly. It is designed for quick and easy setup.

### 4.1 Unpacking

1. Carefully remove all components from the packaging.
2. Inspect all parts for any signs of damage. If any damage is found, do not proceed with assembly and contact customer

support.

3. Retain packaging materials for potential future transport or returns.

## 4.2 Connecting Power

1. Locate the power cord.
2. Plug the power cord into the designated port on the chair's motor mechanism (usually located at the rear or underside of the chair).
3. Plug the other end of the power cord into a grounded electrical outlet.
4. Ensure all connections are secure.

## 5. OPERATING INSTRUCTIONS

---

The chair is operated using the handwand (remote control).

### 5.1 Using the Handwand

The handwand features simple buttons for controlling the chair's movement:

- **Recline Button:** Press and hold to recline the chair backward. Release the button when you reach your desired recline position.
- **Lift Button:** Press and hold to raise the chair to a standing assist position. Release the button when you reach your desired height or the fully lifted position.
- **USB Port:** Located on the handwand, this port allows you to charge your mobile devices (e.g., tablet, phone) while seated.



**Figure 2:** The Mega Motion Trio Piccolo Lift/Recline Chair shown in a reclined position, demonstrating the extended footrest and backrest.

### 5.2 Side Pockets

The chair is equipped with easy-to-reach side pockets on both sides. These are convenient for storing magazines,

newspapers, remote controls, or other small personal items.

## 5.3 Emergency Battery Backup

In the event of a power outage, the integrated emergency battery backup allows for a limited number of cycles to return the chair to a seated position. This feature is for emergency use only and is not intended for regular operation. Ensure the batteries are properly installed and maintained according to manufacturer guidelines (typically located near the motor or battery compartment).

## 6. MAINTENANCE

---

Proper maintenance ensures the longevity and performance of your chair.

### 6.1 Cleaning the Fabric (Popstitch, 100% Polyester-W)

- **Regular Cleaning:** Vacuum regularly to remove dust and loose particles.
- **Spot Cleaning:** For spills, blot immediately with a clean, dry cloth. For stubborn stains, use a mild soap solution and a damp cloth. Test on an inconspicuous area first.
- **Avoid Harsh Chemicals:** Do not use harsh detergents, solvents, or abrasive cleaners, as these can damage the fabric.
- **Drying:** Allow the fabric to air dry completely. Do not use heat.

### 6.2 Mechanical Parts

- Periodically check all bolts and screws to ensure they are tight.
- Do not lubricate any moving parts unless specifically instructed by the manufacturer.

## 7. TROUBLESHOOTING

---

If you encounter issues with your chair, refer to the following troubleshooting guide:

Problem	Possible Cause	Solution
Chair does not move	No power to the chair; Power cord disconnected; Handwand malfunction	Check if the power cord is securely plugged into the wall outlet and the chair. Ensure the outlet is functional. Check for any visible damage to the handwand or cord.
Chair moves slowly or intermittently	Low power; Obstruction in mechanism; Overheating motor	Ensure the chair is plugged directly into a wall outlet (avoid extension cords). Check for any objects obstructing the lift mechanism. Allow the chair to cool down if it has been used continuously.
USB port not charging	Device cable issue; Device not compatible; Chair not powered	Try a different USB cable or device. Ensure the chair is powered on.

If the problem persists after attempting these solutions, please contact customer support.

## 8. SPECIFICATIONS

---

Detailed specifications for the Mega Motion MM-4001-PE Trio Piccolo Lift/Recline Chair:

- **Model Name:** Single Motor 3 Position Lift/Recline Chair

- **Model Number:** TRIO PICCOLO
- **Brand:** Generic (Mega Motion)
- **Color:** Pebble
- **Material:** Popstitch (100% Polyester-W)
- **Product Dimensions (D x W):** 40"D x 32.5"W
- **Weight Capacity:** 375 Lbs
- **Suggested User Height:** 5' to 5'6"
- **Seat-To-Floor Height:** 19"
- **Seat Depth:** 19"
- **Seat Width:** 19"
- **Total Width:** 32.5"
- **Total Depth:** 40"
- **Top Of Back To Seat:** 26"
- **Top Of Back To Floor:** 41.5"
- **Total Recline Depth:** 66"
- **Arm Height From Seat:** 5"
- **Arm Height To Floor:** 24"
- **Extended Ottoman Height:** 21"
- **Extended Back To Ottoman:** 38.5"
- **Overall Lifted Chair Height:** 56.5"
- **Distance From Wall:** 16"
- **Assembly Required:** No

## 9. WARRANTY INFORMATION

---

For specific warranty details, please refer to the warranty documentation provided with your purchase or contact the manufacturer directly. Warranty terms typically cover manufacturing defects for a specified period from the date of purchase.

## 10. CUSTOMER SUPPORT

---

If you have any questions, require assistance, or need to report an issue with your Mega Motion MM-4001-PE Trio Piccolo Lift/Recline Chair, please contact the manufacturer's customer support. Refer to your purchase documentation for contact information.