



Manuals.plus /

- › Comfyer /
- › Comfyer Heated Wearable Breast Pump Model W2 Instruction Manual

## Comfyer W2

# Comfyer Heated Wearable Breast Pump Model W2 User Manual

Model: W2

## 1. INTRODUCTION

---

This manual provides detailed instructions for the safe and effective use of your Comfyer Heated Wearable Breast Pump Model W2. Please read this manual thoroughly before initial use and retain it for future reference. The Comfyer breast pump features industry-leading heating technology, hospital-grade suction, and smart app control for a comfortable and efficient pumping experience.



Image 1.1: Comfyer Heated Wearable Breast Pump and App Interface

## 2. WHAT'S IN THE BOX

Your Comfyer Heated Wearable Breast Pump package includes the following components:

- 1 x Pump Motor
- 1 x Flange
- 1 x Milk Collector
- 1 x Dust Cover
- 1 x Valve
- 1 x Silicone Diaphragm
- 3 x Flange Inserts (17mm, 19mm, 21mm)
- 1 x USB Charging Cable
- 1 x Nipple Size Ruler
- 20 x Breast Milk Storage Bags

## WHAT'S IN THE BOX? (1 PACK)



Image 2.1: Package Contents

### 3. SETUP AND ASSEMBLY

#### 3.1 Initial Sterilization

Before the first use, sterilize all washable components (flange, milk collector, valve, silicone diaphragm, and flange inserts) by boiling or steaming them for 5 minutes. Ensure all parts are completely dry before assembly.

#### 3.2 Choosing the Correct Flange Size

Selecting the correct flange size is crucial for comfort and milk output. Use the included nipple size ruler to measure your nipple diameter after it has been stimulated. Choose the flange insert that provides a comfortable, snug fit without compressing the nipple or areola.

- **17mm Insert:** For nipple measurement 9-11mm
- **19mm Insert:** For nipple measurement 13-15mm
- **21mm Insert:** For nipple measurement 17-20mm

### CHOOSE THE CORRECT FLANGE SIZE

**1**

Stimulate nipple until it perks up

Nipple Measurement	Flange
9-11mm	13mm
11-13mm	15mm
13-15mm	17mm
15-17mm	19mm
17-20mm	21mm
20-23mm	24mm
23-26mm	27mm
26-29mm	30mm
29-32mm	32mm

**2**

Too Small

Correct Fit

Too Large

Image 3.1: Flange Size Selection Guide

#### 3.3 Assembly Steps

1. Insert the chosen flange into the milk collector.
2. Place the silicone diaphragm securely over the milk collector.
3. Attach the valve to the bottom of the milk collector.
4. Align the assembled milk collector with the pump motor and press firmly to ensure a secure, leak-

proof seal.

5. Ensure the duckbill valve, flange tunnel, diaphragm, and pump hub are firmly sealed before use to prevent leaks and ensure optimal suction.

## 4. OPERATING INSTRUCTIONS

### 4.1 Powering On/Off and Mode Selection

The pump features a central power button and mode switch. Press and hold the power button to turn the pump on or off. Short press the power button to cycle through the three available modes:

- **Heating Mode:** Reduces clogs and improves milk flow with gentle warmth.
- **Stimulation Mode:** Mimics a baby's natural sucking to initiate milk let-down.
- **Expression Mode:** Provides steady suction for efficient milk collection.



Image 4.1: Pumping Modes and Levels

### 4.2 Adjusting Suction and Heating Levels

Use the '+' and '-' buttons on the pump motor to adjust the intensity of suction (9 levels) or the heating level, depending on the selected mode. The heating technology operates between 35°C (95°F) and 45°C (113°F) to help prevent clogs and enhance milk flow.



Image 4.2: Heating Technology for Enhanced Milk Flow

### 4.3 App Control

For a hands-free experience and customized settings, control your Comfyer breast pump using the dedicated 'comfyer' app. The app allows you to customize modes, intensity, and pumping duration.

1. Download the 'comfyer' app from the App Store (OS) or Google Play (Android). The app is compatible with iOS 13+ and Android 10+.
2. Create your account and log in.
3. Turn on Bluetooth on your device and power on the breast pump. The device should connect automatically.
4. If the device does not connect automatically, you can manually add it by pressing the '+' button within the app.



Image 4.3: App Connection Guide

## 4.4 Magnetic Charging

The pump features an innovative magnetic charging base that also serves as a storage solution for the pump motor and charging cable. Simply place the pump motor onto the magnetic base to charge. A full charge provides approximately 180 minutes of usage or 8 pumping sessions.



Image 4.4: Magnetic Charging and Battery Life

## 5. MAINTENANCE AND CLEANING

Proper cleaning and maintenance ensure hygiene and prolong the life of your breast pump. All skin-contact parts are BPA-free and food-safe.

### 5.1 Disassembly for Cleaning

To clean, carefully disassemble the pump into its four main washable parts: the flange, diaphragm, valve, and milk collector.



Image 5.1: Parts to Clean

## 5.2 Cleaning Instructions

- Wash all disassembled parts with warm soapy water immediately after each use.
- Rinse thoroughly with clean water.
- Allow parts to air dry completely on a clean surface or use a dedicated drying rack.
- Do not wash the pump motor. Wipe it clean with a damp cloth if necessary.

## 6. TROUBLESHOOTING

---

If you encounter any issues with your Comfyer breast pump, please consider the following:

- **Low Suction:** Ensure all pump components are correctly assembled and securely sealed. Check for any damage to the silicone parts (diaphragm, valve, flange).
- **No Power:** Verify the pump is fully charged. Connect it to the magnetic charging base and ensure the charging indicator light is on.
- **App Connectivity Issues:** Ensure Bluetooth is enabled on your mobile device and the pump is powered on. Try restarting both the app and the pump.
- **Discomfort During Pumping:** Re-evaluate your flange size using the nipple size ruler. Adjust suction levels to a comfortable setting.

For further assistance, please contact Comfyer customer support.

## 7. SPECIFICATIONS

---

Feature	Specification
Model Number	W2
Brand	Comfyer
Material	Polypropylene, Silicone
Noise Level	45 Decibels
Special Feature	Hands-Free Usage, Heating Technology, App Control
Battery Type	2 AAA batteries (included) - <i>Note: Product description mentions 2000mAh battery life, suggesting rechargeable. AAA batteries might be for a remote or specific component, or a data discrepancy. For the main pump, assume rechargeable as per feature bullets.</i>
Package Dimensions	7.17 x 5.94 x 3.39 inches
Item Weight	9.03 ounces

## 8. WARRANTY AND SUPPORT

For warranty information, product support, or to purchase replacement parts, please visit the official Comfyer store on Amazon or contact Comfyer customer service directly. You can find more information at the Comfyer Store.