



Manuals.plus /

› Morelax /

› Morelax ER262 1080P WiFi Ear Wax Remover Otoscope Instruction Manual

Morelax ER262

Morelax ER262 1080P WiFi Ear Wax Remover Otoscope

Model: ER262 | Brand: Morelax

INTRODUCTION

The Morelax ER262 is an advanced 1080P WiFi ear wax remover otoscope designed for visual ear inspection and cleaning. Equipped with a 360° camera and 6 LED lights, it provides a clear view of the ear canal directly on your smartphone or tablet. Its integrated design and IP67 waterproof rating ensure safe and hygienic use for adults, children, and even pets.

PACKAGE CONTENTS

Please verify that all items listed below are included in your package:

- 1 x Otoscope with Light
- 1 x Type-C Charging Cable
- 6 x Silicone Ear Spoons
- 1 x User Manual
- 1 x 8-piece Earwax Removal Tool Kit



Image: All components included in the Morelax ER262 package, including the otoscope, charging cable, ear spoons, and additional cleaning tools.

PRODUCT OVERVIEW

The Morelax ER262 otoscope combines advanced technology with user-friendly design for effective ear care.

High-Definition Camera and Lighting

The device features a built-in 5-million-pixel high-precision lens, offering a 1080P HD resolution. Its 360° wide-angle design, complemented by 6 bright LED lights, ensures a comprehensive and clear view of the ear canal, allowing for precise and hygienic cleaning.



Image: Detailed view of the otoscope's 5-million-pixel 1080P HD camera and 6 LED lights, illustrating its wide-angle and illumination capabilities.

Integrated Ear Spoon Design

The ear spoon is integrated directly with the camera, preventing accidental detachment during use. This design enhances safety and simplifies the cleaning process. The silicone ear covers are soft and skin-friendly for comfort.



Image: Illustration of the integrated ear spoon design, highlighting its secure attachment to the otoscope to prevent accidental detachment.

IP67 Waterproofing

The otoscope lens and ear spoons are IP67 waterproof, allowing for easy cleaning with water or alcohol wipes after each use. **Note: Do not immerse the entire otoscope device in water.**

Q	About cruising ability	▼
A	It can be fully charged for 90 minutes and standby for up to 30 days.	
Q	About waterproof function	▼
A	IP67 waterproof ear lens, earplugs, lenses can be cleaned with water or wiped with alcohol swab after use, but don't put the whole ear lens into the water.	
Q	Will the app that I must download and use reveal my personal information?	▼
A	The app we use only requests access to your camera to view the ear canal on your phone.	
Q	Is the app free to use?	▼
A	The app is completely free to download and use.	
Q	Which smartphones is it compatible with?	▼
A	The application is compatible with any IOS and Android systems.	
Q	About photo storage	▼
A	The photos and videos taken by the device will be stored in the App album.	

Image: The otoscope lens shown with water droplets, demonstrating its IP67 waterproof capability for easy cleaning.

Battery and Compatibility

The device is equipped with a 350mAh rechargeable battery, providing approximately 90 minutes of continuous use on a full charge and up to 30 days of standby time. It is compatible with iOS, iPadOS, and Android devices, connecting via a dedicated mobile application.



MORELAX Ear Cleaning leads in innovative ear care solutions, prioritizing hygiene and comfort. With user-centric design, our gentle, ergonomic tools ensure safety and satisfaction.

Driven by excellence, we explore new technologies to enhance ear health. Experience superior ear care with MORELAX – where precision and advanced technology meet for optimal hygiene and comfort.

Image: The otoscope shown with various mobile devices (iPhone, Android, iPad), illustrating its wide compatibility across different operating systems.

SETUP

1. App Download

Scan the QR code provided in the manual or on the product packaging to download the dedicated application (e.g., "AIR-Look" or "HNDvisual"). Alternatively, search for the app in your device's app store (Apple App Store for iOS/iPadOS, Google Play Store for Android).



Image: A user interacting with a smartphone, showing the four steps for app connection: App Download (via QR code), Power On, Connect WiFi, and Open App to Start. [Download App Here](#)

2. Power On the Device

Press and hold the power button on the otoscope for 2 seconds to turn it on. The indicator light will illuminate.

3. Connect to WiFi

Go to your device's WiFi settings and connect to the network named "JesHome-P6S-XXX" (where XXX is a unique

identifier). No password is required for this connection.

Video: Demonstrates the setup process, including downloading the app, powering on the device, and connecting to its WiFi network. (ACWOO Earwax Remover with Camera)

4. Open the App and Start

Once connected to the WiFi, open the downloaded app and press the "Start" button to begin viewing the live feed from the otoscope camera.

OPERATION

1. Using the App Interface

The app provides various functions such as taking photos, recording videos, and a mirror function to adjust the image orientation for easier use. Photos and videos taken will be stored in the app's gallery.

2. Ear Cleaning Process

Carefully insert the otoscope into the ear canal while observing the live feed on your device. Use the provided silicone ear spoons to gently remove earwax. Always proceed with caution and avoid inserting the device too deeply.

Visual Smart APP Precisely Clear Hard Earwax



Image: A person using the Morelax ER262 otoscope to clean their ear, with the live camera feed visible on a smartphone.

Video: Illustrates the practical use of an ear otoscope for cleaning, showing how to insert and maneuver the device while viewing the ear canal. (Ear Scoop Cleaner Otoscope Cleaning Kit)

3. Safety Guidelines

- Always use the device gently and avoid applying excessive force.
- Do not use the product while it is charging.
- If you experience any pain or discomfort, immediately stop using the device.
- Keep out of reach of small children without adult supervision.

MAINTENANCE

Cleaning the Device

After each use, clean the otoscope lens and silicone ear spoons. The IP67 waterproof lens can be cleaned with water or wiped with alcohol wipes. **Do not submerge the entire otoscope body in water.**

Replacing Ear Spoons

To replace an ear spoon, gently pull off the used silicone cover from the otoscope tip. Select a new silicone ear spoon from the kit and carefully push it onto the tip until it is securely in place.

TROUBLESHOOTING

If you encounter any issues, please refer to the following common solutions:

- **Device refused to connect?**

If the device is connected to another phone, the app cannot use the product. If the phone has set WiFi that cannot connect to the Internet or the phone detects abnormal WiFi, it will automatically replace it with its own WiFi, then the WiFi of the product is disconnected and cannot be used.

- **WiFi connection unstable?**

If the WiFi prompts "connected but no Internet access, whether to switch", please click "No". If you click "Yes", it automatically switches to an Internet network and the product cannot be used. If your own Internet WiFi coverage area is strong, the device WiFi will be automatically replaced. At this point, it is recommended to go to areas with weak or no WiFi to connect the product's own WiFi and use it. If you want to surf the Internet while using this product, you can only use your own traffic to surf the Internet, and the WiFi is still connected to their own products. If you frequently connect to the Internet, you can only delete all your connected Internet WiFi.

- **About suitable parts for earwax removal camera?**

This is for adults, children over 6 years and pets examining ears, nose, teeth and other body parts.

- **Whether the ear spoon can be removed?**

In view of the problem that the ear spoon is easy to fall off on the market, we have optimized and improved the camera with the ear spoon to achieve an integrated design, so it cannot be removed, but the silicone ear spoon can be replaced.

- **About cruising ability?**

It can be fully charged for 90 minutes and standby for up to 30 days.

- **About waterproof function?**

IP67 waterproof ear lens, earplugs, lenses can be cleaned with water or wiped with alcohol swab after use, but don't put the whole ear lens into the water.

- **Will the app that I must download and use reveal my personal information?**

The app we use only requests access to your camera to view the ear canal on your phone.

- **Is the app free to use?**

The app is completely free to download and use.

- **Which smartphones is it compatible with?**

The application is compatible with any iOS and Android systems.

- **About photo storage?**

The photos and videos taken by the device will be stored in the App album.

- **Is a password required for the default WiFi?**

Connect it to your WiFi without a password.

- **The APP connection is unstable.**

Close and restart the product, and then reconnect to the APP.

- **Why does the device ring off?**

Reminds you to use the equipment to dig too deep, to pay attention to safety.

- **Do I need any water for the operation?**

Earwax is easily removed with forceps without repeated irrigation or water accumulation.

- **Is it painful to operate it?**

The product is soothing, non-stimulating, without contact with the ear canal.

- **Not turn on or light up?**

Please charge the item fully first to have a check. Take another charger to have a try.

SPECIFICATIONS

Feature	Detail
Brand	Morelax
Manufacturer	Morelax
Model Number	ER262
Product Dimensions	11 x 10 x 149 cm; 18 g
ASIN	B0DX1M5752
Camera Resolution	1080P HD, 5 Million Pixels
Field of View	360° Wide Angle
Lighting	6 LED Lights
Waterproof Rating	IP67 (Lens and Ear Spoons)
Battery Capacity	350mAh Lithium-ion
Usage Time	Up to 90 minutes
Standby Time	Up to 30 days
Connectivity	WiFi (App-controlled)
Compatibility	iOS, iPadOS, Android
Charging Port	Type-C
Country of Origin	China

WARRANTY AND SUPPORT

For warranty information, technical support, or service inquiries, please refer to the contact details provided on the product packaging or the official Morelax website. Keep your purchase receipt for warranty claims.