

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

Manuals.plus /

- › [Widex](#) /
- › Widex Instant Double Ear-Tip Instruction Manual

Widex Instant Double Ear-Tip

Widex Instant Double Ear-Tip Instruction Manual

Model: Instant Double Ear-Tip (Medium)

INTRODUCTION

The Widex Instant Double Ear-Tip is designed to provide a secure and comfortable fit for your hearing aids, enhancing sound quality and reducing background noise. Its dual-layered design ensures stability and comfort for extended wear. This manual provides essential information on the proper use, maintenance, and care of your ear tips.





Image: Widex Instant Double Ear-Tip, Medium size. This image illustrates the dual-layered silicone design of the ear tip, highlighting its shape and material.

IMPORTANT SAFETY INFORMATION

- **Choking Hazard:** Keep ear tips out of reach of children and pets. Small parts can be swallowed.
- **Proper Fit:** Ensure the ear tip is correctly attached to your hearing aid before insertion. A loose ear tip may become dislodged in the ear canal.
- **Discomfort:** If you experience pain, irritation, or discomfort, remove the ear tip immediately and consult a hearing care professional.
- **Cleaning:** Clean ear tips regularly as instructed to prevent earwax buildup and maintain hygiene.
- **Damage:** Do not use damaged or torn ear tips. Replace them promptly to ensure optimal performance and safety.

PACKAGE CONTENTS

Your Widex Instant Double Ear-Tip package typically includes:

- Widex Instant Double Ear-Tips (quantity as purchased, e.g., 1 pair)

Note: Hearing aids are sold separately and are not included with the ear tips.

SETUP

1. Identifying Your Ear Tip Size

The Widex Instant Double Ear-Tips are available in various sizes (XSmall, Medium, Large) to accommodate different ear canal anatomies. Ensure you are using the correct size for a comfortable and secure fit. The current product is a Medium size ear tip.

2. Attaching the Ear Tip to Your Hearing Aid

1. Gently hold the hearing aid receiver or sound tube.
2. Align the opening of the ear tip with the end of the receiver/sound tube.
3. Push the ear tip firmly onto the receiver/sound tube until it is securely seated. Ensure there are no gaps and the ear tip cannot be easily pulled off.
4. Verify that the ear tip is fully attached and stable before inserting it into your ear.





Image: Attaching the ear tip. This image demonstrates the process of securely attaching the ear tip to the hearing aid's receiver or sound tube.

OPERATING

1. Proper Insertion

1. Hold the hearing aid with the attached ear tip.
2. Gently guide the ear tip into your ear canal until it feels snug and comfortable. Do not force it.
3. Adjust the hearing aid behind your ear until it sits securely.

2. Ensuring a Secure Fit and Optimal Performance

- **Double-Layered Design:** The dual-layered structure of these ear tips is engineered to provide a firm yet comfortable hold within the ear canal, minimizing movement and enhancing sound delivery.
- **Soft & Flexible Material:** Made from skin-friendly material, the ear tips conform to your ear's shape, ensuring all-day comfort and reducing the likelihood of slipping.
- **Enhanced Noise Isolation:** The design helps to reduce ambient background noise, allowing for

improved speech clarity and a better listening experience.

- **Universal Compatibility:** These ear tips are compatible with a variety of hearing aids, including 3.0 Signia hearing aids and most Widex hearing aids, offering a versatile solution for many users.

MAINTENANCE

1. Cleaning Instructions

Regular cleaning of your ear tips is crucial for hygiene and optimal performance.

1. Gently remove the ear tip from the hearing aid receiver/sound tube.
2. Wipe the ear tip with a soft, damp cloth. You may use a mild, non-alcoholic cleaning solution specifically designed for hearing aid accessories.
3. Ensure no moisture enters the hearing aid itself.
4. Allow the ear tip to air dry completely before reattaching it to the hearing aid.
5. Do not use harsh chemicals, abrasive materials, or immerse the hearing aid in water.

2. Replacement Schedule

Ear tips should be replaced regularly to maintain hygiene, comfort, and acoustic performance. The frequency of replacement depends on usage, earwax production, and cleaning habits. It is generally recommended to replace ear tips every 2-4 weeks, or sooner if they show signs of wear, discoloration, or damage.

TROUBLESHOOTING

If you encounter any issues with your ear tips, refer to the following common problems and solutions:

- **Ear tip feels loose or falls out:**
 - Ensure the ear tip is fully seated on the hearing aid receiver/sound tube.
 - You may be using an incorrect size. Consider trying a larger size for a more secure fit.
 - The ear tip might be worn out. Replace it with a new one.
- **Sound quality issues (e.g., muffled sound, feedback):**
 - Check for earwax buildup in or on the ear tip. Clean thoroughly.
 - Ensure the ear tip is properly inserted and creating a good seal in your ear canal.
 - The ear tip might be damaged or torn, affecting sound. Replace it.
- **Discomfort or irritation:**
 - Ensure you are using the correct size ear tip. A tip that is too large or too small can cause discomfort.
 - Verify that the ear tip is clean and free of debris.
 - If discomfort persists, discontinue use and consult a hearing care professional.

SPECIFICATIONS

Product Name	Widex Instant Double Ear-Tip
---------------------	------------------------------

Model	Instant Double Ear-Tip
Size	Medium (also available in XSmall, Large)
Material	Soft, flexible, skin-friendly silicone
Compatibility	Compatible with 3.0 Signia hearing aids and most Widex hearing aids.
Manufacturer	Widex
ASIN	B0DXHSLZM5

WARRANTY AND SUPPORT

Warranty Information

Specific warranty details for the Widex Instant Double Ear-Tip are typically provided by the retailer or manufacturer at the time of purchase. Please refer to your purchase documentation or contact your point of sale for warranty terms and conditions.

Customer Support

For further assistance, questions, or concerns regarding your Widex Instant Double Ear-Tips, please contact your hearing care professional or the retailer from whom you purchased the product. They can provide specialized support and guidance.

© 2025 Widex. All rights reserved.

This manual is for informational purposes only. Specifications are subject to change without notice.