

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Avistar](#) /
- › [Avistar Extra Strong Denture Adhesive Cream User Manual](#)

Avistar AV-DA5

Avistar Extra Strong Denture Adhesive Cream User Manual

Model: AV-DA5

1. PRODUCT OVERVIEW

The Avistar Extra Strong Denture Adhesive Cream is designed to provide a secure and comfortable hold for full and partial dentures. This formula offers an extra-strong, waterproof grip for 8-12 hours, effectively sealing out food particles and reducing gum irritation. It features a refreshing mint flavor and is formulated without zinc or fluoride, prioritizing user health and safety.



Image: Avistar Extra Strong Denture Adhesive Cream tube and its retail packaging, highlighting the product's appearance.

Key Features:

- **Extra Strong Hold:** Provides 8-12 hours of secure, waterproof grip.
- **Superior Sealing:** Effectively seals out food particles, enhancing comfort and reducing irritation.
- **Zinc & Fluoride Free:** Formulated with safe, food-grade ingredients without zinc or fluoride.
- **Easy Application:** Features a precision nozzle for controlled and accurate application.
- **Refreshing Mint Flavor:** Leaves a fresh sensation after application.

2. SAFETY INFORMATION AND PRECAUTIONS

Before using Avistar Denture Adhesive Cream, please read and understand the following safety guidelines:

- Do not use more product than instructed. Excessive use may lead to oozing.
- Do not use if you have a known allergy or sensitivity to any of the ingredients listed.

- Keep out of reach of children.
- If irritation or discomfort occurs, discontinue use and consult a dental professional.
- Do not use on broken or irritated oral tissue.
- Store in a cool, dry place.



NO HARMFUL INGREDIENTS



**ZINC FREE
FLOURIDE FREE**

Image: Visual confirmation that the product is Zinc Free and Fluoride Free, emphasizing its safe formulation.

3. INGREDIENTS AND COMPOSITION

The Avistar Denture Adhesive Cream is formulated with the following ingredients:

- Calcium/sodium PVM/MA copolymer
- Petrolatum
- Cellulose gum
- Mineral oil
- Menthol (for mint flavor)

- Paraffin
- Methyl-P-hydroxybenzoate
- Silica
- Vaseline

100% SAFE INGREDIENTS

 <p>CALCIUM SODIUM PVM/MA COPOLYMER</p> <p>EXTRA STRONG WATERPROOF BOND</p>	 <p>CARBOXYMETHYL CELLULOSE (CMC)</p> <p>COMFORTABLE ADHESION AND BLOCKING FOOD</p>	 <p>PARAFFIN</p> <p>SMOOTH EASY APPLICATION</p>
 <p>METHYL-P-HYDROXYBENZOATE</p> <p>BUILD A PROTECTIVE BARRIER</p>	 <p>SILICA</p> <p>CREATES A COMFORTABLE BARRIER</p>	
 <p>MENTHOL</p> <p>FRESH MINT SENSATION</p>	 <p>VASELINE</p> <p>PROVIDES THICKENING FOR STICKINESS</p>	

 **AVISTAR**

Image: An infographic illustrating the 100% food-grade and safe ingredients used in the Avistar Denture Adhesive Cream, explaining the function of each component.

4. APPLICATION INSTRUCTIONS

Follow these steps for optimal adhesion and comfort:

1. **Clean and Dry Dentures:** Ensure your dentures are thoroughly cleaned and completely dry before applying the adhesive. Any moisture can reduce the effectiveness of the hold.
2. **Apply Adhesive in Thin Lines:** Apply a small amount of adhesive in thin, continuous lines on the denture surface that will come into contact with your gums. For upper dentures, apply along the ridge and in the center. For lower

dentures, apply along the ridge. For partials, apply to the areas that contact your gums. Avoid applying too close to the edge of the denture to prevent oozing.

- 3. Insert Dentures and Hold Briefly:** Place the dentures firmly into your mouth and hold them in place for a few seconds to allow the adhesive to set.

HOW TO USE



1. Clean and Dry Dentures



2. Apply Adhesives in Thin Lines



3. Insert Dentures and Hold Briefly in Place

Image: A step-by-step visual guide on how to apply the denture adhesive, showing cleaning, precise application on different denture types, and insertion.



Image: A close-up view demonstrating the precision nozzle of the Avistar tube, allowing for easy and accurate application of the cream onto a denture.

5. MAINTENANCE AND REMOVAL

Proper maintenance ensures hygiene and extends the life of your dentures.

Removing Dentures:

To remove your dentures, rinse your mouth with warm water. Gently rock the denture back and forth until it loosens. If removal is difficult, a warm saltwater rinse or a denture cleanser soak may help.

Cleaning Dentures and Adhesive Residue:

- After removal, clean your dentures thoroughly with a soft-bristle brush and denture cleanser to remove any adhesive residue and food particles.
- Rinse your mouth with warm water to remove any remaining adhesive from your gums, palate, and tongue. A soft brush or cloth can be used to gently wipe away residue from oral tissues.

- For stubborn adhesive residue on dentures, consider using a specialized denture cleaning tablet or soaking solution.

6. TROUBLESHOOTING

Here are solutions to common issues you might encounter:

Issue	Possible Cause	Solution
Denture not holding securely or hold is short-lived.	<ul style="list-style-type: none"> ◦ Denture not completely dry before application. ◦ Insufficient amount of adhesive applied. ◦ Uneven application. ◦ Denture fit may be poor. 	<ul style="list-style-type: none"> ◦ Ensure dentures are completely dry. ◦ Apply a thin, continuous line of adhesive as instructed. ◦ Reapply if necessary, but avoid excessive use. ◦ Consult your dentist for denture refitting if the problem persists.
Adhesive oozes out from under dentures.	<ul style="list-style-type: none"> ◦ Too much adhesive applied. ◦ Adhesive applied too close to the denture edge. 	<ul style="list-style-type: none"> ◦ Use less adhesive. Start with a small amount and add more if needed. ◦ Apply adhesive slightly inward from the denture edge.
Difficulty removing adhesive residue from gums or dentures.	<ul style="list-style-type: none"> ◦ Adhesive has dried or hardened. 	<ul style="list-style-type: none"> ◦ Rinse mouth with warm water or warm saltwater. ◦ Gently brush gums with a soft toothbrush. ◦ For dentures, soak in a denture cleaning solution or use a denture brush.

7. PRODUCT SPECIFICATIONS

Specification	Detail
Product Name	Avistar Extra Strong Denture Adhesive Cream
Model Number	AV-DA5
Brand	Avistar
Item Form	Cream
Hold Duration	8-12 hours
Key Features	Waterproof, Zinc Free, Fluoride Free, Mint Flavor, Precision Nozzle
Package Dimensions	6.18 x 2.44 x 1.69 inches
Item Weight	3.84 ounces
Country of Origin	China



**EXTRA STRONG,
WATERPROOF GRIP**



**SEALS OUT FOOD,
NO IRRITATION**

Image: Illustrates the product's long-lasting, waterproof grip and its ability to seal out food, enhancing comfort during meals.



Image: Highlights the fresh mint flavor of the adhesive cream, contributing to a pleasant user experience.

8. WARRANTY AND SUPPORT

Avistar stands behind the quality of its products. If you are not completely satisfied with your purchase, Avistar offers a money-back guarantee.

- **Return Policy:** Products are eligible for a refund or replacement within 30 days of purchase.
- **Customer Satisfaction Guarantee:** If for any reason you are not satisfied, please contact Avistar customer service for assistance.
- **Contact Support:** For further inquiries or support, please visit the official Avistar store on Amazon or refer to the contact information provided on the product packaging.

For more information and to explore other Avistar products, visit the [Avistar Store on Amazon](#).

