

[Manuals.plus](#) /

› [Pride Mobility](#) /

› Pride Mobility Go-Go Super Portable Mobility Scooter SC15 User Manual

Pride Mobility SC15

Pride Mobility Go-Go Super Portable Mobility Scooter SC15 User Manual

Model: SC15

INTRODUCTION

This manual provides essential information for the safe and efficient operation, maintenance, and care of your Pride Mobility Go-Go Super Portable Mobility Scooter SC15. Please read this manual thoroughly before using your scooter to ensure proper understanding of its features and functions.



Image: The Pride Mobility Go-Go Super Portable Mobility Scooter SC15 in Mint color, showcasing its compact design and included front basket.

SAFETY INFORMATION

The Pride Mobility Go-Go Super Portable Mobility Scooter SC15 is an FDA Class II Medical Device, designed to assist individuals with mobility impairments. Always prioritize safety by adhering to the following guidelines:

- Ensure the scooter is fully unfolded and locked into position before use.
- Always use the magnetic charging system correctly to avoid damage.
- Be aware of your surroundings and use LED lights for visibility in low-light conditions.
- Do not exceed the maximum weight capacity of 300 lbs.

- Familiarize yourself with the braking system and speed settings before operating in public areas.

UNPACKING AND INITIAL SETUP

1. Unfolding the Scooter

The Go-Go Super Portable Mobility Scooter features a lightweight and foldable design, allowing it to compact to a slim 12-inch depth without removing the seat. It weighs only 41 lbs (without battery).



The Go Go[®] Super Portable folds quickly (without even having to remove the seat), to a depth of only 12 inches.

Image: A sequence demonstrating the quick and easy folding mechanism of the Go-Go Super Portable Mobility Scooter. To unfold, release the locking mechanism and gently raise the handlebars and seat until they click into their upright, locked positions. Built-in supports allow the scooter to stand upright when folded for convenient storage.



Image: The Go-Go Super Portable Mobility Scooter standing upright next to a car, demonstrating its compact storage capability.

2. Battery Installation and Charging

The scooter is equipped with an 11.7 Ah airline-approved lithium-ion battery. It features Charge 360 technology for effortless magnetic charging.

- Locate the battery compartment, typically under the seat or deck.
- Insert the lithium-ion battery securely.
- Connect the magnetic charger to the charging port. The Charge 360 system will illuminate to signal a successful connection.
- Allow the battery to fully charge before initial use.



Our exclusive *Charge360* snaps into place magnetically and illuminates to signal a successful connection.

Image: A close-up of the Charge 360 magnetic charging port, showing the illuminated connection for ease of use.

OPERATING INSTRUCTIONS

1. Controls Overview

The scooter features an intuitive control panel for easy operation. Key controls include:

- **Speed Dial:** Adjusts the maximum speed of the scooter.
- **Throttle Lever:** Engages forward or reverse motion.
- **Headlights/Taillights:** Front and rear LED lights for enhanced visibility.
- **Battery Indicator:** Displays the current battery charge level.

2. Driving and Braking

The Go-Go Super Portable Mobility Scooter offers smooth acceleration and a reliable braking system.

- To move forward, gently press the throttle lever. Release the lever to engage the automatic braking system.
- The scooter utilizes a dual braking system combining electric and regenerative braking for efficient stops.
- Adjust the speed dial to match your comfort level and environment.

3. Navigating Ramps and Thresholds

The scooter is designed to handle various terrains and inclines. Its powerful motor allows for smooth ascent and descent on ramps. The stable 4-wheel design provides confidence when navigating thresholds.

MAINTENANCE AND CARE

1. Battery Maintenance

- Always charge the lithium-ion battery fully after each use.
- Avoid completely draining the battery to prolong its lifespan.
- Store the scooter and battery in a cool, dry place when not in use.

2. Cleaning and General Inspection

- Clean the scooter regularly with a damp cloth. Avoid using harsh chemicals.
- Periodically check all connections, tires, and moving parts for wear or damage.
- Ensure the carbon fiber and magnesium frame components are free from cracks or damage.



41.6 lbs.
without battery

By utilizing **carbon fiber** and magnesium in the frame, the overall weight remains low without sacrificing strength and durability.

Image: A detailed view of the carbon fiber frame, highlighting the lightweight yet durable construction.

SPECIFICATIONS

Feature	Detail
Model Number	SC15
Weight Capacity	300 lbs.
Turning Radius	51.6"
Overall Length	37.2"
Overall Width	21.5"

Maximum Speed	Up to 4.5 mph
Ground Clearance	3.7" at mid-deck, 1.85" at motor
Front Wheels	7" solid
Rear Wheels	8" solid
Braking System	Regenerative and electromechanical
Seat Size	15.7"W x 14.7"D
Seat-to-Deck	16.6"
Per Charge Range (11.7 Ah battery)	Up to 10 miles at 220 lbs. / Up to 8.8 miles at 300 lbs.
Total Weight	41.6 lbs. (without battery) / 49.9 lbs. (with 11.7 Ah battery)
Battery Requirements	11.7Ah Lithium-ion
Battery Charger	Off-board



GO GO
SUPER PORTABLE

UNFOLD YOUR WORLD

Model Number	SC15
Weight Capacity	300 lbs.
Turning Radius	51.6"
Overall Width ³	21.5"
Overall Length ³	37.2"
Maximum Speed ²	Up to 4.5 mph
Ground Clearance ³	3.7" at mid-deck 1.85" at motor
Front Wheels	7" solid
Rear Wheels	8" solid
Braking System	Regenerative and electromechanical
Seat Size	15.7"W x 14.7"D
Seat-to-Deck	16.6"
Per Charge Range ^{3,4} with 11.7 Ah battery	Up to 10 miles at 220 lbs. Up to 8.8 miles at 300 lbs.
Total Weight	41.6 lbs. (without battery) 49.9 lbs. with 11.7 Ah battery
Battery Requirements	11.7Ah Lithium-ion
Battery Charger	Off-board
Warranty	Frame: 3-year limited Electronics: 2-year limited Drive Motors: 2-year limited Batteries: 2-year limited

Image: A visual representation of the scooter's key specifications.

WARRANTY AND SUPPORT

The Pride Mobility Go-Go Super Portable Mobility Scooter includes a comprehensive warranty:

- **Frame:** 3-year limited warranty.
- **Electronic Components:** 2-year limited warranty.
- **Drive Motor:** 2-year limited warranty.
- **Lithium-ion Battery:** 2-year limited warranty.
- **In-Home Service:** 1-year in-home service for professional repairs.

For customer support or warranty claims, please refer to the contact information provided with your purchase documentation or visit the official Pride Mobility website.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question