

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [VERPAL](#) /

› [VERPAL Electric Wheelchair KR8810 User Manual](#)

VERPAL KR8810

VERPAL Electric Wheelchair KR8810 User Manual

Model: KR8810

1. INTRODUCTION

This manual provides essential instructions for the safe and efficient operation, maintenance, and troubleshooting of your VERPAL Electric Wheelchair KR8810. Please read this manual thoroughly before using the wheelchair to ensure proper function and user safety.



Figure 1: VERPAL Electric Wheelchair KR8810. This image displays the electric wheelchair with its green mesh backrest and seat, black frame, and joystick control.

1.1 Safety Warnings

- Always ensure the user is properly secured with the safety belt before operation.
- Do not operate the wheelchair on inclines exceeding the specified limits.
- Regularly inspect all components for wear and tear, especially brakes and tires.
- Keep hands and clothing clear of moving parts.
- Consult a healthcare professional for guidance on safe usage and suitability for individual needs.

2. SETUP AND INITIAL USE

2.1 Unboxing and Assembly

The VERPAL Electric Wheelchair KR8810 comes largely pre-assembled. Carefully remove the wheelchair from its packaging. No complex assembly is required. Ensure all packing materials are removed before use.

2.2 Battery Charging

1. Locate the charging port on the battery unit.
2. Connect the provided charger to the charging port.
3. Plug the charger into a standard wall outlet.
4. Charge the battery for a minimum of 4 hours for a full charge.
5. The battery features an indicator to show the current charge level. Press the button next to the indicator lights to check the power level.

Note: A fully charged battery ensures optimal performance and range.

3. OPERATING INSTRUCTIONS

3.1 Powering On/Off

1. To power on, insert the key into the ignition slot on the battery and turn it clockwise to the red circular symbol. A red indicator light will illuminate.
2. To power off, turn the key counter-clockwise and remove it.

3.2 Basic Movement

The wheelchair is controlled via a joystick located on the armrest. Gently push the joystick in the desired direction of travel (forward, backward, left, right). Releasing the joystick will automatically engage the electromagnetic brakes, bringing the wheelchair to a stop.



Figure 2: Joystick Controller. The joystick allows for intuitive control of the wheelchair's movement.

3.3 Footrest Adjustment

Ensure the user's feet are comfortably and securely placed on the footrests during operation. The footrests can be adjusted or folded up for easier transfers or storage.



Figure 3: Adjustable Footrests. The footrests provide support and can be folded for convenience.

4. MAINTENANCE

4.1 Battery Care

- Charge the battery regularly, even if the wheelchair is not in frequent use, to prolong its lifespan.
- Avoid completely draining the battery before recharging.
- Store the battery in a cool, dry place.

4.2 Cleaning

Clean the wheelchair regularly with a damp cloth and mild detergent. Avoid using abrasive cleaners or solvents. Ensure no water enters electrical components.

4.3 General Inspection

Periodically check the following:

- Tire pressure and condition.

- Tightness of all nuts and bolts.
- Functionality of brakes and controls.
- Condition of the seat and backrest fabric.

5. TROUBLESHOOTING

- **Wheelchair does not power on:** Check if the battery is charged and properly connected. Ensure the key is turned to the 'on' position.
- **Reduced range or power:** The battery may need a full charge or could be nearing the end of its lifespan. Ensure proper battery maintenance.
- **Unusual noises:** Stop using the wheelchair immediately and inspect for loose parts or obstructions. Contact customer support if the issue persists.

6. SPECIFICATIONS

Feature	Detail
Brand	VERPAL
Model	KR8810
Color	Green
UPC	738191924866
Manufacturer	VERPAL
Date First Available	April 15, 2025

Note: Specifications are subject to change without notice.

7. WARRANTY INFORMATION

For detailed warranty information regarding your VERPAL Electric Wheelchair KR8810, please refer to the warranty card included with your purchase or visit the official VERPAL website. Keep your proof of purchase for warranty claims.

8. CUSTOMER SUPPORT

If you have any questions, require technical assistance, or need to order replacement parts, please contact VERPAL customer support. Contact details can be found on the VERPAL website or on the packaging of your product.