

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

Manuals.plus /

› [HOMCOM](#) /

› HOMCOM 3-in-1 Shower Commode Wheelchair User Manual

## HOMCOM 713-190V00GY

# HOMCOM 3-in-1 Shower Commode Wheelchair User Manual

Model: 713-190V00GY

## 1. INTRODUCTION

Thank you for choosing the HOMCOM 3-in-1 Shower Commode Wheelchair. This versatile product is designed to provide assistance with personal hygiene and mobility for individuals requiring support. It functions as a commode toilet, a shower chair, and a transport wheelchair. Please read this manual thoroughly before assembly and use to ensure safe and proper operation.



Image 1.1: The HOMCOM 3-in-1 Shower Commode Wheelchair in its primary configuration.

## 2. SAFETY INFORMATION

- **Weight Capacity:** This chair is designed to support a maximum user weight of 300 lbs (approximately 136 kg). Do not exceed this limit.
- **Assembly:** Ensure all components are correctly assembled and securely fastened before each use. Refer to the assembly instructions provided.
- **Wheel Locks:** Always engage the rear wheel locks when the chair is stationary, especially during transfers, showering, or commode use, to prevent unintended movement.

- **Stability:** Do not lean excessively to one side or the front/back while seated, as this may compromise stability.
- **Transfers:** Exercise caution during transfers to and from the chair. Ensure the chair is stable and locked. Seek assistance if necessary.
- **Environment:** Use the chair on level surfaces. Avoid uneven terrain or obstacles that could cause tipping.
- **Cleaning:** Regularly clean the chair to maintain hygiene. Refer to the 'Maintenance and Cleaning' section.
- **Inspection:** Periodically inspect the chair for any signs of wear, damage, or loose components. Discontinue use if any issues are found and contact customer support.

**Legal Disclaimer:** Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is a mobility and hygiene aid and should be used as directed.

### 3. PRODUCT OVERVIEW AND FEATURES

The HOMCOM 3-in-1 Shower Commode Wheelchair offers multiple functionalities to assist users in various daily activities.

#### Key Features:

- **3-in-1 Functionality:** Serves as a commode toilet, shower chair, and transport wheelchair.
- **Comfort Enhanced:** Features padded armrests, seat, and back with a durable PU leather cover.
- **Foldable Footrests:** Provides support for feet and can be folded up for easier transfers.
- **Waterproof Construction:** Made with waterproof materials suitable for shower use and easy cleaning.
- **Secure Mobility:** Equipped with four lockable wheels for enhanced safety and easy movement.
- **Detachable Bucket:** Includes a removable bucket for commode use.

# 3-IN-1 FUNCTION



Flip-up Footrest



Commode Toilet



Shower Wheelchair



Bedside Toilet

Image 3.1: Visual representation of the chair's three primary functions.



## COMFORTABLE



Padded Armrests



Padded Seat



Padded Back

Image 3.2: Details of the padded components for user comfort.

## 4. SETUP AND ASSEMBLY

The HOMCOM 3-in-1 Shower Commode Wheelchair requires assembly. Please follow the detailed instructions included in the product packaging. A general overview of the assembly process is provided below:

1. **Unpack Components:** Carefully remove all parts from the packaging and verify against the parts list in your instruction manual.
2. **Assemble Frame:** Connect the main frame components as per the diagrams. Ensure all bolts and screws are tightened securely.
3. **Attach Seat and Backrest:** Secure the padded seat and backrest to the frame.
4. **Install Wheels:** Attach the four wheels to the base of the frame. Ensure they rotate freely and the locking mechanisms are functional.
5. **Attach Footrests:** Install the foldable footrests. Confirm they can be easily flipped up and down.

6. **Insert Commode Bucket:** Place the detachable commode bucket into its designated slot under the seat.
7. **Final Check:** Before first use, double-check all connections for tightness and ensure all moving parts operate smoothly. Test the wheel locks.



Image 4.1: The product is delivered in one package, containing all necessary components for assembly.

## 5. OPERATING INSTRUCTIONS

### 5.1. As a Commode Toilet

- **With Detachable Bucket:** Ensure the bucket is securely placed under the commode seat. Engage all wheel locks before use. After use, carefully remove the bucket for emptying and cleaning.
- **Over a Standard Toilet:** Roll the chair over a standard toilet. Ensure the chair is centered and stable. Engage all wheel locks. The commode seat is designed to fit over most standard toilet bowls.



Image 5.1: The commode wheelchair can be used directly over a standard toilet.



Image 5.2: The commode seat and detachable bucket for hygiene needs.

## 5.2. As a Shower Chair

- Roll the chair into the shower area. Engage all wheel locks to prevent movement during showering.
- The waterproof materials allow for safe use in wet environments.
- After showering, disengage the wheel locks and roll the chair out for drying and cleaning.



## DESIGNED FOR SHOWER

Constructed with waterproof materials, ensuring a hygienic surrounding



Easy to clean with just a wipe

Image 5.3: The chair's waterproof design makes it suitable for shower use.

### 5.3. As a Transport Wheelchair

- Ensure the user is seated comfortably and securely.
- Disengage the wheel locks for mobility.
- The four wheels allow for easy maneuverability within various living spaces.
- For transfers, flip up the footrests to clear the path.

# WHEELS FOR MOBILITY

Rear Wheel Locks



Image 5.4: Rear wheel locks provide stability and security during use.



Image 5.5: The chair can be easily positioned next to a bed for convenient access.

## 6. MAINTENANCE AND CLEANING

Regular cleaning and maintenance will extend the life of your HOMCOM 3-in-1 Shower Commode Wheelchair and ensure hygienic use.

- **General Cleaning:** Wipe down the frame, seat, backrest, and armrests with a damp cloth and mild soap. Rinse thoroughly and dry completely to prevent water spots or residue buildup.
- **Commode Bucket:** The detachable bucket should be emptied and cleaned after each use. Use a disinfectant cleaner and rinse thoroughly.
- **Waterproof Materials:** The chair's waterproof construction allows for easy cleaning. Simply wipe with a cloth.
- **Wheel Maintenance:** Keep the wheels free of hair, lint, and debris to ensure smooth operation. Check wheel locks periodically for proper function.
- **Fasteners:** Periodically check all bolts, screws, and fasteners to ensure they are tight. Tighten if necessary.

- **Storage:** Store the chair in a dry, clean environment when not in use.

## 7. TROUBLESHOOTING

If you encounter any issues with your HOMCOM 3-in-1 Shower Commode Wheelchair, please refer to the following common troubleshooting tips:

- **Chair Instability:** Ensure all assembly bolts are tightened. Verify that the wheel locks are fully engaged when the chair needs to be stationary.
- **Difficulty Rolling:** Check wheels for debris (hair, lint) and remove if present. Ensure wheel locks are disengaged.
- **Footrests Not Operating Smoothly:** Check for any obstructions or damage to the hinge mechanism.
- **General Malfunction:** If problems persist after checking the above, discontinue use and contact HOMCOM customer support for assistance.

## 8. SPECIFICATIONS

Feature	Specification
Model Number	713-190V00GY
Product Dimensions (Overall)	22.4" W x 30.7" D x 37" H (56.9 cm W x 78 cm D x 94 cm H)
Seat Dimensions	16.9" W x 16.9" D (42.9 cm W x 42.9 cm D)
Weight Capacity	300 lbs (136 kg)
Item Weight	26 lbs (11.8 kg)
Material	Waterproof materials, PU leather cover (seat/back/armrests)
Color	Gray
Manufacturer	Aosom LLC



Image 8.1: Key dimensions of the HOMCOM 3-in-1 Shower Commode Wheelchair.

## 9. WARRANTY AND SUPPORT

For warranty information or customer support, please refer to the documentation included with your purchase or contact HOMCOM customer service directly. Keep your proof of purchase for any warranty claims.

**Contact Information:** Please visit the official HOMCOM website or refer to your product packaging for the most current customer service contact details.

