

## MODOY ESWT02/ESWT03

# MODOY ESWT02/ESWT03 Shockwave Therapy Machine Replacement Bullets User Manual

Model: ESWT02/ESWT03 | Brand: MODOY

## 1. INTRODUCTION

---

This manual provides essential information for the proper use and maintenance of your MODOY ESWT02/ESWT03 Shockwave Therapy Machine Replacement Bullets. Please read these instructions carefully before use to ensure safe and effective operation. These replacement bullets are designed specifically for the MODOY ESWT02 and ESWT03 shockwave therapy machines.

## 2. SAFETY INFORMATION

---

- Always ensure the shockwave therapy machine is turned off and unplugged before installing or removing replacement bullets.
- Handle bullets with clean hands to prevent contamination.
- Inspect each bullet for any visible damage or wear before use. Do not use damaged bullets.
- Store replacement bullets in a clean, dry place away from direct sunlight and extreme temperatures.
- Keep out of reach of children.
- These bullets are for professional use or under the guidance of a qualified practitioner.

## 3. PACKAGE CONTENTS

---

Your package should contain the following items:

- 5 x MODOY ESWT02/ESWT03 Shockwave Therapy Machine Replacement Bullets



Image 1: Five replacement bullets for MODOY ESWT02/ESWT03 shockwave therapy machines.

## 4. SETUP AND INSTALLATION

---

Follow these steps to properly install the replacement bullets into your MODOY ESWT02/ESWT03 Shockwave Therapy Machine applicator:

1. **Power Off:** Ensure your shockwave therapy machine is completely powered off and unplugged from the electrical outlet.
2. **Remove Old Bullet:** Carefully unscrew or detach the existing bullet from the applicator head of your shockwave therapy machine. Refer to your machine's specific manual for detailed instructions on removing the applicator components.
3. **Clean Applicator:** Clean the applicator head thoroughly to remove any residue or debris before installing a new bullet.
4. **Install New Bullet:** Take one of the new replacement bullets. Align it with the applicator head and gently screw it in or attach it according to your machine's design. Ensure it is securely fastened but do not overtighten.
5. **Verify Installation:** Double-check that the bullet is firmly in place and there is no wobbling or looseness.

## 5. OPERATING INSTRUCTIONS

---

Once the replacement bullet is securely installed, you can proceed with operating your MODOY ESWT02/ESWT03 Shockwave Therapy Machine. Always refer to your main shockwave therapy machine's user manual for comprehensive operating instructions, treatment protocols, and safety guidelines.

- **Pre-Treatment:** Prepare the treatment area as instructed by your machine's manual or your practitioner. This typically involves applying a coupling gel to the skin.
- **Machine Settings:** Power on your shockwave therapy machine and select the appropriate treatment parameters (e.g., frequency, energy level, number of pulses) as recommended for the specific condition being treated.
- **Application:** Place the applicator head with the installed bullet firmly against the skin in the treatment area. Move the applicator slowly and consistently over the target area, maintaining good contact.
- **Post-Treatment:** After the treatment session, turn off the machine, clean the applicator and bullet, and store the equipment properly.

## 6. MAINTENANCE

---

Proper maintenance of your replacement bullets ensures their longevity and optimal performance.

- **Cleaning:** After each use, clean the bullet thoroughly with a soft cloth and a mild disinfectant solution. Ensure all coupling gel residue is removed. Do not use abrasive cleaners or solvents that could damage the bullet material.
- **Drying:** Allow the bullet to air dry completely before storage.
- **Storage:** Store the clean, dry bullets in their original packaging or a designated clean, dry container to protect them from dust and damage.
- **Inspection:** Regularly inspect the bullets for signs of wear, cracks, or deformation. Replace bullets immediately if any damage is observed.

## 7. TROUBLESHOOTING

---

If you encounter issues with your replacement bullets, consider the following common troubleshooting steps:

- **Bullet Not Fitting Securely:** Ensure you are using the correct replacement bullets for your specific MODOY ESWT02 or ESWT03 model. Check for any debris in the applicator's threading or connection point.
- **Reduced Performance:** If the shockwave output feels weaker or inconsistent, inspect the bullet for wear or damage. A worn bullet may need replacement. Also, ensure the bullet is installed correctly and securely.
- **Unusual Noise:** If you hear unusual noises during operation, immediately power off the machine. Check the bullet for proper installation and any damage. Consult your main machine's manual or contact customer support if the issue persists.

For issues related to the shockwave therapy machine itself, please refer to the machine's dedicated user manual or contact MODOY customer support.

## 8. SPECIFICATIONS

---

| Feature             | Detail  |
|---------------------|---|
| Product Name        | Shockwave Therapy Machine Replacement Bullets |
| Model Compatibility | MODOY ESWT02, ESWT03                          |
| Quantity            | 5 pieces per package                          |
| Package Dimensions  | 1.3 x 1.2 x 0.8 inches                        |
| Package Weight      | 8 ounces                                      |
| Manufacturer        | MODOY   |
| ASIN                | B0FFM9P1LW                                    |

## 9. WARRANTY AND SUPPORT

---

For warranty information and customer support regarding your MODOY ESWT02/ESWT03 Shockwave Therapy Machine Replacement Bullets, please refer to the documentation provided with your original shockwave therapy machine or visit the official MODOY website. You may also contact MODOY customer service directly for assistance with product inquiries, technical support, or replacement parts.

**Contact Information:** Please refer to your product packaging or the MODOY official website for the most current contact details.