

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [ZWP](#) /

› [ZWP Replacement Remote Control Instruction Manual for Yireehome TXR-A1-01-CA, TXR-A1-01M-CA, QNR-A1-01, TNR-A1-01, FDR-A1-01 Adjustable Beds](#)

ZWP TXR-A1-01-CA

ZWP Replacement Remote Control User Manual

COMPATIBLE WITH YIREEHOME TXR-A1-01-CA, TXR-A1-01M-CA, YIREE QNR-A1-01, TNR-A1-01, FDR-A1-01 ADJUSTABLE BEDS

1. Product Overview

This manual provides instructions for the ZWP replacement remote control, designed for use with various Yireehome and Yiree adjustable bed models. It offers a comfortable and ergonomic design for easy operation of your adjustable bed's functions.



Image 1: ZWP Replacement Remote Control. This image displays the remote control with its various buttons for adjusting bed positions.

2. Setup

2.1 Battery Installation

The remote control requires two (2) 1.5V AAA batteries (not included). Follow these steps to install them:

1. Locate the battery compartment cover on the back of the remote control.
2. Slide the cover downwards or press the latch to open it.
3. Insert two AAA batteries, ensuring the positive (+) and negative (-) terminals align correctly with the markings inside the compartment.
4. Replace the battery compartment cover until it clicks securely into place.



Image 2: Battery Installation. This image illustrates the open battery compartment on the back of the remote, showing where to insert two AAA batteries.

2.2 Pairing Instructions

To sync the remote control with your adjustable bed, follow these pairing steps:

1. Ensure your adjustable bed is powered on.
2. On the remote control, simultaneously press and hold the **"head down"** button and the **"foot up"**

button.

3. Continue to hold both buttons for approximately 10 seconds.
4. Release the buttons once you hear a series of beeps. This indicates that the remote control has successfully synced with your adjustable bed.



How To Pairing

Pairing:

press and hold the "head down" button and the "foot up" button at the same time for 10 seconds, Once you hear beeps that will indicate that the remote is synced.

Image 3: Pairing Instructions. This image highlights the 'head down' and 'foot up' buttons on the remote, indicating which buttons to press simultaneously for pairing.

3. Operating the Remote Control

The remote control features intuitive buttons for adjusting your adjustable bed to various positions. Familiarize yourself with the button layout and their functions:

- **Head Up/Down Buttons:** Adjusts the head section of the bed upwards or downwards.
- **Foot Up/Down Buttons:** Adjusts the foot section of the bed upwards or downwards.
- **Combined Head/Foot Buttons:** Simultaneously adjusts both head and foot sections for specific comfort positions.
- **FLAT Button:** Returns the bed to a completely flat position.
- **Light Button:** Activates or deactivates any integrated under-bed lighting or remote light feature.

3.1 Pre-set Positions

The remote control may include pre-programmed buttons for common comfort positions:

- **Zero Gravity (40°/15°):** Elevates the head and foot to create a feeling of weightlessness, reducing pressure on the back.
- **Television (45°/15°):** Positions the bed for comfortable viewing of television or reading.
- **Pregnancy Relaxation (15°/20°):** A gentle incline designed for comfort during pregnancy.
- **Sleeping (0°/0°):** The standard flat position for sleeping.
- **Prevent Snoring (15°/0°):** Slightly elevates the head to help reduce snoring.
- **Nursing Mode (50°/15°):** Provides an elevated position suitable for nursing or feeding.



Image 4: Adjustable Bed Positions. This image displays the remote control alongside diagrams illustrating various pre-set bed positions and their corresponding angles.

4. Maintenance

To ensure the longevity and proper functioning of your remote control, follow these maintenance guidelines:

- **Cleaning:** Wipe the remote control with a soft, dry cloth. Avoid using abrasive cleaners, solvents, or harsh chemicals, as these can damage the surface or internal components.
- **Battery Replacement:** Replace batteries promptly when the remote's performance diminishes or if it

stops responding. Always use new, high-quality AAA batteries.

- **Storage:** Store the remote control in a dry place, away from direct sunlight, extreme temperatures, and high humidity.
- **Protection:** Avoid dropping the remote control or exposing it to excessive force, which could cause internal damage.

5. Troubleshooting

If you encounter issues with your remote control, refer to the following troubleshooting tips:

- **Remote Not Responding:**
 - Check if the batteries are correctly installed and have sufficient charge. Replace them if necessary.
 - Ensure there are no obstructions between the remote control and the adjustable bed's receiver.
 - Verify that the adjustable bed is properly plugged into a power outlet and is turned on.
- **Remote Not Pairing:**
 - Re-attempt the pairing process as described in Section 2.2. Ensure you hold the "head down" and "foot up" buttons simultaneously for the full 10 seconds until you hear beeps.
 - Move closer to the adjustable bed during the pairing process.
- **Intermittent Operation:**
 - This could indicate low battery power. Replace the batteries.
 - Ensure there is no strong electromagnetic interference in the vicinity of the bed.

6. Specifications

Feature	Detail
Brand	ZWP
Manufacturer	ZWP
Compatible Devices	Adjustable Bed (Yireehome TXR-A1-01-CA, TXR-A1-01M-CA, Yiree QNR-A1-01, TNR-A1-01, FDR-A1-01, KSR-A1-01)
Max Number of Supported Devices	1
Color	Black
Special Feature	Ergonomic
Item Weight	2.08 ounces
Package Dimensions	7.4 x 2.5 x 1.1 inches

7. Support

For any questions, concerns, or assistance with your ZWP Replacement Remote Control, please do not hesitate to contact our customer support. We are committed to helping you resolve any issues you may encounter.

Please refer to your purchase documentation or the seller's contact information for specific support channels.