

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [VOCIC](#) /

› [VOCIC AX35 Movable Electric Chair Lift Instruction Manual](#)

VOCIC AX35

VOCIC AX35 Movable Electric Chair Lift Instruction Manual

Model: AX35
Brand: VOCIC

INTRODUCTION

This manual provides comprehensive instructions for the safe and effective operation of your VOCIC AX35 Movable Electric Chair Lift. Designed for portable short-distance transport, this chair assists individuals, particularly seniors, with standing from the floor. Please read all instructions carefully before initial use and retain this manual for future reference.

IMPORTANT SAFETY INFORMATION

- Always use the secure wheel locks when the chair is in use.
- Do not attempt to move the chair while the user is seated.
- Always follow the instructions provided to ensure safe operation and to minimize the risk of falls.
- Ensure the chair is on a stable, level surface before use.
- Do not exceed the maximum weight capacity of 350 LBS.

WHAT'S INCLUDED

The VOCIC AX35 package includes all necessary components for assembly and operation. Please verify that all items listed below are present upon unboxing.



Image: All components of the VOCIC AX35 Electric Chair Lift, including the main unit, charger, hand controller, base feet, hex key wrench,

and screws.

- AX35 Main Unit (x1)
- Charger (x1)
- Hand Controller (x1)
- Base Feet (x2)
- Hex Key Wrench (x1)
- Screws (x4)

ASSEMBLY AND SETUP

Follow these steps to assemble your VOCIC AX35 Electric Chair Lift. A hex key wrench is provided for installation.



Image: Visual guide demonstrating the simple one-step installation process for the AX35 chair lift, involving attaching leg screws.

1. **Unboxing:** Carefully remove all components from the packaging.
2. **Armrest Installation:** Attach the armrests to the main unit using the provided screws and hex key wrench. Ensure they are securely fastened.
3. **Caster Installation:** Attach the four wheels (casters) to the base of the main unit. Ensure they are firmly in place for smooth mobility.
4. **Power Connection:** Connect the hand controller to the main unit.
5. **Initial Charge:** Before first use, fully charge the unit using the provided charger. The charging port is typically located near the controller connection.

Video: A quick installation guide demonstrating the unboxing, armrest, and caster installation steps. (Note: This video may feature a similar product model from a related brand, but the assembly steps are generally applicable.)

OPERATING INSTRUCTIONS

The VOCIC AX35 is designed for user-friendly operation. Familiarize yourself with the controls and functions before assisting a user.

1. **Positioning the Chair:** With the user on the floor, roll the chair lift into position, ensuring the base is stable and the seat is aligned with the user.
2. **Securing the User:** Assist the user onto the seat. Fasten the safety belt securely around the user's waist.
3. **Raising the Seat:** Use the hand controller to slowly raise the seat. The lift mechanism operates smoothly to ensure comfort and safety.
4. **Lowering the Seat:** To lower the seat, use the hand controller to gently descend the user to the desired height or to the floor.
5. **Mobility:** Once the user is securely seated and the chair is in an upright position, the 4 wheels allow for short-distance

transfers between rooms. Always engage the wheel locks when the chair is stationary or when assisting a user on or off the seat.

Video: Official VOCIC product video demonstrating the key features and operational aspects of the AX35 Movable Electric Chair Lift, including raising and lowering functions.

Video: A demonstration of how to use the electric chair lift, highlighting its functionality and ease of operation for assisting users.

PRODUCT FEATURES

The VOCIC AX35 Electric Chair Lift incorporates several features designed for user comfort, safety, and ease of use.

Effortless Room-to-room Transfers

Multi-functional Mobility



AX35



Enables effortless movement

AX07



Unable to move using the wheels

Image: A woman gently pushing the VOCIC AX35 Electric Chair Lift across a room, demonstrating its easy maneuverability.

Equipped with 4 smooth rollers that glide over various floor surfaces, this chair lift is ideal for short-distance transfers between rooms without requiring lifting or carrying, thereby reducing the physical burden on both users and caregivers.

Ultra-Compact Design

Breathable and Comfortable Seats



Image: The VOCIC AX35 Electric Chair Lift shown in both its folded, compact state and its extended, operational state, highlighting its adjustable height and space-saving design.

The chair lift features a four-wheel system while maintaining an ultra-lightweight structure, ensuring agile operation and exceptional stability. The seat and armrests fold flat, allowing easy storage in bathrooms or bedrooms without occupying excessive space.

User-Friendly Operation



Image: A man sitting on the floor next to the VOCIC AX35 Electric Chair Lift, demonstrating the chair's low ground clearance of 2.3 inches for easier access.

With a weight capacity of up to 350 LBS, it is suitable for seniors and individuals with limited physical strength seeking both portability and support. Easy-to-grip handles aid with standing, while secure wheel locks maintain chair stability during use. This design aims to reduce the risk of falls and enhance user confidence in daily activities.

Durable Construction

How to choose chair?



AX07



weighs 23 lbs



Available for bathrooms



AX25



weighs 31 lbs



Easy to transfer chair



AX35



weighs 31 lbs



Easy to transfer chair



AX15



weighs 44.53 lbs



Professional and safer rescue lift chair
Recommended by professional paramedics

Image: A close-up view of the VOCIC AX35's breathable Oxford cloth seat, emphasizing its texture and comfort.

Despite its lightweight design, the lift chair maintains high durability. Its legs are constructed from aluminum material for reliable support and stability, ensuring long-lasting performance.

SPECIFICATIONS

Feature	Detail
Brand	VOCIC

Model Name	AX35
Color	Black
Product Dimensions	15.7"D x 14"W x 19.7"H
Maximum Weight Recommendation	350 Pounds
Material	Aluminum
Item Weight	32 pounds
Form Factor	Wheeled, Armless, T-Shape Base
Included Components	Cushion

MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your VOCIC AX35 Electric Chair Lift.

- **Cleaning:** Wipe down the chair with a damp cloth and mild detergent. Avoid abrasive cleaners.
- **Battery Charging:** Charge the battery regularly, even if not in frequent use, to maintain battery health. Refer to the charger's indicator lights for charging status.
- **Inspection:** Periodically check all screws, bolts, and moving parts for tightness and wear. Tighten any loose fasteners.
- **Storage:** When not in use, fold the chair flat and store it in a dry, safe place away from extreme temperatures.

TROUBLESHOOTING

If you encounter any issues with your chair lift, please refer to the following common troubleshooting tips:

- **Chair not operating:** Ensure the hand controller is securely connected and the battery is charged.
- **Unstable movement:** Check that all casters are properly installed and locked when stationary. Verify all assembly screws are tight.
- **Unusual noises:** Discontinue use immediately and contact customer support.

WARRANTY AND SUPPORT

Your VOCIC AX35 Electric Chair Lift is backed by an industry-leading 1-year warranty. For any setup assistance, troubleshooting, or customization advice, our dedicated care team is available 24/7. Please refer to the warranty card included in your package for detailed terms and conditions. For support, visit the official VOCIC website or contact our customer service hotline.