

PAISEEC W-1

PAISEEC W-1 Foldable Electric Wheelchair User Manual

Model: W-1

1. INTRODUCTION

Thank you for choosing the PAISEEC W-1 Foldable Electric Wheelchair. This 2-in-1 device functions as both an electric wheelchair and a rollator walker, designed for adults and seniors seeking enhanced mobility. It features a powerful 400W Hall Sensor Brushless Motor, the PAI Smart Driving System for confident control, and a convenient 2-step folding mechanism for portability. This manual provides essential information for the safe and effective use, maintenance, and care of your PAISEEC W-1.

2. SAFETY INFORMATION

Your safety is paramount. Please read and understand all instructions before operating the PAISEEC W-1. Failure to follow safety guidelines may result in injury or damage to the device.

- **Weight Capacity:** Do not exceed the maximum weight capacity of 300 lbs (136 kg).
- **Terrain:** Operate on stable, level surfaces. Exercise caution on uneven terrain, slopes, and near curbs. The device is designed for all-terrain capability but user awareness is crucial.
- **Battery Safety:** Use only the provided charger. Do not expose the battery to extreme temperatures or puncture it. Ensure the battery is properly installed before use.
- **Folding/Unfolding:** Ensure the device is powered off and stable before attempting to fold or unfold. Keep hands and fingers clear of pinch points.
- **FDA Clearance:** This product is FDA 510(k) cleared, meeting medical-grade safety standards for electrical safety, stable structure, and safe battery operation.

FDA 510(k) CLEARED

Medical Grade Safety Standard

- ✓ **Electrical Safety**
Leakage-free, anti-interference, passes 1500V safety tests.
- ✓ **Stable Structure**
Anti-tipping, impact-resistant, safe for daily movement.
- ✓ **Safe Battery**
Overcharge/short-circuit protection, UN38.3 certified.



The PAISEEC W-1 is FDA 510(k) cleared, ensuring compliance with medical-grade safety standards.

3. PRODUCT OVERVIEW

The PAISEEC W-1 is a versatile mobility solution. Familiarize yourself with its components for optimal use.

Key Components:

- **Main Frame:** Aircraft-grade aluminum alloy for durability and lightness.
- **Joystick Controller:** For electric wheelchair mode operation.
- **Handlebar:** For rollator mode and assisted pushing.
- **Removable Battery:** Lithium-ion battery for power.
- **Footrest:** Provides support for the user's feet.
- **Wheels:** 8-inch swivel casters (front) and 9-inch pneumatic tires (rear).
- **Seat:** Ergonomically designed for comfort.

Packing List:

- PAISEEC W-1 Electric Wheelchair
- Removable Battery
- Charger
- Storage Basket
- User Manual

Product Dimensions



Frame Weight: **42 lbs (19 kg)**

Packing List



Electric Wheelchair*1



Removable Battery*1



Charger*1



Storage Basket*1



User Manual*1

Visual representation of the PAISEEC W-1's dimensions in both unfolded and folded configurations, along with its complete packing list.

4. SETUP

Unfolding the Wheelchair:

1. Place the folded wheelchair on a flat, stable surface.
2. Step on the folding pedal and hold it down.

3. Gently pull the handlebar upwards.
4. Fold the backrest forwards until the chair locks into its open position.
5. Release the folding pedal. The chair is now fully unfolded.



Visual guide for the 3-step folding process, demonstrating how the wheelchair can be easily folded and stored in a car trunk.

Adjusting Handle Height:

The handle height can be adjusted to suit users from 5'1" to 6'1" for comfortable use in rollator mode or when pushing. Rotate the knurled knob on the handlebar to loosen, adjust to the desired height, and then tighten the knob securely.



Turn Radius32 inchesAdditional FeaturesAdjustable, Foldable, Lightweight, PAI Smart Driving System

10. WARRANTY AND SUPPORT

PAISEEC is committed to providing reliable products and excellent customer service.

Warranty Information:

- **Frame Warranty:** 36 months
- **Battery Warranty:** 18 months

This warranty covers manufacturing defects under normal use. Please retain your proof of purchase for warranty claims.

Customer Support:

For any questions, technical assistance, or warranty inquiries, please contact PAISEEC customer service:

- **Phone Support:** 1 (833) 612-2678 (Mon-Sun 9:00 AM - 6:00 PM EST)
- **Email Support:** support@paiseec.com (Replies within 24 hours)



PAISEEC offers comprehensive customer support, including phone and email assistance, along with warranty and return policies.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details.](#)

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question