

Vaginitis Combo Test Kit B0FP96Z2JJ

Generic Vaginitis Combo Test Kit: At-Home Rapid Screening Instruction Manual

Model: Vaginitis Combo Test Kit B0FP96Z2JJ

INTRODUCTION

This manual provides detailed instructions for the proper use of the Generic Vaginitis Combo Test Kit. This at-home self-test is designed for rapid screening of common causes of vaginal discomfort, specifically detecting *Candida albicans*, *Trichomonas vaginalis*, and *Gardnerella vaginalis*. Please read all instructions carefully before performing the test to ensure accurate results.

IMPORTANT SAFETY INFORMATION

Legal Disclaimer:

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This test is for screening purposes only and should not replace a professional medical diagnosis. Consult a healthcare professional for any health concerns or before making any medical decisions.

Keep out of reach of children. Store the kit at room temperature, away from direct sunlight and moisture. Do not use if the packaging is damaged or expired.

KIT CONTENTS

Each Generic Vaginitis Combo Test Kit contains the following items:

- One (1) Vaginitis Combo Test Device
- One (1) Sterile Vaginal Swab for sample collection
- One (1) Sample Diluent Tube
- One (1) Instruction Manual



**Abnormal
discharge**



Discomfort



**Abnormal odor
Itching**

Early signs

**Don't ignore
these symptoms**

Image: Contents of the Vaginitis Combo Test Kit, showing the individually packaged test device, sterile swab, and sample diluent tube.

SETUP AND PREPARATION

1. **Read Instructions:** Carefully read this entire instruction manual before beginning the test.
2. **Gather Materials:** Ensure you have the test kit, a timer, and a clean, flat surface.
3. **Wash Hands:** Thoroughly wash your hands with soap and water before handling any components of the kit.
4. **Open Packaging:** Open the foil pouch containing the test device just before use. Do not touch the reaction areas of the test strip.
5. **Check Expiration Date:** Verify the expiration date on the kit packaging. Do not use if expired.



Image: Visual guide to the testing process, illustrating hand washing, sample collection with a swab, mixing the sample, and applying drops to the test device.

OPERATING INSTRUCTIONS: PERFORMING THE TEST

1. Step 1: Sample Collection (Swab)

Carefully remove the sterile vaginal swab from its packaging. Hold the swab by the handle. Gently insert the swab into the vagina, approximately 2 inches deep. Rotate the swab gently for 10-15 seconds to collect a sufficient sample. Ensure the tip of the swab touches the vaginal walls. Withdraw the swab carefully.

**Take control
of your intimate health
and get comfortable again**



Image: A sterile swab and a sample diluent tube, ready for use in the test.

2. Step 2: Sample Preparation (Mix)

Open the sample diluent tube. Insert the swab with the collected sample into the tube containing the diluent solution. Swirl the swab vigorously in the solution for at least 10 seconds, pressing the swab against the inner wall of the tube to release the sample. Remove the swab, pressing it against the side of the tube to extract as much liquid as possible. Discard the swab.

3. Step 3: Apply Sample to Test Device (Drip)

Place the cap firmly back onto the sample diluent tube. Gently invert the tube several times to ensure thorough mixing. Hold the tube vertically and dispense exactly 3 drops of the mixed solution into each of the three sample wells (marked 'S') on the test device.

4. Step 4: Wait for Results

Start a timer immediately after applying the drops. Read the results at 10-15 minutes. Do not interpret results after 30 minutes, as this may lead to inaccurate readings.

RESULT INTERPRETATION

To interpret the results, observe the lines that appear in the result window of the test device. The test device has three

separate test areas for *Candida albicans* (CA), *Trichomonas vaginalis* (TV), and *Gardnerella vaginalis* (GV).



Trusted Care, Right at Home

**Private.**
Test in the privacy and comfort of your own home.

**Reliable.**
High accuracy for results you can trust.

**Empowering.**
Care for your intimate wellness with confidence.

Image: A visual guide to interpreting the test results, showing the test device with control (C) and test (T) lines for *Candida albicans* (CA), *Trichomonas vaginalis* (TV), and *Gardnerella vaginalis* (GV), alongside a table detailing result patterns and associated symptoms.

Vaginitis Combo Test Result Interpretation Guide

Test Item	Test Result	Interpretation
CA (<i>Candida albicans</i>)	Control line (C) and Test line (T) both appear.	Positive for <i>Candida albicans</i> .
	Only Control line (C) appears.	Negative for <i>Candida albicans</i> .
TV (<i>Trichomonas vaginalis</i>)	Control line (C) and Test line (T) both appear.	Positive for <i>Trichomonas vaginalis</i> .
	Only Control line (C) appears.	Negative for <i>Trichomonas vaginalis</i> .
GV (<i>Gardnerella vaginalis</i>)	Control line (C) and Test line (T) both appear.	Positive for <i>Gardnerella vaginalis</i> .
	Only Control line (C) appears.	Negative for <i>Gardnerella vaginalis</i> .
Invalid Result		No Control line (C) appears, or only a Test line (T) appears. Repeat the test with a new kit.

Note: Results are only valid when all three Control (C) lines are present. If a Control line does not appear for any of the test items, the result for that item is invalid and the test should be repeated.

TROUBLESHOOTING

- **No lines appear or only a Test line (T) appears:**The test is invalid. This could be due to insufficient sample, incorrect procedure, or a faulty test kit. Discard the used kit and repeat the test with a new one, carefully following all instructions.
- **Faint lines:** Even a faint line in the Test (T) region should be considered a positive result, provided the Control (C) line is clearly visible. If unsure, repeat the test or consult a healthcare professional.
- **Symptoms persist despite negative results:** If you continue to experience symptoms of vaginal discomfort, consult a healthcare professional for further evaluation and diagnosis. This test is a screening tool and not a substitute for medical advice.

SPECIFICATIONS

Feature	Detail
Product Name	Vaginitis Combo Test Kit
Model Number	B0FP96Z2JJ
Brand	Generic
Test Type	Colloidal gold immunochromatography for Candida albicans, Trichomonas vaginalis, Gardnerella vaginalis
Result Time	10-15 minutes (Do not read after 30 minutes)
Sample Type	Vaginal Swab
Number of Tests	1 per kit
Package Dimensions	6.54 x 2.95 x 1.34 inches
Weight	0.63 ounces
UPC	700758132472

MAINTENANCE

This is a single-use diagnostic kit and requires no maintenance. Dispose of the used kit components according to local regulations for medical waste after use.

WARRANTY AND SUPPORT

For questions regarding the product, its use, or potential issues, please refer to the retailer's customer service or the manufacturer's contact information provided on the product packaging. As this is a diagnostic product, specific warranty details may vary. Always consult a healthcare professional for medical advice.

