

MAHANCRIS ETHR6801

MAHANCRIS Overbed Table with Wheels ETHR6801 Instruction Manual

Model: ETHR6801

1. INTRODUCTION

This manual provides detailed instructions for the assembly, operation, and maintenance of your MAHANCRIS Overbed Table with Wheels, model ETHR6801. This versatile table is designed to offer a flexible workspace or surface for various activities while in bed or seated on a sofa. It features adjustable height, a tilting desktop, and smooth-rolling casters for easy mobility.

2. SAFETY INFORMATION

- **Maximum Weight Capacity:** Do not exceed the maximum weight capacity of 50 lbs (22.7 kg) on the tabletop. Overloading can cause instability and damage.
- **Stability:** Ensure all locking mechanisms (height adjustment knob, tabletop tilt knobs, and wheel brakes) are securely fastened before use to prevent accidental movement or tilting.
- **Movement:** When moving the table, ensure no items are on the tilting section of the tabletop to prevent them from falling. Unlock the casters before moving and lock them once the desired position is reached.
- **Assembly:** Follow all assembly instructions carefully. Ensure all parts are correctly installed and tightened.
- **Children and Pets:** Keep children and pets away from the table during assembly and ensure they do not climb on or play with the table to prevent injury.
- **Surface Protection:** Avoid placing hot items directly on the tabletop. Use coasters or protective pads.

3. PACKAGE CONTENTS

Verify that all components listed below are present before beginning assembly. If any parts are missing or damaged, please contact customer support.

- 1 x Height Adjustable Laptop Table (disassembled components)

- 1 x Accessory Kit (hardware, tools)
- 1 x Instructions Manual

4. ASSEMBLY INSTRUCTIONS

Assembly of the MAHANCRIS Overbed Table is designed to be straightforward. All necessary hardware and a clear guide are included in the box. Please refer to the included instruction sheet for detailed step-by-step diagrams.

4.1. Overview of Components and Dimensions

Familiarize yourself with the table's structure and dimensions before assembly.



Figure 1: Product Dimensions and Weight Capacity. The table measures 23.6 inches wide by 11.8 inches deep, with an adjustable height from 26 to 33 inches. It supports a maximum weight of 50 lbs.

4.2. Key Features for Assembly

The table includes features that facilitate its functionality and assembly.



Figure 2: Key features include a safety lip on the tilting panel, a height adjustment knob, and smooth casters for mobility.

General Assembly Steps:

1. Attach the casters to the H-shaped base. Ensure they are securely tightened.
2. Connect the vertical support pole to the base.
3. Attach the tabletop panels to the vertical support. The tilting panel will have a separate mechanism for angle adjustment.
4. Secure all connections using the provided hardware and tools. Do not overtighten.

5. OPERATING INSTRUCTIONS

The MAHANCRIS Overbed Table is designed for ease of use and adaptability.

5.1. Height Adjustment

The table's height can be adjusted to four positions, ranging from 26 to 33 inches, to suit various bed or sofa heights and user preferences.



Figure 3: The table offers four height settings: 26, 28, 30, and 33 inches.

To adjust the height:

1. Loosen the height adjustment knob located on the vertical support pole.
2. Carefully raise or lower the tabletop to the desired height.
3. Tighten the height adjustment knob firmly to secure the table at the chosen height.

5.2. Tilting Desktop

One section of the tabletop can be tilted from 0 to 180 degrees, providing an ergonomic angle for reading, writing, or using a laptop. A safety lip is integrated to prevent items from sliding off.



Figure 4: The tabletop can be adjusted to various angles for comfortable use.



Figure 5: Angle adjustment is secured by a locking mechanism.

To adjust the tabletop angle:

1. Loosen the angle adjustment knobs located underneath the tilting tabletop section.
2. Adjust the tabletop to the desired angle.
3. Tighten the adjustment knobs firmly to lock the tabletop in place.

5.3. Mobility

Equipped with four casters, the table can be easily moved. The C-shaped base allows it to slide under most beds or sofas.



Figure 6: The table's C-shaped base allows it to fit under furniture.



Figure 7: The overbed table in use over a bed.

To move the table:

1. Ensure the caster brakes are disengaged.
2. Gently push or pull the table to the desired location.
3. Once in position, engage the caster brakes to prevent unwanted movement.

5.4. Product Demonstration Video

Your browser does not support the video tag.

Video 1: This video demonstrates the MAHANCRIS Overbed Table with Wheels in various settings, highlighting its smooth mobility, adjustable height, and tilting desktop feature for comfortable use in bed or next to a sofa.

6. MAINTENANCE

- **Cleaning:** Wipe surfaces with a soft, damp cloth. Avoid harsh chemicals or abrasive cleaners that could damage the finish.
- **Hardware Check:** Periodically check all screws and bolts to ensure they remain tight. Retighten as necessary.
- **Caster Care:** Keep casters free of debris (hair, dust) to ensure smooth rolling.
- **Storage:** If storing the table for an extended period, ensure it is clean and dry. Store in a cool, dry place.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
Table wobbles or is unstable.	Loose connections; uneven floor; overloaded tabletop.	Ensure all assembly screws are tightened. Move to a level surface. Reduce weight on the tabletop.
Tabletop does not tilt smoothly or stay in position.	Angle adjustment knobs are loose or overtightened; mechanism obstructed.	Loosen the angle adjustment knobs, adjust the tilt, then tighten firmly. Check for any obstructions in the tilting mechanism.
Height adjustment is difficult or does not lock.	Height adjustment knob is too tight or too loose; mechanism obstructed.	Loosen the height adjustment knob, adjust height, then tighten firmly. Ensure the knob aligns with the adjustment holes. Check for obstructions.
Casters do not roll smoothly or lock.	Debris in casters; brakes engaged; casters not properly installed.	Clean any debris from the caster wheels. Disengage brakes. Ensure casters are fully inserted and tightened during assembly.

8. SPECIFICATIONS

Feature	Detail
---------	--------

Feature	Detail
Brand	MAHANCRIS
Model Number	ETHR6801
Product Dimensions (D x W x H)	11.8" x 23.6" x 26"-33" (Adjustable Height)
Item Weight	9 Pounds (4.08 kg)
Maximum Weight Recommendation	50 Pounds (22.7 kg)
Frame Material	Engineered Wood
Top Material Type	Engineered Wood (Particleboard)
Color	Rustic Brown (Top), Black (Frame)
Table Design	Overbed Table (C-Shaped)
Required Assembly	Yes

9. WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the documentation included with your purchase or contact MAHANCRIS customer service directly. Contact details can typically be found on the product packaging or the official MAHANCRIS website.