

SEEDEAR JH-A260

SEEDEAR JH-A260 Rechargeable OTC Hearing Aids Instruction Manual

Brand: SEEDEAR | Model: JH-A260

INTRODUCTION

This manual provides comprehensive instructions for the SEEDEAR JH-A260 Rechargeable OTC Hearing Aids, covering setup, operation, maintenance, and troubleshooting to ensure clear sound quality and comfortable wearing. These hearing aids are designed with advanced features for adults experiencing hearing loss.

Key Features:

- **Clear Voice Technology:** Equipped with an intelligent processing chip and dynamic noise reduction for enhanced speech clarity.
- **Ultra-Discreet Wear:** Micro-design for a virtually invisible fit in the ear canal.
- **All-day Comfort:** Ultra-soft silicone tips conform to the ear canal, eliminating pressure points.
- **Easy to Use:** Automatic power-on/off, 8-second intelligent start-up delay, and large tactile buttons for mode and volume control.
- **Advanced Power Performance:** 2-hour fast charging for up to 48 hours of continuous operation.

WHAT'S IN THE BOX

- Packing Box
- Charging Case
- Hearing Aids x 2
- Eartips x 12
- USB Cable
- Cleaning Brush
- Operating Manual
- Quick Usage Guide



Figure 1: Contents of the SEEDEAR JH-A260 Hearing Aids package.

SETUP

1. Initial Activation and Charging

Before first use, the device needs to be activated. Remove the hearing aids from the charging case. Remove any protective stickers from the charging contacts on the hearing aids. Place them back into the charging case. Connect the charging cable to the charging port of the case. The three LEDs on the case will light up orange, indicating successful activation and charging. Ensure the device is fully charged before initial use.

Your browser does not support the video tag.

Video 1: Instructions on how to activate the SEEDEAR JH-A260 hearing aids for first-time use.

2. Selecting and Attaching Eartips

Choose the appropriate size and style of silicone eartips for a comfortable and secure fit in your ear canal. The correct eartip size is crucial for optimal sound quality and to prevent feedback (whistling sounds).



Figure 2: Selection of various eartips for personalized comfort and fit.

3. Wearing the Hearing Aids

Identify the left (L) and right (R) hearing aids. Gently insert the appropriate hearing aid into your ear canal, ensuring the pull line faces downwards. The eartip should fit snugly without causing pain or discomfort. If you hear a long beeping sound, it may indicate an improper fit; try a different eartip size or adjust the placement.

Your browser does not support the video tag.

Video 2: Demonstration on how to correctly wear the SEEDEAR JH-A260 hearing aids.



Figure 3: Correct wearing position for the hearing aids.

OPERATING INSTRUCTIONS

1. Power On/Off

The hearing aids feature automatic power-on/off. They will turn on automatically when removed from the charging case and turn off when placed back in. There is an 8-second intelligent start-up delay to prevent sudden noise jolts.

To manually power on/off, press and hold the tactile button on the hearing aid for 2 seconds.

2. Adjusting Volume

The device has a single tactile button for volume adjustment. Each short press of the button will adjust the volume in a cycle. You will hear a single "Du" sound with each adjustment. Once the highest volume is reached, you will hear a single "Du" sound again, indicating the cycle has reset to the lowest volume.

Your browser does not support the video tag.

Video 3: Guide on how to adjust the volume of the SEEDEAR JH-A260 hearing aids.

3. Switching Listening Modes

The hearing aids offer two listening modes: Normal Program and Restaurant Program. To switch between programs, quickly press the tactile button twice. Each time the program is switched, you will hear different prompt sounds to indicate the current program:

- "Beep" represents the Normal Program (for daily use).
- "Beep Beep" represents the Restaurant Program (for noisy environments).

Your browser does not support the video tag.


Video 4: Demonstration of how to switch between listening modes on the SEEDEAR JH-A260 hearing aids.





Figure 4: Two available listening modes: Quiet Mode and Noise Reduction Mode.


MAINTENANCE

- **Charging:** For optimal battery health, recharge the charging case every 3 days.
- **Cleaning:** Use the provided cleaning brush to remove wax, dirt, or debris from the eartips and hearing aids as needed. Regular cleaning helps maintain sound quality and device longevity.









Others





	The latest smart chip, more powerful AI intelligent algorithm	Generic chip, backward algorithm
	Ultra-soft silicone ear tips, more comfortable to wear	Ordinary ear tips, wearing uncomfortable
	Intelligent noise reduction, effectively reduce background noise, protect hearing	No effective noise reduction, sound confusion
	Extra large buttons design, easier to use	Difficult to use the mini button

Figure 5: Cleaning and charging tips for the hearing aids.

TROUBLESHOOTING

- **Beeping Sound:** If you hear a long beeping sound while wearing the hearing aid, it may be because they are not fitted correctly. Try changing the eartips to find one that fits your ear canal best, and then refitting them.
- **Unclear or Distant Sounds:** These hearing aids are designed for mild to moderate hearing needs. If sounds feel consistently unclear or distant, consulting a hearing specialist is recommended to optimize your sound experience.

SPECIFICATIONS

Feature	Detail
Model Name	JH-A260
Brand	SEEDDEAR
Color	Light Gray
Material	Silicone
Item Weight	6.4 ounces
Package Dimensions	4.96 x 3.98 x 2.13 inches
Batteries	1 Lithium Polymer batteries required (included)
Battery Weight	43.79 Grams
Manufacturer	SEEDDEAR

WARRANTY & SUPPORT

- **Warranty:** 24 Months Warranty & Lifetime Service.
- **Customer Support:** 24/7 Customer Support available.
- **Money Back Guarantee:** 60-Day Money Back Guarantee.



100,000+ People's choice

We craft hearing aids so those with hearing loss fully experience life's moments-grand-children's laughter, dinner conversations, city sounds-without barriers. No more missing out or feeling isolated. Join the SEEDER family, let's experience this wonderful world together!

Figure 6: Warranty and support information.