

RAYS Scalp Vein Sets



RAYS Scalp Vein Sets Instructions

[Home](#) » [RAYS](#) » RAYS Scalp Vein Sets Instructions 

Contents

- [1 RAYS Scalp Vein Sets](#)
- [2 Product Usage Instructions](#)
- [3 Frequently Asked Questions](#)
- [4 Q: What is the sterilization method used for the Scalp Vein Sets?](#)
- [5 INSTRUCTIONS FOR USE](#)
- [6 DISPOSAL](#)
- [7 Documents / Resources](#)
 - [7.1 References](#)



RAYS Scalp Vein Sets



Product Specifications

- **Purpose of Use:** Sterile, non-toxic, non-pyrogenic, for single use, needle for injections and administration of medicinal solutions.
- **Latex Free:** The device is totally latex free.
- **Sterilization Method:** Sterilized by Ethylene Oxide (EtO), residual amount of gas in accordance with the regulations.
- **Disposal:** Must be made in accordance with applicable government regulations. Incineration or landfill disposal must be carried out under controlled conditions.
- **Conservation:** Store in a cool dry place and far away from any source of heat at a specified temperature.

Product Usage Instructions

1. **Check Expiry Date:** Before use, always check the expiry date of the product.
2. **Inspect Packaging:** Ensure the packaging is intact before use. Do not use if the packaging is damaged.
3. **Prepare Area:** Prepare the injection area carefully and aseptically to prevent infections.
4. **Remove Needle:** Remove the needle from the blister pack just before use.
5. **Therapy Administration:** Proceed with the therapy as prescribed by a healthcare professional.
6. **Proper Disposal:** Dispose of the used needle properly after use to prevent any contamination or injuries.

Frequently Asked Questions

Q: Can the Scalp Vein Sets be reused?

A: No, the Scalp Vein Sets are designed for single use only. Reusing them can lead to infections; always dispose of them after use.

Q: How should I store the Scalp Vein Sets?

A: Store the product in a cool, dry place away from any heat sources to maintain its integrity and sterility.

Q: What is the sterilization method used for the Scalp Vein Sets?

A: The Scalp Vein Sets are sterilized using Ethylene Oxide (EtO) with residual gas amounts complying with regulations to ensure safety.

SCALP VEIN SETS WITHOUT SAFETY MECHANISM

PURPOSE OF USE:

Sterile, non-toxic, non-pyrogenic, for single use, needle for injections and administration of medicinal solutions.

INSTRUCTIONS FOR USE

- Before using check the expire date.
- Check the integrity of the packaging before use
- Open the packaging and use immediately.
- Do not use if packaging is damaged.
- In case of re-use, infections can occur: consult a doctor.
- Prepare the area carefully and aseptically.
- Remove the needle from the blister.
- Go on with the therapy.
- Dispose of it properly after use.

LATEX FREE:

The device is totally latex free.

STERILIZATION METHOD:

Sterilized by Ethylene Oxide (EtO), residual amount of gas in accordance with the regulations.

DISPOSAL

Disposal must be made in accordance with the applicable government regulations. Incineration or landfill disposal must be carried out under controlled conditions. The product that did not come in to contact with infectious agents or other hazardous materials, cannot be considered dangerous waste.

CONSERVATION:

Store in a cool dry place and far away from the any source of heat at a temperature.
These Instructions for Use note must be preserved by the user for the duration of the DM.

RAYS SPA


via Francesco Crispi – 60027 Osimo (AN)

Tel. 071 2868935 – Fax 071 2868911 – www.rays.it – info@rays.it



[RAYS Scalp Vein Sets](#) [pdf] Instructions Scalp Vein Sets, Vein Sets, Sets

References

-  [Home - Rays SpA](#)
- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.