



## Contents [ [hide](#) ]

- [1 Arboleaf CS20N Body Composition Smart Scale](#)
- [2 Specifications](#)
- [3 Measurement Tips](#)
- [4 Trouble Shooting](#)
- [5 Package Contents](#)
- [6 Power Up Scale](#)
- [7 Documents / Resources](#)
  - [7.1 References](#)




## Arboleaf CS20N Body Composition Smart Scale



## Specifications






- **Dimensions:** 11.8" x 11.8" x 1.0" (300mm x 300mm x 25mm)
- **Weight:** 3.61b (1.6kg)
- **Measurement Range:** 0.4- 3961b (0.2-180kg)
- **Measurement Resolution:** 0.21b (0.1kg)
- **Power:** 2×1.5V AAA batteries

## Measurement Tips



Conductive areas

## Scale Display

 Normal working	No display Change battery	 No wireless connection
 Change battery	 Overload	 Calibrating

## Trouble Shooting

### No wireless connection between your scale and the app

1. Reinstall the Arboleaf app. Exit from the app, go to “My Account” -> “Settings” -> “Exit” at the bottom, then delete the app. Go to the Apple App Store or Google Play to install the Arboleaf app again.
2. Turn off Bluetooth, then turn it on again (go to your phone’s Settings).
3. Power on your scale. Pull out the batteries from your scale and put them back in.
4. For Android 6.0 and above, Location needs to be enabled.

### A non-zero number on the scale

Your scale needs another reset.

1. Pull one battery out from your scale and put it back.
2. IMMEDIATELY place the scale on a hard, flat surface.
3. After you see the scale display 0.0, close the battery cover.

## Warranty Registration

<https://arboleaf.com/warranty>

## Measurement Tips

To have accurate, consistent, and comparable measurement results, there are some

DOS and DON'Ts.

## **DO's**

- Place the scale on a flat, hard surface.
- Stand straight on the scale with bare feet, no socks, and no shoes. Both feet uniformly touch the four metal electrodes.
- Do the calibration each time. Wait for "CAL" to return to "0.0" then step on the scale for measurements.
- Use the scale under the same conditions during the same hour each day.
- Empty your bladder before weighing.

## **DON'Ts**

- Do not measure with socks or shoes on.
- Do not measure during or immediately after strenuous exercise.
- Do not measure immediately after shower, bath, or sauna.
- Do not measure after heavy alcohol consumption.
- Do not measure when you have a fever.

## **CAUTION**

- Your body composition may not be accurate after a shower or heavy physical activity.
- This smart scale model is not intended for use by infants or toddlers. Some body composition measurements are not for children under 18.
- Arboleaf smart scale should not be used for treating or diagnosing any medical conditions.
- Always consult your physician or health professional to determine what is suited for you.

## **WARNING**

Please read carefully before use

- Do not use if you have a pacemaker or any other internal medical devices in your body.

- Do not use if you are pregnant.
- Do not use on slippery floors.
- Do not use when your body is wet.
- Do not use rechargeable 1.2V AAA batteries.

## Package Contents

1. Smart Scale
2. 2×1.5V AAA batteries
3. User Guide

## Power Up Scale

Your scale shipment comes with 2 AAA batteries. Put the batteries in the battery compartment. Your scale should display 0.0. The scale is ready for use.

## Initial Setup

- **Step 1:** Download “The Arboleaf” App and create an account

From Apple App Store (iOS) or Google Play (Android), search for “Arboleaf” and download it to your smartphone. You can also scan the QR code to download.



Download “**Arboleaf**” App  
or scan the QR code




Register an account

- **Step 2:** Pair the scale with your smartphone. Make sure Bluetooth is enabled on your phone. For Android 6.0 and above, Location also needs to be enabled. Pairing is done automatically. The Bluetooth icon is displayed on the scale.



Open **Arboleaf** app




Apply pressure to turn on scale

## Take Measurements

Keep your app open on the “Measurements” screen.

- Apply pressure to turn on the scale, wait for “CAL” to return to “0.0”



- Step on the scale barefoot until the numbers are stable



- Check measurements





When the display number flashes twice, it indicates that the weight measurement has finished. When the number 0 stops moving back and forth, the body fat and other body composition measurements have been completed.

## FCC Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for

Compliance could void the user's authority to operate the equipment.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used according to the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.

- Consult the dealer or an experienced radio/TV technician for help.

## **Arboleaf Body Composition Smart Scale User Guide**

Arboleaf@ Body Composition Smart Scale measures weight, body fat percentage, BMI, and other body composition data using bioelectrical impedance analysis (BIA) technology. This smart scale is intended for home use only.

### **E – CrossStu – GmbH**

Mainzer Landstr. 69

60329 Frankfurt am Main

**Mail:** [E-CrossStu@web.de](mailto:E-CrossStu@web.de)

### **Arboleaf Corporation**

5465 Legacy Drive, Suite 650 Plano, TX 75024

1-[800-658-1148](tel:800-658-1148)

[support@arboleaf.com](mailto:support@arboleaf.com)

Custom Support Service (EU)

[support.eu@arboleaf.com](mailto:support.eu@arboleaf.com)

[www.arboleaf.com](http://www.arboleaf.com)

**Company:** YH Consulting Limited

**Address:** C/O YH Consulting Limited

Office 147, Centurion House, London

Road Staines-upon-Thames,

sSurrey TW18 4AX

**Telephones:** +44 07514-677868

**E-Mail:** [H2YHUK@gmail.com](mailto:H2YHUK@gmail.com)

**Manu acted as Shenzhen Yolan Technology Co., Ltd**


**Address:** Room 201, J Jinfula Building, No.49-1, Dabao Road, Dalang Area, Xinan


Street, Baoan, Shenzhen, China

[service@yolanda.hk](mailto:service@yolanda.hk)

## **Documents / Resources**



  
arboleaf  
www.arboleaf.com

  
arboleaf Body Composition Smart Scale  
User Guide

Arboleaf® Body Composition Smart Scale measures weight, body fat percentage, body water and other body composition data using bioelectrical impedance analysis (BIA) technology. This smart scale is designed for home use only.

**CS20N**

• Capacity: 300 lbs  
• Accuracy: ±0.5 lbs  
• Weight: 1.5 lbs  
• Dimensions: 11.5" x 11.5" x 1.5"

**CS20N**

• Capacity: 300 lbs  
• Accuracy: ±0.5 lbs  
• Weight: 1.5 lbs  
• Dimensions: 11.5" x 11.5" x 1.5"

Arboleaf Body Composition Smart Scale measures weight, body fat percentage, body water and other body composition data using bioelectrical impedance analysis (BIA) technology. This smart scale is designed for home use only.


Arboleaf Body Composition Smart Scale measures weight, body fat percentage, body water and other body composition data using bioelectrical impedance analysis (BIA) technology. This smart scale is designed for home use only.

[arboleaf CS20N Body Composition Smart Scale \[pdf\]](#) User Guide

CS20N, CS20N Body Composition Smart Scale, Body Composition Smart Scale, Smart Scale, Scale

## References

- [User Manual](#)

 arboleaf, Body Composition Smart Scale, CS20N, CS20N Body Composition Smart Scale, Scale, Smart Scale

## Leave a comment

Your email address will not be published. Required fields are marked \*

Comment \*

Name

Email

Website

☐ Save my name, email, and website in this browser for the next time I comment.

Post Comment

## Search:

e.g. whirlpool wrf535swhz

Search

[Manuals+](#) | [Upload](#) | [Deep Search](#) | [Privacy Policy](#) | [@manuals.plus](#) | [YouTube](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.