

Homedics PAR-200

HoMedics PAR200 ParaSpa Deluxe Paraffin Bath Instruction Manual

Model: PAR-200

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, basic safety precautions should always be followed, including the following:

- **READ ALL INSTRUCTIONS BEFORE USING.**
- **DANGER:** To reduce the risk of electrocution: Always unplug the appliance from the electrical outlet immediately after using and before cleaning. Do not reach for an appliance that has fallen into water. Unplug it immediately.
- **WARNING:** To reduce the risk of burns, fire, electric shock, or injury to persons: Close supervision is necessary when this appliance is used by, on, or near children, invalids, or disabled persons.
- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by HoMedics.
- Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return it to a HoMedics Service Center for examination and repair.
- Keep cord away from heated surfaces.
- Never use while sleeping or drowsy.
- Never drop or insert any object into any opening.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Do not use outdoors.
- Do not immerse the unit in water or any other liquid.
- This appliance is intended for household use only.
- Do not attempt to repair the appliance. There are no user-serviceable parts.

SAVE THESE INSTRUCTIONS.

PRODUCT FEATURES

- Paraffin heat treatment for softening and smoothing hands, elbows, and feet.
- Heat therapy helps soothe sore muscles and aching joints.
- Adjustable temperature control for safe and comfortable use.
- Cool-touch tub design, keeping the inner tub cooler to the touch.
- Variable control dial to select desired temperature.
- Power on/off LED indicator.
- Heater on/off indicator.
- Includes 6 pounds of pure, hypo-allergenic paraffin wax (no scents or dyes).
- Includes 60 plastic liners to extend heat treatment benefits.
- Attractive transparent lid.
- Removable grate.

SETUP INSTRUCTIONS

1. Place the ParaSpa unit on a flat, sturdy, heat-resistant surface.
2. Remove the lid and place it aside.
3. Ensure the removable grate is at the bottom of the tub.
4. Open the bags of paraffin wax and break the wax into smaller pieces if necessary. Place the wax into the tub. The unit requires a minimum of 4 lbs of wax to cover the heating elements. The unit can hold up to 6 lbs of wax.
- 5.



This image illustrates the process of adding the included paraffin wax into the main basin of the HoMedics ParaSpa Deluxe Paraffin Bath. The unit is white with a control dial visible on the front.

6. Plug the unit into a 120V AC electrical outlet. The red "Power" indicator light will illuminate.
7. Turn the control dial to the "Hot" setting. The green "Heater On" indicator light will illuminate.
8. Allow approximately 2-2.5 hours for the wax to melt completely. The green "Heater On" light will turn off once the wax reaches the maximum temperature.
9. Once the wax is melted, turn the control dial to the "Low" or "Med" setting for a comfortable temperature. Wait approximately 15-20 minutes for the wax to cool to a safe and comfortable

temperature for dipping.

10. Before dipping, test the wax temperature on a small area of your wrist to ensure it is comfortable and not too hot.

OPERATING INSTRUCTIONS

1. Wash the area to be treated (hands, elbows, or feet) thoroughly and dry completely.
2. Slowly immerse the area into the melted paraffin wax. Remove immediately. Wait a few seconds for the wax to solidify.
3. Repeat this dipping process 3-5 times, allowing each layer to harden slightly before re-dipping. This builds up a thick wax layer.
4. For extended heat therapy, place a plastic liner (provided) over the waxed area. This helps retain heat and moisture.
5. Allow the wax to remain on the skin for 10-20 minutes.
6. To remove the wax, peel it off the skin. The wax should come off easily. Discard the used wax. Do not reuse wax that has been removed from the skin.
7. After treatment, apply a moisturizer to the treated area if desired.

Note: Always ensure the wax is at a comfortable temperature before dipping. Never leave the unit unattended while operating.

MAINTENANCE AND CLEANING

Cleaning the Unit:

- Always unplug the unit from the electrical outlet and allow it to cool completely before cleaning.
- To clean the exterior, wipe with a soft, damp cloth. Do not use abrasive cleaners or harsh chemicals.
- Never immerse the unit in water or any other liquid.

Replacing Wax:

- The paraffin wax can be reused for personal use. However, it is recommended to replace the wax after approximately 20-30 uses, or if it appears dirty or contaminated.
- To remove old wax, ensure the unit is unplugged and cooled. The wax will solidify. You can then gently pry out the solid block of wax.
- Discard the old wax and follow the "Setup Instructions" to add new paraffin wax.

Storage:

- Store the unit in a cool, dry place.
- Do not wrap the power cord around the unit.

TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Unit does not turn on.	Not plugged in; power outlet malfunction.	Ensure unit is securely plugged into a working electrical outlet. Check circuit breaker.
Wax is not melting.	Temperature dial not set to "Hot"; insufficient time for melting.	Ensure dial is set to "Hot" and allow 2-2.5 hours for initial melting. Verify green "Heater On" light is illuminated.
Wax is too hot.	Temperature dial set too high; insufficient cooling time.	After wax melts, turn dial to "Low" or "Med" and wait 15-20 minutes for wax to cool to a comfortable temperature. Always test on wrist.
Wax is too cool.	Temperature dial set too low.	Adjust dial to a higher setting (e.g., "Med") and allow time for wax to reheat.

SPECIFICATIONS

- **Model:** PAR-200
- **Product Dimensions:** 14.25 x 7 x 7.5 inches (or 14"L x 8"W x 8.5"H)
- **Item Weight:** 11 Pounds
- **Material:** Plastic, Stainless Steel
- **Wax Capacity:** 4 to 6 pounds
- **Power:** 120V AC (standard household outlet)

LIMITED WARRANTY

HoMedics sells its products with the intent that they are free from defects in manufacture and workmanship for a period of one year from the date of original purchase, except as noted below. HoMedics warrants that its products will be free of defects in material and workmanship under normal use and service. This warranty extends only to consumers and does not extend to Retailers.

To obtain warranty service on your HoMedics product, contact a Consumer Relations Representative by telephone at 1-800-466-3342 for assistance. Please make sure to have the model number of the product available.

HoMedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate HoMedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism, or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of HoMedics.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY

OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of HoMedics.

This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

CONTACT AND SUPPORT

For service or support, please contact HoMedics Consumer Relations:

- **Phone:** 1-800-466-3342
- **Website:** www.homedics.com
- **Address:** HoMedics Service Center, 3000 Pontiac Trail, Commerce Township, MI 48390

When contacting support, please have your product model number (PAR-200) available.