

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [TANITA](#) /

› [Tanita 2001W-B Ultimate Body Fat Bathroom Scale User Manual](#)

TANITA 2001W-B

Tanita 2001W-B Ultimate Body Fat Bathroom Scale User Manual

Model: 2001W-B

INTRODUCTION

The Tanita 2001W-B Ultimate Body Fat Bathroom Scale is designed to provide both accurate weight measurements and body fat analysis. Utilizing Bioelectrical Impedance Analysis (BIA) technology, this scale offers a comprehensive view of your body composition. It features modes for adult and child measurement, making it suitable for various users. Please read this manual thoroughly before use to ensure proper operation and to maximize the benefits of your scale.



This image displays the Tanita 2001W-B Ultimate Body Fat Bathroom Scale. It is white with a central gray control panel featuring an LCD screen and two circular buttons (orange and blue). Four metallic electrodes are visible on the scale's surface, two rectangular at the front and two circular at the back, designed for body fat analysis.

SAFETY INFORMATION

- **Do not use if you have a pacemaker or other internal medical device.** The bioelectrical impedance analysis (BIA) feature sends a low-level electrical signal through the body, which may interfere with such devices.
- **Do not use if you are pregnant.** The BIA readings may not be accurate and the effects of the electrical current on a fetus are unknown.
- Always place the scale on a hard, flat, and stable surface. Using it on carpet or uneven flooring may result in inaccurate readings.
- Ensure your feet are clean and dry before stepping on the scale for body fat analysis.
- Do not use the scale in a slippery environment, such as a wet bathroom floor, to prevent falls.
- Keep the scale away from water and extreme temperatures.
- This scale is intended for household use only.

PACKAGE CONTENTS

Verify that all items are present in the package:

- Tanita 2001W-B Ultimate Body Fat Bathroom Scale
- AA Batteries (4 included)
- Instruction Manual (this document)

SETUP

1. Battery Installation

1. Locate the battery compartment on the underside of the scale.
2. Open the battery compartment cover.
3. Insert the four AA batteries, ensuring correct polarity (+/-) as indicated inside the compartment.
4. Close the battery compartment cover securely.

2. Scale Placement

Place the scale on a hard, flat, and stable surface. Avoid placing it on carpets, rugs, or uneven flooring, as this can affect measurement accuracy.

OPERATING INSTRUCTIONS

1. Weight Measurement Only

For a quick weight reading without body fat analysis:

1. Gently tap the scale with your foot to turn it on. The display will show '0.0'.
2. Step onto the scale with bare feet, standing still and centered.
3. Your weight will be displayed. The scale will automatically turn off after a few seconds.

2. Body Fat Analysis

To measure body fat percentage, you will need to input personal data. The scale features easy-to-program individual information buttons.

1. Gently tap the scale to turn it on.
2. Press the 'SET' button (orange button) to begin entering your personal data.
3. Use the arrow buttons (blue button) to adjust values for parameters such as height, age, and gender. Press 'SET' to confirm each entry and move to the next.
4. Once all data is entered, the display will show '0.0'.
5. Step onto the scale with bare feet, ensuring your feet are positioned on the metallic electrodes. Stand still and centered.
6. The scale will first display your weight, then your body fat percentage (displayed in 1-percent increments).
7. The scale will automatically turn off after a few seconds. Note that personal data must be re-entered for each body fat reading as this model does not store user profiles.

3. Adult and Child Mode

The Tanita 2001W-B includes distinct modes for adult and child measurements to ensure appropriate body fat analysis. The specific method for switching between these modes is typically integrated into the personal data setup process. During the initial setup or when re-entering your data for body fat analysis, you will be prompted to select either 'Adult' or 'Child' based on the user's age. Ensure the correct mode is selected for accurate readings.

MAINTENANCE

Cleaning

- Wipe the scale surface with a soft, damp cloth. Do not use abrasive cleaners or immerse the scale in water.
- Ensure the metallic electrodes are clean and free of debris for accurate body fat readings.

Storage

- Store the scale in a cool, dry place, away from direct sunlight and extreme temperatures.
- Avoid placing heavy objects on the scale when not in use.

TROUBLESHOOTING

- **Inaccurate Weight Readings:** Ensure the scale is on a hard, flat surface. Avoid moving the scale frequently. Recalibrate by stepping on and off the scale once before taking a measurement.
- **No Display/Low Battery Indicator:** Replace all four AA batteries. Ensure batteries are inserted with correct polarity.
- **Body Fat Reading Errors:** Ensure feet are bare, clean, and dry, and positioned correctly on the electrodes. Verify that personal data (height, age, gender) has been entered correctly.
- **Scale Does Not Turn On:** Check battery installation. If batteries are new and correctly installed, contact customer support.

SPECIFICATIONS

Model Number	2001W-B
Weight Capacity Maximum	300 Pounds
Power Source	4 AA batteries (included)
Display	1-inch LCD screen
Body Fat Analysis Increments	1-percent
Color	White
Manufacturer	Tanita
UPC	742496200143

WARRANTY INFORMATION

The Tanita 2001W-B Ultimate Body Fat Bathroom Scale comes with a **1-year limited warranty** from the date of purchase. This warranty covers defects in materials and workmanship under normal use. It does not cover damage caused by misuse, accidents, unauthorized modifications, or failure to follow the instructions in this manual. Please retain your proof of purchase for warranty claims.

SUPPORT

For further assistance, technical support, or warranty inquiries, please contact Tanita customer service. Refer to the official Tanita website or product packaging for the most current contact information.