

[manuals.plus](#) /

› [Waterpik](#) /

› [Waterpik PS-100E Plaque Seeker Replacement Tips Instruction Manual](#)

## Waterpik PS-100E

# Waterpik PS-100E Plaque Seeker Replacement Tips

## INSTRUCTION MANUAL

### Introduction

---

Thank you for choosing Waterpik Plaque Seeker Replacement Tips. These tips are specifically designed to enhance your oral hygiene, particularly for individuals with dental implants, crowns, bridges, or permanent retainers. This manual provides detailed instructions for proper use, maintenance, and care to ensure optimal performance and longevity of your tips.

### Product Overview

---

The Waterpik Plaque Seeker Tip (PS-100E) features a unique design with three thin tufts of bristles. This specialized configuration is engineered to effectively access and remove stubborn plaque from areas that are prone to accumulation, such as around dental implants, crowns, bridges, and other restorations. Each package contains two color-coded tips, allowing multiple users to easily identify their personal tip for hygienic use.



Image: Two Waterpik Plaque Seeker Tips (PS-100E) in their packaging, illustrating the product as received.



Image: A close-up view of the Plaque Seeker Tip, highlighting the three thin tufts of bristles designed for targeted plaque removal.

## Compatibility

---

These Plaque Seeker Tips are compatible with most Waterpik Water Flosser models. For optimal performance, ensure your Waterpik Water Flosser is not among the incompatible models listed below.

### Compatible Models:

For use with any Waterpik water flosser **except** the models listed under incompatible models.

### Incompatible Models:

- Waterpik Classic (WP 60/65/72)
- Waterpik Whitening (WF-05, WF-06)
- Waterpik Sidekick (WF-04)
- Waterpik Sonic Fusion (SF-01, SF-02, SF-03, SF-04)

## Installation and Setup

---

1. **Select a Tip:** Choose a Plaque Seeker Tip. If multiple users share a device, use the color-coded

base to identify your personal tip.

2. **Insert the Tip:** Firmly press the tip into the handle of your Waterpik Water Flosser until it clicks into place. Ensure it is securely seated.
3. **Remove the Tip:** To remove, press the tip eject button (if available on your model) and pull the tip from the handle.

## Operating Instructions

---

Using the Plaque Seeker Tip is similar to using other Waterpik tips, with a focus on targeting specific areas.

1. **Fill Reservoir:** Fill the Waterpik Water Flosser reservoir with lukewarm water.
2. **Adjust Pressure:** Set the pressure control to the lowest setting for first-time use, then gradually increase to a comfortable level.
3. **Position the Tip:** Place the Plaque Seeker Tip with the bristles gently touching the tooth surface and the gumline, specifically targeting areas around implants, crowns, or bridges.
4. **Start Flossing:** Turn on the Water Flosser. Allow the water to flow through the bristles, gently massaging the gumline and cleaning around the dental work.
5. **Glide Along:** Move the tip slowly along the gumline, pausing briefly between teeth and around any restorations to allow the water and bristles to clean effectively.
6. **Maintain Posture:** Lean over the sink and allow water to flow from your mouth into the sink to avoid splashing.

Your browser does not support the video tag.

Video: A short demonstration by Water Pik Inc. showing the Plaque Seeker Tip in action, illustrating how it cleans around dental structures.



Image: An animated illustration depicting the Waterpik Plaque Seeker Tip effectively cleaning around a dental implant and crown, demonstrating its targeted action.

## Maintenance

---

Proper maintenance ensures the effectiveness and hygiene of your Plaque Seeker Tips.

- **Daily Cleaning:** After each use, rinse the tip thoroughly under running water.
- **Deep Cleaning:** Periodically, you may soak the tip in a solution of 2 parts water to 1 part white vinegar for 5-7 minutes. Rinse thoroughly afterward.
- **Storage:** Store tips in a clean, dry place.

## Troubleshooting

---

If you encounter any issues with your Plaque Seeker Tips, consider the following:

- **Tip Not Seating Properly:** Ensure the tip is fully inserted into the handle until it clicks. Check for any debris obstructing the connection.
- **Reduced Water Flow:** This may indicate a clog. Try soaking the tip in vinegar solution as described in the Maintenance section. If the issue persists, consult your Waterpik Water Flosser manual for device-specific troubleshooting.

- **Bristles Wearing Out:** Bristle wear is normal over time. Refer to the Replacement Schedule section for recommended replacement frequency.

## Specifications

<b>Model Number</b>	PS-100E
<b>Unit Count</b>	2 Count
<b>Material Type</b>	Plastic
<b>Product Benefits</b>	Gum Health, Removes Plaque
<b>Item Dimensions</b>	3.13 x 0.75 x 6.75 inches
<b>Item Weight</b>	0.02 Grams

## Replacement Schedule

For optimal hygiene and performance, it is recommended to replace your Waterpik Plaque Seeker Tips every **3 months**. Regular replacement ensures the bristles remain effective at removing plaque and maintaining gum health.

## Warranty and Support

Waterpik stands by the quality of its products. For any product questions or needs, Waterpik's dedicated support team is available to assist you. Please refer to the official Waterpik website or your Waterpik Water Flosser device manual for specific warranty information and contact details.

**Customer Support:** Waterpik's Colorado support team is available to help with any product questions or needs.