

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [NuFACE](#) /

› [NuFACE Trinity Microcurrent Facial Toning Device Instruction Manual](#)

## NuFACE Trinity

# NuFACE Trinity Microcurrent Facial Toning Device Instruction Manual

MODEL: TRINITY

## Product Overview

NuFACE Trinity is a handheld microcurrent facial toning device designed to improve facial contour, skin tone, and reduce the appearance of fine lines and wrinkles. It is clinically tested and delivers professional-level results. The included Hydrating Leave-On Gel Primer is formulated without alcohol, fragrance, parabens, or oil.

## Key Features

- **NuFACE Trinity Starter Kit:** Enjoy a 5-Minute Facial-Lift with this pro-level microcurrent facial toning device.
- **Microcurrent Technology:** Delivers low-level electrical current with gentle waves to facial muscles, mimicking the body's natural current to tone, lift, and contour.
- **Targeted Treatment:** Stimulates face and neck to tighten skin and improve the look of fine lines and wrinkles. Recommended for 5 minutes a day, 5 days a week.
- **Hydrating Leave-On Gel Primer:** Essential for device use, ensuring proper microcurrent delivery.
- **Clinically Proven:** Products are clinically shown to provide instant results and cumulative benefits in treated areas.

## Important Safety Information

**The NuFACE Trinity, TRINITY+, Mini, MINI+ and FIX facial trainer devices are contraindicated against and should not be used if you have epilepsy or seizures, or if you have an electronic implanted device such as a defibrillator, pacemaker or ECG monitor. The devices should also not be used if you are pregnant, have active cancer or are under the age of 18. If you have any medical concerns, please consult your doctor before using a device. We recommend waiting to use your device 2 weeks after receiving a hydrafacial, botox or fillers, and waiting 1 week after a microneedling treatment.**

## Initial Setup and Charging

---

Before first use, fully charge your NuFACE Trinity device. Connect the charging cable to the device and a power source. The indicator lights will show charging progress.

## Preparing for Treatment

---

Cleanse and dry your face and neck thoroughly. Apply a mask-like layer of the Hydrating Leave-On Gel Primer to the treatment area. Ensure the skin is well-covered with primer for optimal microcurrent conductivity.

## Using Your NuFACE Trinity Device

---

The NuFACE Trinity device is designed for a 5-minute facial-lift.

### Step-by-Step Guide:

- 1. Power On:** Press the power button to turn on the device. Select your desired intensity level (the device typically starts at the lowest level).
- 2. Glide and Hold:** Perform glides and holds using the device across targeted zones of your face and neck. Follow the recommended treatment areas and techniques for best results.
- 3. Automatic Shut-off:** The device is programmed to automatically shut off after 5 minutes of use. You may also manually turn it off.

### Recommended Usage:

- Use for 5 minutes a day, 5 days a week.
- For detailed guidance on specific glides and holds for different facial areas, refer to the official NuFACE app or website.

## Care and Cleaning

---

After each use, wipe the device clean with a damp, soft cloth. Do not submerge the device in water. Store in a cool, dry place.

## Troubleshooting Tips

---

- **No Sensation:** Ensure sufficient gel primer is applied to the skin. The device requires proper conductivity.
- **Device Not Turning On:** Check if the device is fully charged. Recharge if necessary.
- **Irritation:** If irritation occurs, reduce the intensity level or discontinue use and consult a healthcare professional.

## Product Specifications

---

<b>Brand</b>	NuFACE
<b>Model Name</b>	Trinity

<b>Item Model Number</b>	Trinity Starter Kit
<b>Product Dimensions</b>	9.2 x 5.6 x 5.5 inches
<b>Item Weight</b>	1 Pounds
<b>Batteries</b>	1 Lithium Metal batteries required.
<b>UPC</b>	685239988888, 858271002250
<b>Manufacturer</b>	Carol Cole Company

## Legal Disclaimer

This product is intended for cosmetic use only. It is not intended to diagnose, treat, cure or prevent any medical condition or disease. Please read all instructions and warnings before use. Use only as directed. Always consult with a healthcare professional if you have any medical concerns or conditions before using the device.

## Additional Information

**Verified by Transparency:** Each item has a unique code that we verify before shipping.

## Product Visuals and Demonstrations

### Product Images



NuFACE Trinity Facial Toning Device in use on a woman's face.

## NuFACE Trinity® Facial Toning Device



Improve **facial contour, tone** and appearance of **fine lines and wrinkles**

- **5-Minute Facial-Lift**

- **5 minutes** a day, **5 days** a week

- **Interchangeable attachments** for personalized routine

- **FDA-Cleared, Aesthetician Approved**

NuFACE Trinity Facial Toning Device with Hydrating Leave-On Gel Primer.

## NuFACE Trinity® FACIAL TONING DEVICE



**85%** of users experienced improved facial contour

Unretouched. Clinical study, a single center, single arm, non-randomized trial to evaluate the efficacy of the NuFACE Trinity device with Facial Trainer Attachment.

Before and after comparison of facial contour after 60 days of NuFACE Trinity use.



NuFACE Trinity device with various interchangeable attachments for targeted treatments.

	 <b>Trinity®</b> <i>Pro-Level Targeted Treatments</i>	 <b>Mini</b> <i>On-The-Go Treatments</i>	 <b>FIX®</b> <i>Instant, Quick Fix</i>
<b>WHAT IT DOES</b>	<ul style="list-style-type: none"> <li>Tones, lifts, and contours face and neck</li> <li>Reduces appearance of fine lines and wrinkles</li> <li>Delivers the microcurrent down to the muscle</li> </ul>		<ul style="list-style-type: none"> <li>Helps smoothes lines around eyes, mouth, and forehead</li> <li>Fuller-looking lips</li> <li>Delivers the microcurrent along surface layer of the skin</li> </ul>
<b>TREATMENT LENGTHS</b>	5-15 minutes		3 minutes
<b>TREATMENT AREAS</b>	Neck & Jowls ✓ Cheeks ✓ Forehead ✓ Eyes & Lips With attachment ✓ Nasolabial folds With attachment ✓	✓ ✓ ✓	✓ ✓
<b>ATTACHMENTS</b>	<ul style="list-style-type: none"> <li>Effective Lip &amp; Eye (ELE)</li> <li>Trinity Wrinkle Reducer (TWR)</li> </ul> <small>*Attachments sold separately.</small>		

Comparison chart detailing features, treatment lengths, and areas for NuFACE Trinity, Mini, and FIX devices.



NuFACE Trinity device being used on a woman's cheek to tone and lift.



Packaging of the NuFACE Trinity Advanced Facial Toning Device.

## Official Product Videos

Your browser does not support the video tag.

An official video demonstrating the components and basic usage of the NuFACE Trinity Starter Kit.