

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [DonJoy](#) /

› [DonJoy IceMan Universal Shoulder Wrap-On Pad for Cold Therapy Instruction Manual](#)

DonJoy 11-1334-9-00000

DonJoy IceMan Universal Shoulder Wrap-On Pad Instruction Manual

Model: 11-1334-9-00000

INTRODUCTION

The DonJoy IceMan Universal Shoulder Wrap-On Pad is designed to provide effective cold therapy for the shoulder and other body areas. This pad is an accessory intended for use with DonJoy IceMan cold therapy units, specifically the DonJoy IceMan Clear3 and DonJoy IceMan Classic models, which are sold separately. It facilitates continuous cold water circulation to help manage pain and swelling.

PRODUCT OVERVIEW



Pad Design and Features

The Universal Shoulder Wrap-On Pad features a large surface area, ensuring comprehensive cooling coverage. Its design allows for efficient water flow, reaching the very edges of the pad for consistent cold therapy. The pad is constructed from Polyurethane (PU) and measures approximately 14 x 23 x 12.5 inches, weighing about 12.63 ounces.

One side of the pad is designed to be very cold for direct application, while the other side is well-insulated to prevent excessive cold transfer to unintended areas or for user comfort. The pad is versatile and can be adapted for use on other body areas, such as the knee, in addition to the shoulder.

Compatibility

This pad is compatible with all DonJoy IceMan units, including the IceMan Clear3 and IceMan Classic models. Please note that the cold therapy unit itself is sold separately.

SETUP INSTRUCTIONS

- 1. Prepare the IceMan Unit:** Ensure your DonJoy IceMan Clear3 or Classic cold therapy unit is clean and ready for use.
- 2. Fill the Unit:** Fill the IceMan unit with ice and water according to the instructions provided with your specific IceMan model.
- 3. Connect the Pad:** Connect the tubing from the Universal Shoulder Wrap-On Pad to the appropriate ports on your IceMan unit. Ensure connections are secure to prevent leaks.
- 4. Prime the System:** Follow your IceMan unit's instructions to prime the system, allowing water to circulate through the pad.

OPERATING INSTRUCTIONS

- 1. Position the Pad:** Place the cold therapy pad directly onto the desired area (e.g., shoulder, knee). Ensure the side designed for cold application is against the skin or a thin protective layer as advised by a medical professional.
- 2. Secure the Pad:** Use the integrated straps or wraps to secure the pad firmly but comfortably in place. The pad is designed to wrap around the shoulder completely, and its versatility allows for secure placement on other areas like the knee.
- 3. Initiate Therapy:** Turn on your DonJoy IceMan unit to begin the circulation of cold water through the pad.
- 4. Monitor Comfort:** Periodically check the skin under the pad for any signs of discomfort or excessive cold. Adjust placement or duration as needed. The pad delivers consistent coldness effectively.
- 5. Duration of Use:** Follow your medical professional's recommendations for the duration and frequency of cold therapy sessions.

MAINTENANCE

- Cleaning:** The pad and tubing are made of excellent quality materials that are easy to clean. After each use, drain any remaining water from the pad and tubing. Clean the pad with a mild soap and water solution, then rinse thoroughly. Allow to air dry completely before storage.
- Storage:** Store the pad in a clean, dry place away from direct sunlight and extreme temperatures.
- Single Patient Use:** For safety and hygiene reasons, this cold therapy pad is intended for use by a single patient only. After its first use, water remains in the pad. It should not be returned or transferred to another person.

TROUBLESHOOTING

No specific troubleshooting steps are provided for the pad itself. For issues related to the cold therapy unit (e.g., water circulation, leaks from the unit), please refer to the instruction manual for your specific DonJoy IceMan unit.

If the pad is not delivering sufficient coldness, ensure the IceMan unit is properly filled with ice and water and that the system is primed correctly. Check for any kinks in the tubing that might restrict water flow.

SPECIFICATIONS

Product Name	DonJoy IceMan Universal Shoulder Wrap-On Pad for Cold Therapy
Model Number	11-1334-9-00000
Brand	DonJoy
Material	Polyurethane (PU)

Product Dimensions (L x W x H)	14 x 23 x 12.5 inches
Item Weight	12.63 ounces
Manufacturer	DonJoy Braces
Date First Available	August 29, 2014

IMPORTANT INFORMATION

Safety Information

For use by a single patient. After first use, water remains in the cold therapy pad. For safety reasons, the cold pad should not be returned or transferred to another person.

Indications

This product is indicated for cold therapy applications, including but not limited to: Knee Injury, Knee Sprain, Knee Inflammation, Knee Surgery, Post-Op Knee Recovery, Knee Replacement, Knee Ligament Injury. While specifically designed for shoulder, its universal design allows for effective use on knees as well.

Legal Disclaimer

Item may be returned only if unused due to hygiene. If unused, then normal return policy applies.

WARRANTY AND SUPPORT

Specific warranty details for the DonJoy IceMan Universal Shoulder Wrap-On Pad are not provided in this manual. Please refer to the official DonJoy website or contact DonJoy customer support for detailed warranty information and product support. For further assistance or inquiries, please visit the [DonJoy Store on Amazon](#) or consult your medical professional.