

## Best Choice Products SKY2095

# Shiatsu Foot Massager Instruction Manual

Brand: **Best Choice Products** | Model: **SKY2095**

## INTRODUCTION

---

Thank you for choosing the Best Choice Products Shiatsu Foot Massager. This device is designed to provide therapeutic relief for your feet, targeting pain, plantar fasciitis, neuropathy, and improving blood circulation. With its compact design and remote control functionality, it offers a convenient and effective way to relax and rejuvenate your feet. Please read this manual thoroughly before use to ensure safe and optimal operation.



Figure 1: Best Choice Products Shiatsu Foot Massager with remote control.

## SAFETY INFORMATION

---

- Avoid using the massager while sleeping or with wet feet.
- Always follow electrical safety guidelines. Ensure the power cord is not damaged and is properly plugged into a suitable outlet.
- Be cautious when using high intensity rollers and heating functionality to prevent discomfort or injury. Start with lower intensity settings and gradually increase as comfortable.
- Always follow the instructions provided with the product. Do not attempt to disassemble or repair the unit yourself.
- This product is intended for adult use.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. Consult a healthcare professional if you have pre-existing medical conditions or concerns.

## SETUP

---

1. **Unpack:** Carefully remove the foot massager and all accessories from its packaging.
2. **Placement:** Place the massager on a stable, flat surface. Ensure there is enough space around the unit for proper ventilation.
3. **Power Connection:** Plug the power cord into a standard electrical outlet. The unit is now ready for use.
4. **Removable Sleeves:** The massager features removable foot sleeves for hygiene. Ensure they are properly zipped into place before use.

## ADDITIONAL FEATURES

This foot massager is designed to help provide both physical pain relief and peace of mind in the form of convenience

REMOVABLE &  
WASHABLE SLEEVE

BRIGHT  
LCD SCREEN

NON-SLIP  
BOTTOM

ANTI-SKID  
PAD



Figure 2: Removable and washable foot sleeves for easy cleaning.

## OPERATING INSTRUCTIONS

The Best Choice Products Shiatsu Foot Massager offers various settings for a personalized experience. You can control the unit directly via its LCD touch panel or using the included remote control.

### Remote Control & LCD Panel Functionality:

- **Power (ON/OFF):** Press the power button to turn the massager on or off.
- **Auto Mode:** Select from 3 different automatic modes that target the toe, arch, and sole of your feet for a deep, therapeutic massage.
- **Custom/Manual Mode:** Allows you to manually select specific massage techniques and areas.
- **Speed:** Adjust the intensity of the massage rollers.
- **Direction:** Change the direction of the massage rollers.
- **Time:** Set the duration of your massage session.

# SCREEN AND REMOTE FUNCTIONALITY

Operate the massager manually on its LCD touch panel, or use the handy remote control



Figure 3: LCD screen and remote control for easy operation.

## Massage Techniques:

The massager utilizes a variety of techniques to provide comprehensive relief:

- **Pulsing:** Gentle, rhythmic pressure to stimulate circulation.
- **Kneading:** Deep, circular motions to relieve muscle tension.
- **Rolling:** Rollers target specific pressure points on the soles of your feet.





Figure 4: Visual representation of various massage techniques.

## MAINTENANCE

---

Regular maintenance ensures the longevity and hygiene of your foot massager.

- **Cleaning Sleeves:** The foot sleeves are removable and designed for easy cleaning. Unzip them from the unit and wash according to fabric care instructions. Allow them to dry completely before reattaching.
- **Exterior Cleaning:** Wipe the exterior of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water.
- **Storage:** When not in use, store the massager in a cool, dry place away from direct sunlight. Its compact structure allows for effortless storage.

## ADDITIONAL FEATURES

This foot massager is designed to help provide both physical pain relief and peace of mind in the form of convenience

REMOVABLE &  
WASHABLE SLEEVE

BRIGHT  
LCD SCREEN

NON-SLIP  
BOTTOM

ANTI-SKID  
PAD



Figure 5: The removable sleeves can be unzipped for cleaning.

## TROUBLESHOOTING

Problem	Possible Cause / Solution
Unit does not power on.	<ul style="list-style-type: none"><li>• Ensure the power cord is securely plugged into both the massager and a working electrical outlet.</li><li>• Check if the power outlet is functional by plugging in another device.</li></ul>
Remote control is not working.	<ul style="list-style-type: none"><li>• Check and replace the remote control batteries if necessary.</li><li>• Ensure there are no obstructions between the remote and the massager's sensor.</li></ul>
Massage feels too intense.	<ul style="list-style-type: none"><li>• Reduce the speed/intensity setting using the control panel or remote.</li><li>• Wear thicker socks to cushion the massage.</li><li>• Adjust your foot position within the massager to find a more comfortable pressure point.</li></ul>
Unusual noise during operation.	<ul style="list-style-type: none"><li>• Ensure no foreign objects are inside the foot chambers.</li><li>• If the noise persists, discontinue use and contact customer support.</li></ul>

## SPECIFICATIONS

- **Product Dimensions:** 23 x 13 x 10 inches
- **Item Weight:** 13.5 Pounds
- **Material:** Plastic
- **Power Source:** Plug-in
- **Model Number:** SKY2095
- **Manufacturer:** Best Choice Products
- **First Available Date:** October 23, 2015

## PRODUCT DIMENSIONS



Figure 6: Product dimensions for easy reference.

## WARRANTY & SUPPORT

For any questions, concerns, or support regarding your Best Choice Products Shiatsu Foot Massager, please refer to the official user guide or contact Best Choice Products customer service. The official Instructions for Use (IFU) document can be found [here](#).

