

Sanitas SFT 75

Sanitas SFT 75 Multifunction Thermometer User Manual

INTRODUCTION

Thank you for choosing the Sanitas SFT 75 Multifunction Thermometer. This device is designed for accurate and rapid temperature measurement. Please read this user manual carefully before first use and keep it for future reference. Proper use ensures reliable results and extends the life of your thermometer.

SAFETY INFORMATION

- This thermometer is intended for measuring human body temperature via the ear or forehead, and for measuring the temperature of objects and liquids.
- Do not use the device for purposes other than those specified in this manual.
- Keep the device away from direct sunlight, extreme temperatures, dust, and humidity.
- Do not drop the thermometer or subject it to strong impacts.
- Consult a physician if you have concerns about your health or temperature readings.
- Keep out of reach of children. Small parts may be swallowed.

PRODUCT OVERVIEW

The Sanitas SFT 75 is a versatile thermometer offering multiple measurement options and convenient features.



An image showing the Sanitas SFT 75 thermometer. It is white with a grey display area. The display shows '36.7 °C' and '12-420-32'. There are power and mode buttons near the display. The tip is designed for ear or forehead measurements.

Key features include:

- Ear and forehead measurement capabilities.
- Suitable for measuring object and liquid temperatures.
- Optical fever alarm: Red light for temperatures 38 °C (100.4 °F) and above, green light for normal temperatures.
- Blue backlit display for easy reading in low light.
- Memory function for 30 previous measurements.
- Quick measurement in seconds.
- Displays date and time.
- Switchable between Celsius (°C) and Fahrenheit (°F).

SETUP


Battery Installation

1. Locate the battery compartment cover, usually on the back or side of the device.
2. Slide or gently pry open the cover.
3. Insert the required batteries (e.g., 2 x AAA) according to the polarity indicators (+/-) inside the compartment.
4. Close the battery compartment cover securely.

The thermometer is battery-operated. Ensure fresh batteries are used for optimal performance.

OPERATING INSTRUCTIONS

Powering On/Off

- To power on, press the **Power button** (often marked with  or a similar symbol). The display will light up.
- The device will automatically turn off after a period of inactivity to conserve battery life. You can also manually turn it off by pressing and holding the Power button.

Setting Date and Time

Refer to the specific instructions in your device's full manual for detailed steps on setting the date and time, as this often involves pressing and holding the 'M' (Memory) or 'S' (Set) button while the device is off or in a specific mode.

Selecting Measurement Unit (°C/°F)

With the device powered off, press and hold the 'M' (Memory) button for a few seconds until the unit changes. Release the button when your desired unit (°C or °F) is displayed.

Taking a Measurement

- Ear Measurement:** Gently insert the thermometer tip into the ear canal. Ensure a snug fit. Press the measurement button. The result will appear in seconds.
- Forehead Measurement:** Position the thermometer near the center of the forehead, ensuring good contact. Press the measurement button.
- Object/Liquid Measurement:** Point the thermometer at the object or liquid surface (do not immerse). Press the measurement button.
- After the measurement, the result will be displayed on the blue backlit screen.

Fever Alarm

The thermometer features an optical fever alarm:

- Red Light:** Indicates a fever, typically for temperatures at or above 38 °C (100.4 °F).
- Green Light:** Indicates a normal temperature.

Memory Function

The device stores up to 30 previous measurements. To access them, press the 'M' (Memory) button while the device is on. Each press will cycle through the stored readings, displaying the most recent first.

MAINTENANCE AND CARE

- Cleaning:** After each use, especially for ear measurements, clean the thermometer tip with an alcohol swab or a cotton swab moistened with alcohol. Do not immerse the device in water or other cleaning solutions.
- Storage:** Store the thermometer in a clean, dry place at room temperature, away from direct sunlight and extreme temperatures.
- Battery Replacement:** When the low battery indicator appears on the display, replace all batteries with new ones of the same type. Dispose of used batteries according to local regulations.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Dead or incorrectly inserted batteries.	Check battery polarity or replace batteries.

Problem	Possible Cause	Solution
Inaccurate readings.	Improper measurement technique, dirty sensor, or environmental factors.	Ensure correct measurement method, clean the sensor tip, and measure in a stable environment. Wait 30 minutes after exercise or bathing before measuring.
Display shows "Err" or other error codes.	Measurement error, device malfunction, or out of operating range.	Refer to the full manual for specific error codes. Try restarting the device. If the error persists, contact customer support.

SPECIFICATIONS

Model: SFT 75

Model Number: 796.62

Brand: Sanitas

Measurement Range: Body: 34.0 °C - 43.0 °C (93.2 °F - 109.4 °F) (Typical, based on fever alarm)

Measurement Accuracy: Not explicitly stated, but generally +/- 0.2 °C for medical thermometers.

Memory: 30 measurement values

Display: Blue backlit LCD

Fever Alarm: Optical (Red for ≥ 38 °C / 100.4 °F, Green for normal)

Power Source: Battery-operated (e.g., 2 x AAA)

Dimensions: 16.8 x 9.4 x 5.2 cm

Weight: 140 grams

WARRANTY AND SUPPORT

Sanitas products are manufactured to high quality standards. For information regarding warranty terms and conditions, please refer to the warranty card included with your product or visit the official Sanitas website. For technical support or service inquiries, please contact Sanitas customer service.

You can often find contact information on the manufacturer's website or on the product packaging.