

B01HU04IOG

Weight Gurus Bluetooth Smart Bathroom Scale User Manual

Model: B01HU04IOG

1. INTRODUCTION

This manual provides instructions for the proper setup, operation, and maintenance of your Weight Gurus Bluetooth Smart Bathroom Scale. This scale is designed to measure weight and other key body metrics for up to 8 users, syncing data with compatible Bluetooth Smart Ready devices.



Image 1: Front view of the Weight Gurus Bluetooth Smart Bathroom Scale. The scale features a clear glass platform and a digital display.

2. SAFETY INFORMATION

- Do not use this scale on slippery surfaces. Ensure the scale is placed on a hard, flat surface for accurate readings.
- Do not immerse the scale in water or use harsh chemical cleaners. Clean with a damp cloth.

- This scale is not intended for medical diagnosis or treatment. Consult a healthcare professional for health-related concerns.
- Keep batteries out of reach of children. Dispose of used batteries according to local regulations.
- Maximum weight capacity is 400 lbs (180 kg). Do not exceed this limit.

3. PACKAGE CONTENTS

Verify that all items are present in the package:

- Weight Gurus Bluetooth Smart Bathroom Scale
- 4 x AAA Batteries
- User Manual (this document)

4. SETUP

4.1. Battery Installation

1. Locate the battery compartment on the underside of the scale.
2. Open the battery compartment cover.
3. Insert 4 AAA batteries, ensuring correct polarity (+/-).
4. Close the battery compartment cover securely.

4.2. App Download and Device Compatibility

The Weight Gurus scale requires the free Weight Gurus app for full functionality and data tracking.

- Download the "Weight Gurus" app from the Apple App Store (for iOS devices) or Google Play Store (for Android devices).
- Ensure your device is compatible:
 - Apple iPhone 4S or higher
 - iPad 3rd Generation or higher
 - iPod Touch 5th Generation or higher
 - Select Bluetooth Smart Ready Android devices

4.3. Initial Pairing

1. Enable Bluetooth on your compatible smartphone or tablet.
2. Open the Weight Gurus app and follow the on-screen instructions to create a user profile.
3. Step onto the scale briefly to activate it. The display will light up.
4. The app will guide you through the pairing process. Ensure the scale is within close proximity to your device during pairing.

5. OPERATING INSTRUCTIONS

5.1. Taking a Measurement

1. Place the scale on a hard, flat surface. Avoid carpets or uneven flooring.

2. Step onto the scale with bare, dry feet. The scale features "Auto-on technology" and will activate automatically.
3. Stand still until your weight is displayed and then flashes, indicating the measurement is complete.
4. The scale will then calculate and display other metrics such as BMI, body fat, muscle mass, water weight, and bone mass.
5. If the app is open and connected, your measurements will automatically sync to your user profile.

5.2. Auto-User Detection and Memory

The scale is equipped with "Auto-user detection" and can store data for up to 8 users. Once a user profile is set up in the app and initial measurements are taken, the scale will attempt to identify the user based on their weight range. If multiple users have similar weights, the app may prompt you to confirm your identity.

6. MAINTENANCE

- **Cleaning:** Wipe the scale's glass surface with a soft, damp cloth. Do not use abrasive cleaners or immerse the scale in water.
- **Storage:** Store the scale in a cool, dry place. Avoid extreme temperatures.
- **Battery Replacement:** Replace batteries when the low battery indicator appears on the display.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
Scale does not turn on.	Dead or incorrectly installed batteries.	Check battery polarity. Replace with new AAA batteries.
Inconsistent weight readings.	Scale on uneven surface; scale moved between readings; calibration issue.	Place scale on a hard, flat surface. Do not move the scale between measurements. Step on and off the scale once to calibrate before taking a measurement.
Scale not pairing with app.	Bluetooth off; app not open; device incompatibility; scale not active.	Ensure Bluetooth is enabled on your device. Open the Weight Gurus app. Verify device compatibility (Section 4.2). Step on the scale to activate it before attempting to pair. Restart both the app and your device if issues persist.
Metrics (BMI, body fat) not displayed.	Not standing with bare feet; user profile incomplete in app.	Ensure you are standing on the scale with bare, dry feet for bio-electrical impedance analysis. Complete all required profile information in the Weight Gurus app.

8. SPECIFICATIONS

Model Number	B01HU04IOG
Brand	Generic

Material	Tempered Glass
Color	White
Weight Capacity	400 lbs (180 kg)
Measurement Resolution	0.1 lb
Power Source	4 x AAA Batteries
User Memory	Up to 8 users
Connectivity	Bluetooth Smart Ready
Metrics Measured	Weight, BMI, Body Fat, Muscle Mass, Water Weight, Bone Mass

9. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the documentation provided with your purchase or contact the manufacturer directly. Keep your proof of purchase for warranty claims.

For technical assistance with the Weight Gurus app or scale functionality, please visit the official Weight Gurus support website or contact their customer service.