

Leaone BF8028-B

Leaone BF8028-B Smart Bluetooth Body Weight Scale User Manual

Model: BF8028-B

INTRODUCTION

The Leaone BF8028-B Smart Bluetooth Body Weight Scale is designed to provide comprehensive body composition analysis. This digital scale utilizes advanced sensors to measure not only your weight but also 11 additional body data points, including body fat, water percentage, muscle mass, bone mass, and more. Data is seamlessly synced to a dedicated smartphone application for iOS and Android devices, allowing you to track your health and fitness progress over time.

This manual provides instructions for setting up, operating, and maintaining your Leaone Smart Bluetooth Body Weight Scale, along with troubleshooting tips and product specifications.

SETUP GUIDE

Follow these steps to set up your Leaone BF8028-B Smart Bluetooth Body Weight Scale:

1. **Unpack the Scale:** Carefully remove the scale from its packaging. Ensure all components, including any included batteries, are present.
2. **Install Batteries:** Open the battery compartment on the underside of the scale. Insert the required batteries, ensuring correct polarity.
3. **Place the Scale:** Position the scale on a hard, flat surface. Avoid carpets or uneven flooring, as this can affect measurement accuracy.
4. **Download the App:** Scan the QR code provided in the product's instruction booklet with your smartphone to download the official Leaone health tracking application. Alternatively, search for "Leaone" in your device's app store (iOS App Store or Google Play Store).
5. **Create an Account:** Open the app and follow the on-screen prompts to create a user account.

6. **Pair with Bluetooth:** Enable Bluetooth on your smartphone. Step onto the scale to activate it. Within the app, navigate to the device pairing section and select your Leaone scale from the list of available Bluetooth devices. The app will guide you through the pairing process.
7. **Complete Profile:** Enter your personal information (age, height, gender) into the app. This data is crucial for accurate body composition calculations.



Image: The Leaone BF8028-B Smart Bluetooth Body Weight Scale, showing its sleek design and display.

OPERATING INSTRUCTIONS

Using your Leaone Smart Bluetooth Body Weight Scale is straightforward:

1. Taking a Measurement:

- Ensure your feet are clean and dry.
- Step onto the scale barefoot. The scale will automatically turn on and begin measuring.
- Stand still on the scale until your weight is displayed and stabilized. For body composition analysis, ensure your feet are in contact with the two highly-sensitive electrodes on the scale's surface.
- The scale will display your weight. If the app is open and connected via Bluetooth, all 12 body composition data points will sync automatically to your smartphone. If the app is not open, the weight data will be stored on the scale and sync the next time the app is connected.

2. **Changing Units:** The scale supports units in pounds (LB), kilograms (KG), and stones (ST). You can typically change the unit setting within the Leaone app. Some models may have a physical button on the underside of the scale for this purpose.

3. **Multiple Users:** The Leaone app allows for the creation of multiple user profiles. When different users step on the scale, the app will intelligently identify the user based on their weight history or prompt you to select the correct user.



Image: The Leaone Smart Scale displaying a weight measurement on its digital screen.



Image: A screenshot of the Leaone smartphone application interface, showing various body composition metrics.

MAINTENANCE

Proper care ensures the longevity and accuracy of your scale:

- **Cleaning:** Wipe the scale's tempered glass surface with a soft, damp cloth. Avoid abrasive cleaners or harsh chemicals, as these can damage the surface or electrodes.
- **Storage:** Store the scale in a cool, dry place. Avoid extreme temperatures or direct sunlight.
- **Battery Replacement:** When the low battery indicator appears on the display or in the app, replace all batteries with new ones of the specified type. Ensure correct polarity during installation.
- **Avoid Impact:** While the scale features durable tempered glass, avoid dropping it or subjecting it to heavy impact, which could cause damage.

TROUBLESHOOTING

If you encounter issues with your Leaone Smart Bluetooth Body Weight Scale, refer to the following common problems and solutions:

- **Scale Does Not Turn On:**
 - Check if batteries are installed correctly.
 - Replace old batteries with new ones.
- **"Err" Message on Display:**

- This indicates an overload. The scale's maximum capacity is 400 lbs (180 kg). Remove excess weight.
- **Inaccurate Readings:**
 - Ensure the scale is on a hard, flat, and level surface.
 - Perform a recalibration by stepping on the scale briefly, then stepping off and waiting for it to display "0.0" before stepping on again for a measurement.
 - Ensure your feet are clean and dry and making full contact with the electrodes for body composition analysis.
- **Bluetooth Connection Issues / Data Not Syncing:**
 - Ensure Bluetooth is enabled on your smartphone.
 - Make sure the Leaoe app is open and running in the foreground or background.
 - Try restarting your phone's Bluetooth or the phone itself.
 - Re-pair the scale with the app if necessary (refer to Setup Guide).
 - Check if the scale's batteries are low.
- **Scale Displays Metric Units (KG) Instead of Imperial (LB):**
 - Open the Leaoe app and navigate to the settings or profile section to change the unit preference to LB.
 - Some models may have a unit switch button on the underside of the scale.

SPECIFICATIONS

Feature	Detail
Model Number	BF8028-B
Weight Capacity	3 kg - 180 kg (6.6 lbs - 400 lbs)
Measurement Division	100g / 0.2 lb
Units	KG / LB / ST
Display	Digital LCD/LED (depending on model)
Body Composition Metrics	12 data points: Fat, Water, Visceral Fat, Muscle, Bone, BNR, BMI, Weight, Protein, Obesity, Body Age, LBM
Connectivity	Bluetooth (iOS and Android App)
Operating Temperature	-10°C to 60°C (14°F to 140°F)
Material	Tempered Glass
Dimensions	Approximately 10 x 8 x 15 inches
Item Weight	Approximately 2.8 Pounds
Features	Auto zero, automatic switch-off, low voltage prompt, overload prompt

WARRANTY AND SUPPORT

Specific warranty details are not provided in the product information. For warranty inquiries, technical support, or customer service, please contact Leaone directly through their official website or the contact information provided in the product packaging.

Always refer to the official Leaone website or your purchase documentation for the most current warranty terms and support options.