

B1007

Professional Hair Care Brush Velform Power Grow Laser Hair Comb Massager

Instruction Manual

[Introduction](#) [Safety Information](#) [Package Contents](#) [Setup](#) [Operating Instructions](#) [Maintenance](#)
[Troubleshooting](#) [Specifications](#) [Warranty & Support](#)

1. INTRODUCTION

Thank you for purchasing the Velform Power Grow Laser Hair Comb. This device is designed for hair care and massage, aiming to provide a relaxing experience. Please read this manual carefully before use to ensure proper operation and maintenance.

2. SAFETY INFORMATION

Important Safeguards:

- Do not immerse the device in water or any other liquid.
- Keep out of reach of children.
- Do not use if the device appears damaged.
- Consult a healthcare professional before use if you have any pre-existing medical conditions, especially scalp conditions or light sensitivity.
- Use only as directed in this manual.
- This device is intended for personal use only and not for medical treatment.

3. PACKAGE CONTENTS

Please check the package to ensure all items are present and in good condition:

- Velform Power Grow Laser Hair Comb (Main Unit)
- Standard Hair Brush
- Cleaning Brush

- Instruction Manual (this document)



Image 1: The Velform Power Grow Laser Hair Comb main unit and its detachable comb head, showing the red and blue light emitters.



Image 2: Contents of the package, including the Velform Power Grow Laser Hair Comb, a traditional hair brush, and a small cleaning brush.



Image 3: The Velform Power Grow Laser Hair Comb and its accessories as they appear packaged in the box.

4. SETUP

1. **Unpacking:** Carefully remove all components from the packaging.
2. **Battery Installation:** Locate the battery compartment (usually at the handle). Insert the required batteries (not included) according to the polarity markings.
3. **Attachment:** Ensure the comb head is securely attached to the main unit.

5. OPERATING INSTRUCTIONS

The Velform Power Grow Laser Hair Comb is designed for ease of use. Follow these steps for optimal results:

1. **Preparation:** Ensure your hair is clean and dry before use.
2. **Power On:** Press the "POWER" button located on the device handle. The indicator lights (red and blue) on the comb head should illuminate.
3. **Application:** Gently comb through your hair, moving the device slowly from the front of your scalp to the back. Ensure the comb teeth make contact with your scalp.
4. **Duration:** Use for approximately 10-15 minutes per session, 3-4 times a week. Do not exceed recommended usage.
5. **Power Off:** Press the "POWER" button again to turn off the device after your session.

Note: The device is intended for massage and general hair care. Individual results may vary.

6. MAINTENANCE

- **Cleaning:** After each use, turn off the device and use the provided cleaning brush to remove any hair or debris from the comb teeth and light emitters.
- **Storage:** Store the device in a cool, dry place, away from direct sunlight and moisture.
- **Battery Replacement:** If the device's power weakens or lights dim, replace the batteries as per the setup instructions.
- **Do Not Disassemble:** Do not attempt to disassemble or repair the device yourself. This will void any potential warranty.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Device does not turn on.	Low or dead batteries; Incorrect battery installation.	Replace batteries; Check battery polarity.
Lights are dim.	Low battery power.	Replace batteries.
Comb head feels loose.	Improper attachment.	Ensure comb head is securely clicked into place.

8. SPECIFICATIONS

- **Model Number:** B1007
- **Item Type:** Massage & Relaxation
- **Application:** Head
- **Material:** Non-woven
- **Size:** Medium
- **Power Source:** (Not specified, typically batteries)

Note: Specifications are subject to change without prior notice for product improvement.

9. WARRANTY & SUPPORT

This product is manufactured by Generic. For any support inquiries, troubleshooting assistance beyond this manual, or warranty information, please contact your retailer or the manufacturer directly.

Please retain your purchase receipt as proof of purchase for any warranty claims.

Contact Information:

- **Brand:** Generic
- **ASIN:** B0727LGG2X
- For customer service, please refer to the contact details provided by your point of purchase.

