

[Manuals.plus](#) /

› [NORLANYA BEAUTY EQUIPMENT](#) /

› NORLANYA BEAUTY EQUIPMENT 3MHz Ultra Photon Rejuvenation 3-Color Light Facial Care Massager Instruction Manual

NORLANYA BEAUTY EQUIPMENT 3MHz Ultra Photon Rejuvenation 3-Color Light Facial Care Massager

NORLANYA BEAUTY EQUIPMENT 3MHz Ultra Photon Rejuvenation 3-Color Light Facial Care Massager Instruction Manual

Model: 3MHz Ultra Photon Rejuvenation 3-Color Light Facial Care Massager

1. IMPORTANT SAFETY INFORMATION

Please read all instructions carefully before using the device. Keep this manual for future reference.

Contraindications:

- Do not use if pregnant.
- Do not use if suffering from serious health or skin diseases.
- Avoid use on skin areas that have undergone surgery, implants, or injections within the last 3 months.
- Do not use if hypersensitive to photon light, metal, or heating.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This device is intended for cosmetic use only.

2. PRODUCT OVERVIEW

The NORLANYA BEAUTY EQUIPMENT 3MHz Ultra Photon Rejuvenation 3-Color Light Facial Care Massager is a skin toning device designed for facial beautification. It utilizes natural light waves transmitted by LEDs into the skin, combined with ultrasonic technology, to support skin renewal and rejuvenation. This device aims to address concerns such as wrinkles, spots, pigment irregularities, and sagging skin.



Figure 2.1: Front view of the NORLANYA 3MHz Ultra Photon Rejuvenation Device. The device is white with an orange accent on the side and a clear treatment head with a metallic center.

3. PACKAGE CONTENTS

Verify that all items are present in the package:

- 1 x 3MHz Ultra Photon Rejuvenation Device
- 1 x Pair of Goggles
- 1 x Power Adapter
- 1 x English Instruction Manual
- 1 x Gift Box (packaging)



Figure 3.1: The device, goggles in a blue case, and power adapter neatly arranged within the product's gift box.

4. SPECIFICATIONS

Feature	Detail
Brand	NORLANYA BEAUTY EQUIPMENT
Model	3MHz Ultra Photon Rejuvenation 3-Color Light Facial Care Massager
Product Dimensions	1.77 x 1.97 x 6.5 inches
Item Weight	4.94 ounces (0.14 Kilograms)
Color Options	Red, Green, Blue (LED lights)
Power Source	Power Adapter
Red Light Wavelength	625 ±3nm
Green Light Wavelength	525 ±3nm
Blue Light Wavelength	465 ±3nm

5. SETUP

- Unpack the Device:** Carefully remove the device and all accessories from the packaging.
- Inspect for Damage:** Check the device and power adapter for any visible damage. If any damage is found, do not use the device and contact customer support.
- Connect Power:** Plug the power adapter into the device and then into a suitable electrical outlet. The device requires continuous power from the adapter during operation.
- Wear Goggles:** Always wear the provided goggles to protect your eyes from the LED light during use.

6. OPERATING INSTRUCTIONS

Follow these steps for safe and effective use of your Photon Rejuvenation Device:

- Prepare Skin:** Cleanse your face thoroughly to remove makeup, dirt, and oils. Pat dry.
- Apply Product (Optional):** For enhanced results, apply a thin layer of your preferred serum or gel to the treatment area.
- Power On:** Press the **ON/OFF** button to turn on the device.
- Select Mode:**
 - Press the **Sonic** button to activate the ultrasonic function.
 - Press the **Photon** button to cycle through the different LED light colors (Red, Green, Blue). Refer to Section 7 for details on each light mode.
- Adjust Intensity:** Use the 'High', 'Middle', 'Low' settings to adjust the intensity of the treatment as desired.
- Begin Treatment:** Gently move the treatment head over your skin in small, circular motions. Avoid prolonged contact with one area.
- Treatment Duration:** A typical session lasts 10-15 minutes per area. Do not exceed 20 minutes per session.
- Power Off:** After your session, press the **ON/OFF** button to turn off the device.

9. **Clean Device:** Wipe the treatment head with a soft, damp cloth. Do not immerse the device in water.



Figure 6.1: The control panel showing 'High', 'Middle', 'Low' indicators, 'Sonic' and 'Photon' buttons, and the 'ON/OFF' button.



Figure 6.2: A detailed view of the device's treatment head, featuring multiple LED lights and a central metallic contact point.

7. LIGHT THERAPY MODES

The device features three distinct LED light modes, each targeting specific skin concerns:

Red Light (625 ±3nm)

The red light mode is designed to activate skin cells, promote the formulation of collagen, and help reduce the appearance

of wrinkles, particularly around the eyes and in eye pouches.



Figure 7.1: The device's treatment head illuminated with red LED light, indicating the red light therapy mode is active.

Green Light (525 ±3nm)

The green light mode assists in balancing water and oil levels in the skin, making it suitable for oily and sensitive skin types. It aims to soothe and harmonize the skin's condition.

Blue Light (465 ±3nm)

The blue light mode works to activate skin cells, reduce melanin production, and contribute to skin brightening and whitening. It also helps to improve skin moisture levels.



Figure 7.2: The device's treatment head illuminated with blue LED light, indicating the blue light therapy mode is active.

8. MAINTENANCE

1. **Cleaning:** After each use, wipe the treatment head with a soft, slightly damp cloth. Do not use abrasive cleaners or immerse the device in water. Ensure the device is completely dry before storage.
2. **Storage:** Store the device in its original gift box or a clean, dry place away from direct sunlight and extreme temperatures. Keep out of reach of children.

9. TROUBLESHOOTING

If you encounter any issues with your device, please refer to the following common solutions:

- **Device Not Turning On:**

- Ensure the power adapter is securely connected to both the device and a working electrical outlet.
- Verify that the power outlet is functional by testing it with another appliance.

- **Lights Not Emitting:**

- Confirm the device is powered on.
- Press the 'Photon' button repeatedly to cycle through the light modes.

- **No Ultrasonic Sensation:**

- Ensure the 'Sonic' button has been pressed to activate the ultrasonic function.
- Ultrasonic vibrations are often subtle and may not be strongly felt. Apply a small amount of water to the treatment head; if it atomizes, the function is working.

If the problem persists after attempting these solutions, please contact customer support.

10. WARRANTY AND SUPPORT

Warranty Information:

This device comes with a one-year warranty from the date of purchase. Please retain your proof of purchase for warranty claims.

Return Policy:

A 30-day refund/replacement return policy is typically offered from the date of purchase. Please check with your retailer for specific return terms and conditions.

Customer Support:

For technical assistance, warranty claims, or any other inquiries, please contact NORLANYA BEAUTY EQUIPMENT customer support through your purchase platform or the contact information provided with your device.