

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [FOREO](#) /
- › [FOREO ISSA 2 Sonic Electric Toothbrush User Manual \(Model F3616\)](#)

## FOREO ISSA 2

# FOREO ISSA 2 Sonic Electric Toothbrush User Manual

Model: ISSA 2 (F3616) | Brand: FOREO

## INTRODUCTION

---

Welcome to the user manual for your FOREO ISSA 2 Sonic Electric Toothbrush. This device is designed to provide an advanced oral care experience, combining enhanced sonic pulse technology with a unique brush head. Please read these instructions carefully before first use to ensure optimal performance and longevity of your device.





Image: Front view of the FOREO ISSA 2 Sonic Electric Toothbrush in Mint color.



Image: A collection of FOREO ISSA toothbrushes in various colors, showcasing the design.

## SAFETY INFORMATION

---

- This device is intended for oral hygiene as described in this manual. Any other use may result in damage to the device or injury.
- Consult a dentist or physician before use if you have any oral health concerns or medical conditions.
- Do not use the device if it appears damaged in any way.
- Keep out of reach of children. This device is not a toy.
- Avoid using abrasive cleaners or chemicals on the device.

## WHAT'S IN THE BOX

---

Your FOREO ISSA 2 package includes the following items:

- FOREO ISSA 2 Sonic Electric Toothbrush Unit
- USB Charging Cable

**365** DAYS  
JOURS  
WITH ONLY ONE CHARGE  
AVEC UNE SEULE CHARGE

FOREO

ISSA™ 2

WORLD'S 1ST  
SILICONE  
**SONIC**  
TOOTHBRUSH

**LA 1ÈRE**  
BROSSE À DENTS  
**ÉLECTRIQUE**  
EN SILICONE  
**AU MONDE**



Image: The FOREO ISSA 2 Sonic Electric Toothbrush displayed in its original retail packaging.

# SETUP

---

## Charging the Device

Before first use, fully charge your ISSA 2 device. A full charge provides extended usage time.

1. Locate the charging port at the base of the toothbrush.
2. Connect the provided USB charging cable to the device.
3. Plug the other end of the USB cable into any compatible USB power source (e.g., wall adapter, computer USB port).
4. The indicator light on the device will show charging status. Refer to the device's specific light patterns for full charge indication.

The device requires 1 Lithium Polymer battery, which is included and integrated.

## First Use

After charging, your ISSA 2 is ready for use. Ensure the brush head is securely attached.

## OPERATING INSTRUCTIONS

---

### Brushing Your Teeth

1. Dampen the brush head and apply your preferred toothpaste.
2. Press the power button to activate the device.
3. Brush your teeth for two minutes using circular motions, similar to how you would use a manual toothbrush. The device may have a built-in timer to guide you.
4. After brushing, press the power button again to turn off the device.
5. Rinse your mouth and the ISSA 2 with water. If using fluoride toothpaste, you may choose to avoid rinsing your mouth immediately after brushing to maximize fluoride benefits, but always rinse the device.



Image: A user demonstrating the brushing motion with the FOREO ISSA 2 toothbrush.

## Using the Tongue Cleaner

The ISSA 2 features an integrated tongue cleaner on the back of the brush head to help complete your oral care routine.

1. After brushing your teeth, turn the device over so the back of the brush head faces your tongue.
2. Gently glide the tongue cleaner across your tongue from back to front.
3. Rinse your mouth and the device thoroughly after use.





Image: Back view of the FOREO ISSA 2, highlighting the integrated tongue cleaner.

## MAINTENANCE

---

### Cleaning the Device

To maintain hygiene and performance, clean your ISSA 2 after each use.

1. Rinse the brush head and the entire device under running water.
2. Wipe the body of the device with a clean, dry cloth.
3. Avoid using alcohol, petrol, or acetone-based cleaning products, as these may damage the silicone material.



Image: The FOREO ISSA 2 toothbrush being rinsed under running water for cleaning.

### Replacing the Brush Head

For optimal oral hygiene and performance, it is recommended to replace the brush head every six months.

1. Grip the handle of the ISSA 2 firmly.
2. Gently pull the current brush head straight off the handle.
3. Align the new brush head with the handle and push it down until it clicks securely into place.

## TROUBLESHOOTING

---

If you encounter any issues with your FOREO ISSA 2, please refer to the following common solutions:

- **Device not turning on:** Ensure the device is fully charged. If the issue persists, try charging it with a different USB port or adapter.
- **Device not charging:** Check that the USB cable is securely connected to both the device and the power source. Verify the power source is active.
- **Reduced performance:** The brush head may need replacement. Ensure the brush head is securely attached.

For further assistance, please contact FOREO customer support.

## SPECIFICATIONS

---

<b>Brand</b>	FOREO
<b>Model Number</b>	2028686 (F3616)
<b>ASIN</b>	B077V4P76W
<b>Age Range Description</b>	Adult
<b>Power Source</b>	Battery Powered
<b>Battery Type</b>	1 Lithium Polymer battery (included)
<b>Product Dimensions</b>	2.79 x 2.79 x 15.24 cm
<b>Item Weight</b>	77.11 grams
<b>Special Feature</b>	Integrated Tongue Scraper
<b>Product Firmness</b>	Soft
<b>Color</b>	Mint
<b>Brush Head Material</b>	PBT Polymer and Silicone
<b>Toothbrush Head Shape</b>	Round
<b>First Available Date</b>	March 30, 2018

## WARRANTY INFORMATION

---

FOREO products are manufactured with the highest quality standards. For detailed information regarding the warranty coverage for your ISSA 2 device, please refer to the official FOREO website or the warranty card included with your purchase. Keep your proof of purchase for warranty claims.

## CUSTOMER SUPPORT

---

If you have any questions, concerns, or require technical assistance with your FOREO ISSA 2, please contact FOREO customer support. Contact details can typically be found on the official FOREO website or within your product packaging.