

[Manuals.plus](#) /

› [Philips Sonicare](#) /

› Philips Sonicare C2 Plaque Control Replacement Brush Heads HX9023/65 - Instruction Manual

Philips Sonicare HX9023/65

Philips Sonicare C2 Plaque Control Replacement Brush Heads

Model: HX9023/65

INTRODUCTION

This manual provides instructions for the Philips Sonicare C2 Plaque Control replacement brush heads. These brush heads are designed for effective plaque removal and are compatible with Philips Sonicare click-on electric toothbrushes. For optimal oral hygiene, it is recommended to replace brush heads every three months.



Image: Philips Sonicare C2 Plaque Control replacement brush heads, showing a 3-pack in its retail packaging and three individual brush heads.

PRODUCT FEATURES

- **Effective Plaque Removal:** Densely packed, high-quality tufts and a specially curved power tip are designed to remove up to 700% more plaque than a manual toothbrush, targeting hard-to-reach areas.
- **Sonicare Fluid Action:** The brush head works with Sonicare technology to drive fluid deep between teeth and along the gumline, assisting in thorough cleaning.
- **Clinically Tested:** All Philips Sonicare brush heads undergo rigorous clinical testing to ensure

high standards of quality, safety, and performance.

- **BrushSync Compatibility:** For compatible handles, BrushSync technology tracks usage and reminds you when to replace the brush head. Blue indicator bristles also fade to white, signaling replacement is due.
- **Sustainable Design:** The brush head is made with 70% bio-based plastic, and the packaging is paper-based and recyclable.

Brush head's power tip is designed to target plaque in hard-to-reach areas



Image: A close-up view of the Philips Sonicare C2 brush head, illustrating its power tip and bristles effectively targeting plaque on teeth.

Sonicare Fluid Action

cleans deep between the teeth
and along the gumline



Image: An illustration depicting Sonicare Fluid Action, showing water and toothpaste bubbles reaching deep between teeth and along the gumline for comprehensive cleaning.

SETUP AND INSTALLATION

1. **Remove Old Brush Head:** Grip the old brush head and pull it straight off the toothbrush handle.
2. **Attach New Brush Head:** Align the new C2 Plaque Control brush head with the metal shaft of your Philips Sonicare handle. Push the brush head down firmly until it clicks into place. Ensure it is securely attached.

These brush heads are compatible with all Philips Sonicare click-on electric toothbrush handles, with the exception of Philips One and e-series handles.

All genuine Philips Sonicare brush heads fit all Philips Sonicare click-on handles⁵



Image: A diagram illustrating how Philips Sonicare brush heads, including the C2 Plaque Control, click onto various compatible Philips Sonicare electric toothbrush handles for easy replacement.

OPERATING INSTRUCTIONS

Once the brush head is securely attached, apply a small amount of toothpaste to the bristles. Moisten the brush head. Place the bristles against your teeth at a slight angle towards the gumline. Turn on your Philips Sonicare toothbrush and gently move the brush head slowly across each tooth surface. Allow the sonic vibrations to do the cleaning. Refer to your specific Philips Sonicare toothbrush handle's instruction manual for detailed brushing techniques and modes.



The brush head for thorough cleaning

Image: A person demonstrating the use of a Philips Sonicare electric toothbrush with a C2 Plaque Control brush head, showing the brush head positioned on the teeth for thorough cleaning.

MAINTENANCE

- **Rinsing:** After each use, rinse the brush head thoroughly under running water to remove any toothpaste residue and debris.
- **Air Drying:** Allow the brush head to air dry in an upright position.
- **Replacement Frequency:** Dental professionals recommend replacing your brush head every 3 months.
- **Replacement Indicators:**
 - **BrushSync Reminder:** If your toothbrush handle is BrushSync-compatible, a light on the handle will illuminate when it's time for replacement.
 - **Blue Bristle Fade:** The blue indicator bristles on the brush head will fade to white over time, signaling that it's time to replace the brush head.



Image: A Philips Sonicare C2 brush head with an icon indicating the recommendation to replace the brush head every 3 months.

COMPATIBILITY

The Philips Sonicare C2 Plaque Control brush heads are designed to fit all Philips Sonicare click-on electric toothbrush handles. This includes most models within the Sonicare range. Please note that these brush heads are **not compatible** with Philips One and e-series handles.

PHILIPS sonicare							
Find your brush head:	C2 Plaque Control	W Diamond Clean	S2 Sensitive	G3 Two-in-One Plaque Control + Gum Care	C3 Two-in-One Plaque Control + Whitening	A3 Premium All-in-one	
Plaque removal vs. manual toothbrush	700%	700%	1000%	1000%	1000%	2000%	
Focus on healthier gums				✓		✓	
Focus on whiter teeth		✓			✓	✓	
Focus on sensitive teeth			✓				
Bristle stiffness	Soft	Medium soft	Extra soft	Soft	Medium soft	Soft	
Flexible rubber back for more surface contact				✓	✓	✓	
Colors	White	White/Black	White	White/Black	White/Black	White/Black	
Compatible with all Philips Sonicare handles**	✓	✓	✓	✓	✓	✓	
70% bio-based plastic***	✓	✓	✓	✓	✓	✓	
BrushSync	✓	✓	✓	✓	✓	✓	

*In hard-to-reach areas / **Except for Philips One / ***allocated to plastic of brush head on mass balance basis

Image: A comparison table detailing various Philips Sonicare brush head types, including C2 Plaque Control, W DiamondClean, S2 Sensitive, G3 Two-in-One Plaque Control + Gum Care, C3 Two-in-One Plaque Control + Whitening, and A3 Premium All-in-one, highlighting their features and compatibility.

SPECIFICATIONS

Attribute	Detail
Model Number	HX9023/65
Brand	Philips Sonicare
Color	White
Size	3 Count (Pack of 1)
Material	Nylon
Product Dimensions	0.98 x 2.36 x 8.7 inches
Item Weight	1.41 ounces
Manufacturer	Philips

SAFETY INFORMATION

These brush heads are crafted to be safe and gentle, ensuring effective cleaning without causing harm to teeth or gums. If you experience any discomfort or irritation, discontinue use and consult a dental professional.

LEGAL DISCLAIMER

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is for oral hygiene purposes only.

WARRANTY AND SUPPORT

For information regarding product warranty, returns, or customer support, please refer to the official Philips Sonicare website or the documentation included with your Philips Sonicare electric toothbrush handle. You may also contact Philips customer service directly for assistance.

© 2025 Philips Sonicare. All rights reserved.