

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [ProfiCare](#) /

› [ProfiCare PC-PW 3008 BT 9-in-1 Diagnostic Scale User Manual](#)

ProfiCare PC-PW 3008 BT

ProfiCare PC-PW 3008 BT 9-in-1 Diagnostic Scale User Manual

Model: PC-PW 3008 BT (330080)

1. INTRODUCTION

The ProfiCare PC-PW 3008 BT is a sophisticated 9-in-1 diagnostic scale designed to provide comprehensive insights into your body composition. Utilizing modern Bluetooth Smart technology, this scale measures not only your weight but also calculates body fat percentage, water content, muscle mass, bone mass, organ fat, calorie needs, Body Mass Index (BMI), and biological age. All data is transmitted directly to a dedicated smartphone application for continuous tracking and evaluation.



Image 1.1: ProfiCare PC-PW 3008 BT Diagnostic Scale with key features such as easy-to-read display, 180kg measurement range, step-on function, 10-user memory, low battery indicator, and Bluetooth smart technology.

2. SAFETY INFORMATION

- This device is intended for personal use only and not for medical diagnosis or treatment.
- Do not use this device if you have a pacemaker or other internal medical device.
- Keep the scale away from water and extreme temperatures.
- The scale surface can be slippery when wet. Step on and off carefully.
- Do not jump on the scale.
- Keep out of reach of children.
- Dispose of batteries according to local regulations.

3. PACKAGE CONTENTS

Please check that all items are present and in good condition:

- ProfiCare PC-PW 3008 BT Diagnostic Scale
- User Manual

Note: 3 x 1.5V AAA batteries are not included and must be purchased separately.

4. SETUP

4.1 Battery Installation

1. Turn the scale over to access the battery compartment on the underside.
2. Open the battery compartment cover.
3. Insert three (3) 1.5V AAA batteries, ensuring correct polarity (+/-).
4. Close the battery compartment cover securely.

4.2 App Installation

1. Search for the official ProfiCare app on your smartphone's app store (e.g., Apple App Store or Google Play Store).
2. Download and install the free application.
3. Follow the on-screen instructions to create a user profile within the app.

4.3 Bluetooth Pairing

1. Ensure Bluetooth is enabled on your smartphone or tablet.
2. Open the ProfiCare app.
3. Step onto the scale to activate it. The display will light up.
4. The app should automatically detect and connect to the scale. If not, follow the pairing instructions provided within the app.

5. OPERATING INSTRUCTIONS

5.1 Taking a Measurement

1. Place the scale on a hard, flat surface. Avoid carpets or uneven floors.
2. Ensure your feet are clean and dry.
3. Step onto the scale barefoot. Stand still with your weight evenly distributed.
4. The scale will automatically turn on (Step-On function) and display your weight.
5. After a few seconds, the scale will perform the full diagnostic analysis.
6. Your measurements will be displayed on the LCD and automatically transmitted to the ProfiCare app via Bluetooth.

5.2 Switching Measurement Units

The scale supports measurement in kilograms (kg), pounds (lb), and stone:pounds (st:lb). To switch units, locate the unit switch button on the underside of the scale (usually near the battery compartment) and press it to cycle through the available units.

5.3 User Profiles

The scale can store data for up to 10 different users. User profiles are managed through the ProfiCare app. Ensure you select the correct user profile in the app before taking a measurement to ensure data is assigned accurately.



Easy-to-read
high quality
multi-function



Measurement
range



Step-on
automatic
function



Memory for up
to 10 different
users



Low battery
indicator



Bluetooth
smart
technology



Image 5.1: Close-up view of the ProfiCare PC-PW 3008 BT scale, showing the stainless steel surface, LCD display, and touch controls for navigation.

6. DATA INTERPRETATION

The ProfiCare app provides detailed explanations and historical data for each measurement. Consult the app for personalized insights and trends.

- **Weight:** Your total body mass.
- **Body Fat Percentage:** The proportion of fat to total body weight.
- **Water Content:** The percentage of water in your body.
- **Muscle Mass:** The total weight of muscle in your body.
- **Bone Mass:** The total weight of bone in your body.
- **Organ Fat (Visceral Fat):** Fat stored around the organs in the abdominal cavity.
- **Calorie Needs (BMR/TDEE):** An estimate of the calories your body needs daily.
- **BMI (Body Mass Index):** A measure of body fat based on height and weight.
- **Biological Age:** An estimate of your body's physiological age compared to your chronological age.

7. MAINTENANCE AND CARE

- Clean the scale surface with a soft, damp cloth. Do not use abrasive cleaners or immerse the scale in water.
- The stainless steel surface is designed to be hygienic and easy to clean.
- If the scale will not be used for an extended period, remove the batteries to prevent leakage.
- Store the scale in a cool, dry place.
- The scale features an automatic shut-off function to conserve battery life. A low battery indicator will appear on the display when batteries need replacement.

8. TROUBLESHOOTING

Problem	Possible Cause	Solution
Scale does not turn on.	Batteries are dead or incorrectly inserted.	Replace batteries (3x AAA) or check polarity.
Inaccurate weight readings.	Scale on uneven surface; scale moved; calibration needed.	Place on a hard, flat surface. Step on and off to re-calibrate to zero.
Bluetooth connection issues.	Bluetooth off on device; app not open; scale not active.	Ensure Bluetooth is on. Open app. Step on scale to activate. Restart app/phone if needed.
Diagnostic data not displayed/transmitted.	Not standing barefoot; feet not making proper contact.	Ensure you are barefoot and standing still on the electrodes.

9. SPECIFICATIONS

- **Model:** PC-PW 3008 BT
- **Model Number:** 330080
- **EAN/UPC:** 04006160300802
- **Measurement Range:** 5 kg to 180 kg (11 lb to 396 lb)
- **Measurement Increments:** 100 g (0.2 lb)
- **Measurement Units:** kg, lb, st:lb (switchable)
- **Display:** High-quality multifunction LCD (black/white)
- **User Memory:** Up to 10 different users
- **Power Supply:** 3 x 1.5V AAA batteries (not included)
- **Surface Material:** Extra-large, high-quality stainless steel (approx. 30 x 30 cm)
- **Features:** Modern touch control, automatic Step-On function, anti-slip feet, automatic shut-off, low battery indicator.
- **Product Weight:** 1.7 kg
- **Country of Origin:** China



Image 9.1: The ProfiCare PC-PW 3008 BT scale being held, showcasing its sleek design and compact form factor.

10. WARRANTY AND SUPPORT

This product comes with a 1-year availability for spare parts. For warranty claims, technical support, or further assistance, please contact ProfiCare customer service. Refer to the official ProfiCare website or your purchase documentation for contact details.