

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Microderm GLO](#) /
- › [Microderm GLO MINI Premium Diamond Microdermabrasion Tips Instruction Manual](#)

### Microderm GLO MDG-TIP-PR1804

# Microderm GLO MINI Premium Diamond Microdermabrasion Tips Instruction Manual

Brand: Microderm GLO | Model: MDG-TIP-PR1804

## 1. INTRODUCTION

The Microderm GLO MINI Premium Diamond Microdermabrasion Tips are medical-grade stainless steel accessories designed for use with the Microderm GLO MINI System. These tips feature patented Safe3D Technology, ensuring safe and effective treatments for all skin types.

Regular replacement of these diamond tips is crucial for maintaining professional treatment results and optimizing your skincare routine.

# Introducing Microderm GLO MINI Premium Diamond Tip

Brightens skin tone.  
Rebuilds collagen & elastin.

Medical-grade diamond  
exfoliation on the go.



\*Compatible with the  
Microderm GLO MINI system

*Image: Close-up of the Microderm GLO MINI Premium Diamond Microdermabrasion Tip. This tip is essential for effective exfoliation and skin renewal.*

## 2. COMPATIBILITY

These Premium Diamond Microdermabrasion Tips are exclusively designed for and compatible with the **Microderm GLO MINI System**. They are not intended for use with the original Microderm GLO device or any other microdermabrasion systems.

## 3. SETUP AND INSTALLATION

To ensure proper function and hygiene, follow these steps for installing a new diamond tip:

1. **Remove Used Tip:** Gently pull off the old diamond tip from your Microderm GLO MINI device.
2. **Insert New Filter:** Install a new filter disc into the replacement diamond tip. Ensure it is seated correctly.
3. **Attach New Tip:** Insert the new diamond tip into the Microderm GLO MINI by gently pushing it back into the unit until it is secure.



Visible results  
after one use.

Guaranteed.

*Image: The Microderm GLO MINI device held in hand, demonstrating its ergonomic design. The tip is visible at the top.*

#### 4. OPERATING INSTRUCTIONS

For optimal results and to ensure skin safety, adhere to the following operating procedures:

1. **Prepare Skin:** Wash your face thoroughly and pat it dry before starting the treatment.
2. **Position Device:** Brace the skin with your free hand. Place the diamond tip flat against the skin surface.
3. **Apply Strokes:** Using short, even strokes, gently move the diamond tip across the skin. Ensure the tip remains flush with the skin at all times. Avoid lingering in one spot for too long to prevent irritation.

4. **Post-Treatment Care:** After completing the microdermabrasion treatment, apply a Peptide Complex Serum to hydrate, repair, and moisturize your skin. For enhanced deep cleaning, consider using a facial steamer before treatments.



*Image: A woman demonstrating the proper technique for using the Microderm GLO device on her face, moving the tip gently across the skin.*

## 5. MAINTENANCE

To maintain the effectiveness and hygiene of your Microderm GLO MINI, regular maintenance is essential:

- **Tip Replacement:** Diamond tips should be replaced approximately every 3 months with normal use to ensure optimal exfoliation and treatment quality.
- **Filter Replacement:** Always use a new filter disc with each new diamond tip. Discard used filter discs after each treatment.
- **Cleaning:** After each use, clean the diamond tip and the device according to the main Microderm GLO MINI System manual to prevent buildup and maintain performance.

## 6. TROUBLESHOOTING

If you encounter any issues with your Microderm GLO MINI Premium Diamond Microdermabrasion Tips, consider the following:

- **Reduced Effectiveness:** If you notice a decrease in exfoliation or suction, it may be time to replace your diamond tip and filter disc.
- **Skin Irritation:** Ensure you are not applying too much pressure or lingering in one spot. Always move the device in short, even strokes. If irritation persists, discontinue use and consult a skincare professional.

- **Device Malfunction:** For issues related to the Microderm GLO MINI System itself, refer to the main device's instruction manual or contact Microderm GLO customer support.

## 7. SPECIFICATIONS

Feature	Detail
Item Model Number	MDG-TIP-PR1804
Manufacturer	Microderm GLO
Material	Medical-Grade Stainless Steel
Technology	Patented Safe3D Technology
Compatibility	Exclusively with Microderm GLO MINI System
Recommended Replacement	Every 3 months with normal use

## 8. WARRANTY AND SUPPORT

For warranty information, product support, or to purchase replacement tips and filters, please visit the official Microderm GLO website or contact their customer service directly. Refer to your original Microderm GLO MINI System manual for specific warranty details.

You can also visit the [Microderm GLO Store on Amazon](#) for additional products and support.

## 9. PRODUCT VIDEOS

### Microderm GLO Brand Showcase

Your browser does not support the video tag.

*Video: An official brand showcase video from Microderm GLO, highlighting the benefits and features of their advanced professional skincare products.*