

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [TOUCHBeauty](#) /
- › [TOUCHBeauty 3D Facial and Body Roller Massager AG-1613A Instruction Manual](#)

TOUCHBeauty AG-1613A

TOUCHBeauty 3D Facial and Body Roller Massager

MODEL: AG-1613A

1. Product Overview

The TOUCHBeauty 3D Facial and Body Roller Massager is designed to provide a comprehensive massage experience for both the face and body. Its unique V-shaped and 360-degree rotating rollers are engineered to gently lift and firm the skin, promote blood circulation, and help reduce the appearance of fine lines and wrinkles.



Image 1.1: The TOUCHBeauty 3D Facial and Body Roller Massager.

Key Features:

- **70-Degree V-Shaped Lifting:** Unique rollers designed to massage various body parts, aiding in skin firming and lifting.
- **Enhanced Skin Elasticity:** Helps improve skin elasticity and cellular activity, reducing fine lines and wrinkles.
- **3D Roller Design:** Features 360-degree rotating rollers for a deep, double-contact massage that tightens muscles and restores skin firmness.
- **Versatile Use:** Suitable for face, neck, arms, waist, and legs to relax muscles and relieve fatigue.
- **Ergonomic and Portable:** Compact and lightweight design for comfortable handling and easy portability.

2. What's in the Box

Upon opening your product packaging, you will find the following components:

- 1 x TOUCHBeauty 3D Facial and Body Roller Massager
- 1 x Storage Bag
- 1 x Cleaning Cloth

- 1 x Instruction Manual

Video 2.1: Unboxing and close-up view of the TOUCHBeauty 3D Facial and Body Roller Massager, including the instruction manual and cleaning cloth.

3. Setup

The TOUCHBeauty 3D Facial and Body Roller Massager requires minimal setup. It is a manual device and does not require batteries or charging.

1. Remove the massager from its packaging.
2. Before first use, gently wipe the rollers with the provided cleaning cloth to ensure they are clean.
3. Ensure your skin is clean and free of heavy makeup or oils that could hinder the roller's movement. For best results, apply a thin layer of your preferred skincare product (serum, oil, or moisturizer) to the area you wish to massage.

4. Operating Instructions

The 3D Facial and Body Roller Massager is designed for easy and effective use. Follow these guidelines for optimal results:

General Usage Tips:

- Hold the massager firmly but comfortably.
- Apply gentle pressure as you roll the device over your skin. The rollers are designed to gently grasp and lift the skin.
- Always roll in an upward or outward direction to promote lifting and contouring.
- Use daily, either in the morning or evening, as part of your skincare routine.

3D Facial & Body Roller Massager



Image 4.1: Recommended directions for facial and neck massage.

Specific Application Areas:

- **Face:** Start from the chin and roll upwards towards the ears. For the cheeks, roll from the corners of the mouth towards the temples. For the forehead, roll from the eyebrows upwards to the hairline.
- **Neck:** Roll from the base of the neck upwards towards the jawline.
- **Arms:** Roll from the elbow towards the shoulder, focusing on areas needing firming.
- **Waist/Abdomen:** Roll in upward motions or circular motions across the abdomen.
- **Legs:** Roll from the ankle upwards towards the knee and thigh.

For Face and Body Care



Face



Arms



Waist



Legs

Image 4.2: Various body parts where the massager can be used.

Video 4.3: Demonstration of the TOUCHBeauty 3D roller massager on the face, neck, arms, and other body areas.

5. Maintenance

Proper maintenance ensures the longevity and hygiene of your TOUCHBeauty 3D Facial and Body Roller Massager.

- **Cleaning:** After each use, wipe the rollers with the provided cleaning cloth or a soft, damp cloth. Ensure no skincare product residue remains on the rollers.
- **Drying:** Allow the massager to air dry completely before storing.
- **Storage:** Store the massager in the provided storage bag to protect it from dust and scratches. Keep it in a cool, dry place.
- **Avoid Harsh Chemicals:** Do not use abrasive cleaners, alcohol, or harsh chemicals, as these can damage the product's finish.

6. Troubleshooting

The TOUCHBeauty 3D Facial and Body Roller Massager is a simple, non-electronic device. Most issues can be resolved with basic checks.

- **Rollers Not Moving Smoothly:** Ensure the rollers are clean and free of any debris or dried skincare product. Clean thoroughly with the provided cloth.
- **Discomfort During Use:** Reduce the pressure applied. The massage should be comfortable and gentle. Ensure your skin is adequately moisturized or lubricated with a serum or oil.
- **Product Appears Damaged:** If the product shows any signs of damage, discontinue use immediately. Do not attempt to repair it yourself.

7. Specifications

Feature	Detail
Brand	TOUCHBeauty
Model Number	AG-1613A
Product Dimensions	15.7 x 8.7 x 5.4 cm
Weight	240 grams
Material	Plastic
Style	Compact
Power Source	Battery powered (Note: This device is manual and does not require batteries for operation.)
Included Components	Massager, Storage Bag, Cleaning Cloth, Instruction Manual
Country of Origin	China

8. Safety Information

Please read and understand all safety instructions before using the TOUCHBeauty 3D Facial and Body Roller Massager. Keep this manual for future reference.

Warnings:

- Do not use this product if you are pregnant or breastfeeding.
- Avoid use on areas with allergies, chapped skin, lumps, severe skin diseases, or any other skin conditions.
- Inspect the product for any damage before each use. Do not use if it appears damaged.
- Discontinue use immediately and consult a doctor if you experience any discomfort or feel unwell during or after use.
- Do not attempt to modify, disassemble, or repair the product yourself.
- This product is for personal use only. Do not share it with others to maintain hygiene.
- Keep out of reach of children.

9. Warranty and Support

For warranty information, product support, or any inquiries regarding your TOUCHBeauty 3D Facial and Body Roller Massager, please refer to the official TOUCHBeauty website or contact their customer service directly.

You can visit the TOUCHBeauty Store for more information: [TOUCHBeauty Store](#)