

[Manuals.plus](#) /

› [Tao Clean](#) /

› Tao Clean UV Sanitizing Sonic Toothbrush and Cleaning Station User Manual (Model BA-0130-BLK-US)

Tao Clean BA-0130-BLK-US

Tao Clean UV Sanitizing Sonic Toothbrush and Cleaning Station User Manual

MODEL: BA-0130-BLK-US

Brand: Tao Clean

1. Introduction

This manual provides detailed instructions for the proper use and maintenance of your Tao Clean UV Sanitizing Sonic Toothbrush and Cleaning Station. This system is designed to provide an effective oral cleaning experience combined with hygienic brush head sanitization.



Image: The Tao Clean UV Sanitizing Sonic Toothbrush and Cleaning Station in black, showcasing both the toothbrush handle and its accompanying base station.

2. Safety Information

WARNING: [GENERAL SAFETY WARNINGS APPLY]

Please read all instructions carefully before using this product. Keep this manual for future reference. Ensure the power adapter and cable are not damaged. Do not immerse the main unit or charging base in water. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.

3. Package Contents

- 1 x Tao Clean Sonic Toothbrush Handle
- 1 x Tao Clean Base Station (UV Sanitizing and Charging)

- 1 x Hygienically Sealed Toothbrush Head
- 1 x Power Adaptor and Cable

4. Product Overview

4.1 Components

- **Sonic Toothbrush Handle:** Ergonomically designed for comfortable grip, housing the motor and control buttons.
- **Toothbrush Head:** Replaceable brush head with bristles designed for effective cleaning.
- **Base Station:** Serves as a charging dock, UV sanitizing unit, and drying station for the toothbrush head.
- **Power Adaptor:** Connects the base station to an electrical outlet for power.

4.2 Key Features

- **Sonic Technology:** Delivers 40,000 brush strokes per minute for efficient plaque removal and polishing.
- **UV-C Sanitization:** The patented docking station uses UV-C rays to kill 99.9% of bacteria on the brush head.
- **Built-in Dryer:** Wicks away moisture from the brush head after use, preventing germ growth.
- **Dual Speed Settings:** Offers 'Super Clean' for thorough cleaning and 'Sensitive Clean' for gentle brushing.
- **2-Minute Smart Timer:** Ensures optimal brushing time, with 30-second pauses to indicate when to move to the next quadrant of your mouth.

Find your perfect routine
with two speed settings:
Super Clean & Sensitive Clean



Image: Close-up of the Tao Clean sonic toothbrush handle, illustrating the power button and the two speed setting controls.



A patented dock uses **UV-C** light to sanitize the brush, ensuring you have a clean brush **every time**.

Image: An internal view of the Tao Clean base station, highlighting the UV-C light actively sanitizing the toothbrush head.

5. Setup

1. **Unpack:** Carefully remove the toothbrush handle, base station, toothbrush head, and power adapter from their packaging.
2. **Attach Brush Head:** Securely attach one hygienically sealed toothbrush head to the toothbrush handle. Ensure it clicks firmly into place.
3. **Connect Power:** Plug the power adapter into the base station, then plug the adapter into a suitable electrical outlet.
4. **Initial Charge:** Place the toothbrush handle onto the base station. The base station will begin charging the toothbrush and initiate the UV sanitizing cycle. Allow for a full initial charge before first use (approximately 24 hours).

6. Operation

6.1 Brushing Your Teeth

1. **Wet Brush Head:** Wet the toothbrush head and apply a small amount of toothpaste.
2. **Power On:** Press the power button on the toothbrush handle to turn it on.
3. **Select Mode:** Press the mode button to cycle between 'Super Clean' and 'Sensitive Clean' settings. Choose the mode suitable for your preference.
4. **Brush:** Gently guide the toothbrush along your gum line, allowing the sonic vibrations to clean your teeth. Do not scrub aggressively.
5. **Quadrant Timer:** The toothbrush will briefly pause every 30 seconds to remind you to move to the next quadrant of your mouth. Brush each of the four quadrants (upper right, upper left, lower right, lower left) for 30 seconds.
6. **Auto Shut-off:** The toothbrush will automatically turn off after 2 minutes, indicating the completion of the recommended brushing cycle.

UPGRADE

your daily brushing routine with
powerful sonic cleaning.



Image: A person brushing their teeth with the Tao Clean sonic toothbrush, illustrating the proper technique for effective cleaning.

6.2 After Brushing

After brushing, rinse the toothbrush head thoroughly under running water. Then, place the damp toothbrush handle

upside down into the base station. The base station will automatically initiate a UV sanitizing and drying cycle to ensure your brush head is clean and dry for the next use.

7. Maintenance

7.1 Cleaning the Toothbrush Handle

Wipe the toothbrush handle with a damp cloth as needed. Do not immerse the handle in water.

7.2 Cleaning the Base Station

Unplug the base station before cleaning. Wipe the exterior with a damp cloth. Do not use abrasive cleaners or immerse the base station in water.

7.3 Replacing the Brush Head

It is recommended to replace the toothbrush head every 3 months, or sooner if the bristles become frayed. To replace, firmly pull the old brush head straight off the handle and push a new one on until it clicks into place.

8. Troubleshooting

- **Toothbrush not turning on:** Ensure the toothbrush is fully charged. Place it on the base station for several hours.
- **UV light not activating:** Check that the toothbrush is correctly seated in the base station and the base station is powered on.
- **Reduced brushing power:** The battery may be low. Charge the toothbrush fully. If the issue persists, the brush head may need replacement.
- **Unusual noise:** Ensure the brush head is securely attached. If noise persists, contact customer support.

9. Specifications

Feature	Detail
Model Number	BA-0130-BLK-US
Product Dimensions	4 x 4 x 8 inches
Item Weight	1.26 Pounds
Power Source	Corded Electric
Battery Type	1 Lithium Ion battery (included)
Special Feature	UV sanitizing docking base station
UPC	859332005524

10. Warranty and Support

Tao Clean products come with a **two-year warranty** covering manufacturer-related defects. For warranty claims or technical support, please refer to the contact information provided on the product packaging or visit the official Tao Clean website.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

Why does the light on the front of my Tao Clean BA-0130 UV Sanitizing Station keep flashing white and won't turn off unless unplugged?

A flashing white light on the front of the Tao Clean BA-0130 base station typically indicates that the unit is in its active cleaning and drying cycle. Under normal operation, the base station initiates...

[Read full answer](#)