

Vitaciti 600685985144

Vitaciti Ionic Detox Foot Bath Machine User Manual

Model: 600685985144

1. INTRODUCTION

Thank you for choosing the Vitaciti Ionic Detox Foot Bath Machine. This system is designed to provide a relaxing foot spa experience. This manual provides essential information for the safe and effective operation, maintenance, and troubleshooting of your device. Please read it thoroughly before first use and retain it for future reference.



Image: A person relaxing while using the Vitaciti Ionic Detox Foot Bath Machine.

2. IMPORTANT SAFETY INFORMATION

Please observe the following safety precautions to prevent injury or damage to the device:

- Keep out of reach of children.
- Do not use near water or touch the device with wet hands.

- Do not use while operating a vehicle or heavy machinery.
- Seek medical help immediately if any part of the device is ingested.
- Avoid contact with skin and eyes if any internal components are exposed.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.
- Ensure the output voltage is DC 15V for safe and effective operation.

3. PACKAGE CONTENTS

Verify that all items are present in your package:

- 1x Professional Ionic Foot Bath Machine (Main Unit)
- 1x Professional Ergonomic Foot Massage Basin
- 2x Detox Arrays
- 1x Wrist Strap
- 1x Far Infrared Waist Belt
- 1x Power Adapter
- 1x Instruction Manual (this document)

DO DETOX FOOT SPA

30 MINUTES

TWICE A WEEK



Image: All components included in the Vitaciti Ionic Detox Foot Bath Machine package, laid out for inspection.

4. SETUP INSTRUCTIONS

Follow these steps for initial setup:

1. Place the foot massage basin on a flat, stable surface.
2. Insert a disposable liner into the basin for hygiene.
3. Fill the basin with warm water, ensuring it covers your feet comfortably.
4. Place one detox array into the water, ensuring at least two-thirds of it is submerged.
5. Connect the detox array cable to the main unit.
6. Connect the wrist strap to the main unit and secure it comfortably around your wrist.
7. If using, connect the Far Infrared Waist Belt to the main unit and position it around your waist.
8. Plug the power adapter into the main unit and then into a power outlet.

CLEAR FUNCTIONAL DIVISION



Image: A close-up of a hand connecting the wrist strap to the main control unit.

HEATING RELAXES THE WAIST



Image: A woman using the foot bath, with the infrared waist belt wrapped around her waist.

5. OPERATING INSTRUCTIONS

Once the device is set up, follow these steps to begin your foot detox session:

1. Turn on the main unit using the power button.
2. Immerse your feet into the water in the basin.
3. Add approximately 1/8 teaspoon (about 0.7 oz) of salt to the water. This helps to improve conductivity.
4. Adjust the mode and foot bath time according to your preference using the control panel. A typical session lasts 30 minutes.
5. Relax and enjoy your foot spa session. The unit will automatically shut off after the set time.
6. After the session, remove your feet, discard the liner, and clean the basin and array.



Image: A detailed view of the main unit's control panel, showing the LED display and various buttons for settings.



Image: A woman's feet submerged in the foot bath, with the detox array visible in the center.

Official Operation Video

Your browser does not support the video tag.

Video: An official demonstration of how to operate the Vitaciti Ionic Detox Foot Bath Machine, showing the setup and usage process.

Your browser does not support the video tag.

Video: A brief official video showcasing the Vitaciti Ionic Detox Foot Bath Machine in action, highlighting its ease of use.

6. MAINTENANCE

Cleaning

- After each use, discard the disposable liner from the foot basin.
- Rinse the foot basin with clean water and wipe it dry.
- Clean the detox array with a brush and warm water to remove any residue. Ensure it is completely dry before storage.
- Wipe the main unit, wrist strap, and waist belt with a soft, damp cloth. Do not submerge the main unit in water.

Storage

- Store the device and all accessories in a cool, dry place, away from direct sunlight and extreme temperatures.
- Ensure all components are clean and dry before storing to prevent mold or damage.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not power on	Power adapter not connected or faulty; power outlet issue.	Ensure power adapter is securely connected. Try a different power outlet.
"Current error" displayed on screen or flickering screen	Poor connection of array/wrist strap; insufficient salt in water; faulty array.	Check all connections. Add more salt to the water (up to 1/8 tsp). Ensure the array is properly submerged. If the issue persists, the array may need replacement.
Water does not change color or array does not fizz	Insufficient salt in water; faulty array; poor connection.	Add 1/8 tsp of salt to the water. Ensure the array is correctly connected and submerged. Test with a different array if available.
Wrist strap is finicky	Loose connection; strap not making good contact with skin.	Ensure the wrist strap is securely plugged into the unit and fits snugly on your wrist.

8. SPECIFICATIONS

Feature	Detail
Model Name	Vitaciti Ionic Detox Foot Bath System with Infrared
Model Number	600685985144
Product Dimensions	15.55 x 5.51 x 19.69 inches
Item Weight	8.82 Pounds (4 Kilograms)
Power Source	Corded Electric
Output Voltage	DC 15V
Material Type	Plastic

Feature	Detail
Tub Capacity	3.2 gallons (12 liters)
Tub Dimensions	20"L x 16"W x 6"H
Water Resistance Level	Not Water Resistant

9. WARRANTY AND SUPPORT

The Vitaciti Ionic Detox Foot Bath Machine comes with **aone-year warranty**. For any technical support or warranty claims, please contact VitaCiti Store directly. Free technical support is available to assist you with any questions or issues you may encounter.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question