



Manuals.plus /

› iRestore /

› iRestore Professional Laser Red Light Therapy for Hair Growth System - Instruction Manual

iRestore B07M5X4D3P

iRestore Professional Laser Hair Growth System

Model: B07M5X4D3P | Brand: iRestore

INTRODUCTION

This instruction manual provides essential information for the safe and effective use of your iRestore Professional Laser Hair Growth System. Please read this manual thoroughly before operating the device and retain it for future reference. The iRestore Professional is an FDA-cleared device designed to treat androgenic alopecia and promote hair growth in men and women.

SAFETY INFORMATION

The iRestore device is FDA 510(k) cleared for laser hair regrowth. Users must be at least 18 years of age. Please refer to the Norwood-Hamilton Scale (approved: Class IIa to V) / Ludwig-Savin Scale (approved: Class I to II) and Fitzpatrick Skin Type Scale (approved: Type I to IV) for compatibility. This device is not indicated for users with severe hair loss. iRestore has been clinically tested for males between the ages of 18 and 48 and females between the ages of 18 and 60. While not clinically tested outside this range, it does not indicate ineffectiveness. If you are taking medications that affect skin sensitivity to light or have a medical condition that makes you extremely sensitive to light, consult your doctor before use. This device is indicated to treat androgenic alopecia.

PRODUCT OVERVIEW

The iRestore Professional Laser Hair Growth System is a helmet-style device designed for convenient at-home use. It utilizes a combination of 282 medical-grade lasers and LEDs to deliver low-level light therapy (LLLT) to the scalp.



Image: The iRestore Professional helmet and its wired remote control.

Key Features:

- **Full Scalp Coverage:** Features 282 lasers and LEDs for comprehensive treatment of the hairline, sides, and crown.
- **Hands-Free Operation:** Designed for ease of use, allowing users to perform other activities during sessions.
- **Comfortable Design:** Lightweight, ventilated design with medical-grade interior padding for a secure and comfortable fit.
- **Clinically Studied:** A drug-free solution for supporting thicker-looking hair.
- **Flexible Integration:** Can be used independently or in conjunction with other hair care products.

SETUP

1. **Unboxing:** Carefully remove the iRestore Professional device and all accessories from its packaging.
2. **Inspect Components:** Verify that all components are present and undamaged: the helmet, power

adapter, and user manual.

3. **Connect Power:** Plug the power adapter into the helmet's port and then into a standard electrical outlet. Ensure a secure connection.
4. **Optional Battery Pack:** If using the optional rechargeable battery pack (sold separately), connect it to the device as per its specific instructions for cordless operation.



Image: The iRestore Professional helmet with its power cable connected, ready for use.

OPERATING INSTRUCTIONS

1. **Position the Device:** Place the iRestore Professional helmet comfortably on your head, ensuring it covers the entire scalp area requiring treatment. The medical-grade interior padding should provide a snug fit.
2. **Start Treatment:** Press the 'Power' button on the remote control. The device will automatically begin a 25-minute treatment session. The timer will be displayed on the remote.
3. **During Treatment:** The device is hands-free, allowing you to relax, read, or perform light activities.
4. **Session Completion:** The device will automatically turn off after 25 minutes.

5. **Frequency:** Use the device for 25 minutes, 2-3 times per week. Consistency is key for optimal results.

Increase hair count by 43.2% in just 4 months!*



5.7% PLACEBO SYSTEM

43.2% REAL SYSTEM

* BASED ON A DOUBLE-BLIND CLINICAL STUDY CONDUCTED BY DR. BODIAN, A BOARD-CERTIFIED DERMATOLOGIST AND LASER SURGEON, WITH 40 SUBJECTS TO EVALUATE THE SAFETY AND EFFICACY OF THE iRESTORE PROFESSIONAL. INDIVIDUAL RESULTS MAY VARY.

Image: A user wearing the iRestore Professional helmet, illustrating its hands-free design.

Expected Results:

Visible results can typically be observed within 3 to 6 months of consistent use. Individual results may vary.

MAINTENANCE

- **Cleaning:** Gently wipe the exterior of the helmet with a soft, dry cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in its original packaging or a safe, dry place away from direct sunlight and extreme temperatures when not in use.
- **Cable Care:** Avoid kinking or bending the power cable excessively to prevent damage.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Power cable not properly connected or no power from outlet.	Ensure the power cable is securely plugged into both the device and a working electrical outlet. Check if the outlet is functional.
Red lights are not visible.	Device may not be properly seated on the head, or safety sensors are not activated.	Adjust the helmet to ensure it sits firmly on your head, activating the internal safety sensors. The device will not emit light if sensors are not engaged.
Remote control display is blank.	Power issue or remote control malfunction.	Verify power connection to the helmet. If the issue persists, contact customer support.

SPECIFICATIONS

- **Model:** B07M5X4D3P
- **Brand:** iRestore
- **Light Source:** 282 Medical-Grade Lasers & LEDs
- **Power Output:** 1,410 mW Total Output
- **Treatment Time:** 25 minutes per session
- **Recommended Use:** 2-3 times per week
- **Product Dimensions:** 10.8 x 9.9 x 6.6 inches
- **Weight:** 2.8 Pounds
- **UPC:** 859690003675
- **Manufacturer:** iRestore

OFFICIAL PRODUCT VIDEOS

iRestore Professional Hair Growth System - How It Works

Your browser does not support the video tag.

Video: An official iRestore video explaining the operational principles and benefits of the Professional Hair Growth System.

The iRestore Professional Laser Hair Growth System

Your browser does not support the video tag.

Video: An official iRestore overview of the Professional Laser Hair Growth System, highlighting its features and ease of use.

WARRANTY AND SUPPORT

The iRestore device is protected by a 1-year manufacturer's warranty against defects. For warranty claims, technical assistance, or any questions regarding your iRestore Professional Laser Hair Growth System, please contact iRestore customer support. Refer to the product packaging or the official iRestore website for the most current contact information.

Additionally, iRestore offers a 6-month money-back guarantee. If you are not satisfied for any reason, you may return the device within 6 months of the purchase date for a refund. (Service fees may apply, refer to terms and conditions).